

NITRYLEX® BASIC - NITRILE GLOVES POWDER FREE

Non-sterile, powder free, examination and protective glove for single use. Optimal sensitivity to touch. **EN 455, EN 374**

GIMA code	NITRILE GLOVES NON STERILE	Size	Colour	Minimum order
21801	Small	6-7	navy blue	box of 100
21802	Medium	7-8	navy blue	box of 100
21803	Large	8-9	navy blue	box of 100
21804	Extra large	9-10	navy blue	box of 100

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.



21802

NITRYLEX® CLASSIC - NITRILE GLOVES POWDER FREE

Non-sterile, powder free, examination and protective glove for single use. Greater sensitivity to touch. **EN 455, EN 374**

GIMA code	NITRILE GLOVES NON STERILE	Size	Colour	Minimum order
21806	Small	6-7	blue	box of 100
21807	Medium	7-8	blue	box of 100
21808	Large	8-9	blue	box of 100
21809	Extra large	9-10	blue	box of 100

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.



21807

SAFE PACK SYSTEM FOR NITRYLEX® CLASSIC ONE BY ONE

Safe Pack dispenser system made of plastic, for the distribution of gloves. It can be used only with Nitrylex® Classic One by One, powder free examination gloves (box of 200). **EN 455, EN 374**

GIMA code	NITRILE GLOVES NON STERILE	Size	Colour	Minimum order
21811	Small	6-7	blue	box of 200
21812	Medium	7-8	blue	box of 200
21813	Large	8-9	blue	box of 200
21814	Extra large	9-10	blue	box of 200

21815 Safe Pack wall dispenser system - plastic

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.

MINIMIZE RISK OF CROSS-CONTAMINATION

There is no contact between the user's hand with the working surface of the glove or with the box.



21815



Special folding method allows extraction of a single glove from the cuff

NITRYLEX® BLACK - NITRILE GLOVES POWDER FREE

Non-sterile, powder free, examination and protective glove for single use. Opaque finish. **EN 455, EN 374**

GIMA code	NITRILE GLOVES NON STERILE	Size	Colour	Minimum order
21816	Small	6-7	black	box of 100
21817	Medium	7-8	black	box of 100
21818	Large	8-9	black	box of 100
21819	Extra large	9-10	black	box of 100

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.



21817

NITRYLEX® AMBULANCE HIGH RISK - NITRILE GLOVES POWDER FREE

Non-sterile, powder free, examination and protective glove for single use. Thicker for extra safety in demanding work conditions. **EN 455, EN 374**

GIMA code	NITRILE GLOVES NON STERILE	Size	Colour	Minimum order
21821	Small	6-7	orange	box of 100
21822	Medium	7-8	orange	box of 100
21823	Large	8-9	orange	box of 100
21824	Extra large	9-10	orange	box of 100

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.



Extended cuff for additional protection

EN 455, EN 374



21822

TECHNICAL SPECIFICATIONS NITRYLEX®

	Basic	Classic	Classic One by One	Black	Orange
External surface:	fingertips textured	micro + fingertips textured, polymerized	micro + fingertips textured, polymerized	micro + fingertips textured, polymerized	micro + fingertips textured, polymerized
Powder content:	< 2mg/glove	< 2mg/glove	< 2mg/glove	< 2mg/glove	< 2mg/glove
Cuff:	beaded	beaded	beaded	beaded	beaded
AQL:	≤ 1.5	≤ 1.0	≤ 1.0	≤ 1.0	≤ 1.0
MD (Medical Device):	Class I	Class I	Class I	Class I	Class I
PPE (Personal Protective Equipment):	Category III	Category III	Category III	Category III	Category III