

X-LIFE AUTOMATIC DIGITAL B.P. MONITORS - ESH VALIDATED

Size: 140x104x50 mm
Weight: 247 g


LCD size 60x51 mm

Cuff size for arm
circumferences
23-33 cm



32924

Supplied with
carrying bag



- **32924 X-LIFE B.P. MONITOR - arm**
 - **32926 X-LIFE B.P. MONITOR - wrist**
- Accurate and precise automatic digital blood pressure meters for smooth and fast measurements with innovative functions:
- **IHB:** irregular heartbeat detector
 - **WHO risk category indicator** (classifies your measurements in 6 stages of risk, from optimal to severe hypertension)
 - proven clinically to meet European Society of Hypertension ESH protocol
 - easy to read, wide angle display
 - 120 memories for 3 users: 3 zones of 40 memories each
 - start/stop unique button, with blue backlight



European Society
Hypertension

Size: 74.5x74x30 mm

Weight: 95 g

LCD size 34.5x42 mm

Cuff size for wrist
circumferences
13.5-19.5 cm



32926



Liquid crystal
display

- multilingual manual and box: GB, FR, IT, ES
- memory average: automatically average last 3 memories
- **32925 SPARE CUFF for 32924 - adult**
- **32927 AC ADAPTOR for 32924**

TECHNICAL SPECIFICATIONS

Measuring method:	oscillometric
Measuring range:	pressure: 0-300 mmHg pulse: 40-199 bpm
Accuracy:	pressure: ± 3 mmHg pulse: $\pm 5\%$
Power supply:	arm type: 4xAAA alkaline batteries wrist type: 2xAAA alkaline batteries
Battery life:	arm type: 300 measures wrist type: 250 measures

JOLLY DIGITAL B.P. MONITOR



European Society
Hypertension



32777



- **32777 JOLLY DIGITAL BLOOD PRESSURE MONITOR**
- Professional Blood Pressure monitor for medical use.
- large easy-to-read illuminated LCD display 86x74 mm
 - WHO classification
 - arrhythmia detection (IHB)
 - 4-person memory zone - 30 memory space each
 - average of all saved measured values and morning and evening blood pressure of the last 7 days
 - time and date
 - automatic switch-off
- Multilingual manual and box:
GB, FR, IT, ES, PT, DE, PL, GR.
- **32922 CUFF - adult - 22-30 cm (8-11") - black - spare**
 - **32923 CUFF - adult large - 30-42 cm (11-16") - grey**
 - **32920 CUFF - adult extra large - 42-48 cm (16-18.2") - grey**



TECHNICAL SPECIFICATIONS

Measurement system:	Oscillometric method
Measurement range:	Pressure: 0-300 mmHg ± 3 mmHg Pulse: 40-180 pulse/min. $\pm 5\%$
Power source:	4x1.5 V AAA batteries included
Outside dimension:	150x95xh 65 mm
Weight:	approx 235 g (without batteries)

DOMINO DIGITAL B.P. MONITOR

Automatic inflation
deflation



32803



- **32803 DOMINO DIGITAL BLOOD PRESSURE MONITOR**
 - **32804 DOMINO DIGITAL B.P.M. on trolley with basket**
- Professional Blood Pressure monitor for hospital and medical use.
- extra large LCD display with innovative design
 - date/time/ambient temperature ($^{\circ}\text{C}/^{\circ}\text{F}$) display
 - IHB Irregular heartbeat detector
 - WHO blood pressure classification indicator
 - 120 memories (60 sets for 2 users)
 - alarm clock reminding function
- Multilingual manual and box:
GB, FR, IT, ES, PT, DE, GR, Arabic
- **32802 AC ADAPTOR 230 V - 50 Hz**
 - **32809 ADULT CUFF - spare**
 - **32865 COILED TUBING**
 - EXTENSION 3 m (42/45 spirals)**



32804

Measuring range:	Pressure - 40-260 mmHg Pulse: 40-160 beats/min.
Accuracy:	Pressure ± 3 mmHg Pulse: $\pm 5\%$
Power supply:	4x "AA" batteries or AC adaptor (option)
External dimension:	158x120xh 127 mm
Weight:	490 g

