

LATEX GLOVES - PRE POWDERED OR POWDER FREE



LATEX GLOVES - ambidextrous - with hem
 Pure latex gloves non sterile.
 Sensitive and resistant, ensure optimal grip and protection.
 Suitable for dressing in hospital, ambulatory and laboratory.
 Packed in dispenser of 100. Medical device Class I (all gloves).
 PPE cat III (only powder free gloves).
EN455-EN374-AQL 1.5

Pre powdered	Powder free	LATEX GLOVES NON STERILE (box of 100)
-	25409	Extra small size - non sterile
25706	25711	Small size - non sterile
25707	25712	Medium size - non sterile
25708	25713	Large size - non sterile
25408	25410	Extra large size - non sterile

LATEX SURGICAL GLOVES - STERILE



LATEX STERILE SURGICAL GLOVES
 The micro textured surface assures a very good grip in any use. High protection and anatomic shape assure a very good comfort. Sterilized by EtO.
 Length 270 mm. Conform to USP: international pharmacopoeia.
 Packed in pairs, dispenser of 50 pairs.
 MD Class II A. PPE cat III.
EN455-EN 374-AQL 0.65

Pre powdered	Powder free	LATEX SURGICAL GLOVES STERILE - box of 50 pairs
25685	25695	Size 6.5 - sterile
25686	25696	Size 7 - sterile
25687	25697	Size 7.5 - sterile
25688	25698	Size 8 - sterile
25689	25699	Size 8.5 - sterile

COPOLYMER GLOVES



COPOLYMER GLOVES ON PAPER
 Examination gloves, elasticized, for gynaecological use. MD class 1. AQL 1.5
 Sterilized by EtO (only 25709, 25710)

Non sterile box of 1,000*	Sterile** box of 100	COPOLYMER GLOVES
25714	25709	Medium
25720	25710	Large

*10 packs of 100 **individually packed

SANTEX® AND DERMAGEL® LATEX GLOVES

SANTEX® POWDERED - ambidextrous
 Non-sterile, powdered, examination and protective glove for single use. Made of natural rubber latex. Smooth external surface guarantees a tactile sensitivity during the work. They minimize hand sweating. Enhanced tactile sensitivity. **AQL≤1.0**



Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.

DERMAGEL® COATED - ambidextrous
 Non-sterile, powder free, examination and protective gloves for single use. Made of natural rubber latex. The polymer layer inside ensures easy donning even on wet hands. **AQL≤1.0**



Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.

TECHNICAL SPECIFICATIONS	SANTEX®	DERMAGEL®
	External surface:	Smooth
Cuff:	— beaded —	—
MD (Medical Device):	— Class I —	—
PPE (Personal protective equipment):	— Category III —	—

Santex® Pre powdered	Dermagel® Powder free	LATEX GLOVES NON STERILE (box of 100)
21851	21856	Small size - non sterile
21852	21857	Medium size - non sterile
21853	21858	Large size - non sterile
21854	21859	Extra large size - non sterile