

## SIRIO SPHYGMO AND 3 MULTICUFF KIT WITH BAG

CE 0476



32729

- **32729 SIRIO SPHYGMO - adult**  
Classic sphygmomanometer with fabric cuff with velcro closure (green/grey colour). Reliable large dial manometer (60 mm) for easy readings, fixed on bulb. Scale 0-300 mmHg, graduated every 2 mmHg. Supplied in carrying case and box, with multilingual user manual: GB, FR, IT, ES.



- **32705 SIRIO PEDIATRIC KIT - sphygmo with 3 pediatric cuffs**

Sirio with 3 pediatric cuffs: minicuff (32812) 4 months - 3 years, minicuff (32813) 3-7 years, minicuff (32814) 7-10 years in a nylon bag.

- **27177 SIRIO MULTICUFF PEDIATRIC BAG - empty**



32705

- **32706 SIRIO KIT 3 - sphygmo with 3 one piece cuffs**

Sirio with 3 one piece cuffs: child (32910), adult (32911) and adult large (32912) in a nylon bag.

- **27178 SIRIO MULTICUFF BAG 3 - empty**



32706

- **32707 SIRIO KIT 5 - sphygmo with 5 cuffs**

Sirio with 5 cuffs: infant cuff (32886), child cuff (32888), adult cuff (32830), adult large cuff (32834) thigh cuff (32837) in a nylon bag.

- **27179 SIRIO MULTICUFF BAG 5 - empty**



32707

## LONDON: PROFESSIONAL SPHYGMO

CE 0476



32725

- **32725 LONDON SPHYGMO - adult - black cuff**
  - **32726 LONDON SPHYGMO - adult - red cuff**
- Professional, stylish, reliable, lightweight and latex free sphygmomanometer with ABS manometer, integral PVC bulb, high quality precision mechanism, velcro closure calibrated nylon cuff and one tube TPU bladder. Dial Ø 55 mm, with specially designed eye-catching colour printing and easy to read. Special outer ring protects the manometer from shocks and the glass from scratches. Supplied in a nice colour co-



32726

ordinated carrying bag with multilingual box and user manual: GB, FR, IT, ES, DE, PT, GR.



32726

