

RIESTER ORIGINAL CUFFS - LATEX FREE

RIESTER VELCRO CUFFS WITH BLADDER 1 OR 2 TUBES

Child, adult, thigh cuffs made of high quality materials.

| Child, adult, thigh cuffs made of high quality materials. | | | | | | |
|---|--------------|-------------------------|------------|-----------|------------|---|
| GIMA code | RIESTER code | RIESTER CUFFS 1 TUBE | Colour | Size cm | Circum. cm | Compatibility |
| 49811 | 129 | Child velcro cuff | light blue | 35.5x10 | 13 - 17 | All Riester single tube manual aneroid sphygmomanometers* |
| 49810 | 120 | Adult small velcro cuff | black | 42x13 | 17 - 24 | |
| 49812 | 107 | Adult velcro cuff | black | 54.5x14.5 | 24 - 34 | |
| 49813 | 122 | Adult large velcro cuff | black | 70x15 | 34 - 44 | |
| 49814 | 125 | Thigh velcro cuff | black | 70x22 | 42 - 50 | |
| RIESTER CUFFS 2 TUBES | | | | | | |
| 49821 | 130 | Child velcro cuff | light blue | 35.5x10 | 13 - 17 | All Riester twin tube sphygmomanometers including BigBen |
| 49820 | 121 | Adult small velcro cuff | black | 42x13 | 17 - 24 | |
| 49822 | 108 | Adult velcro cuff | black | 54.5x14.5 | 24 - 34 | |
| 49823 | 123 | Adult large velcro cuff | black | 70x15 | 34 - 44 | |
| 49824 | 126 | Thigh velcro cuff | black | 70x22 | 42 - 50 | |

*Metal valve (included in Minimus II, E-Mega, Babyphon cuff) not included with spare cuff, can use valve from old cuff

RIESTER ONE PIECE CUFFS 1 OR 2 TUBES

They prevent cross contaminations. Comfortable soft surface, latex free. Washable at 60°C and disinfected in liquid, reusable up to 10,000 times. ISO 105 E01 tested.

| GIMA code | RIESTER code | ONE PIECE RIESTER CUFFS 1 TUBE - length 50 cm | Colour | Size cm | Circum. cm | Compatibility |
|--|--------------|---|------------|-----------|------------|-----------------------|
| 49816 | 154 | Child one-piece cuff | dark green | 34.5x10 | 13 - 17 | Same as cuffs 1 tube |
| 49817 | 150 | Adult one-piece cuff | dark green | 54x14.5 | 24 - 32 | |
| 49818 | 152 | Adult large one-piece cuff | dark green | 63x16 | 32 - 41 | |
| ONE PIECE RIESTER CUFFS 2 TUBES - length 50 cm | | | | | | |
| 49827 | 151 | Adult one-piece cuff | dark green | 54.5x14.5 | 24 - 32 | Same as cuffs 2 tubes |
| 49828 | 153 | Adult large one-piece cuff | dark green | 63x16 | 32 - 41 | |




HEINE ORIGINAL CUFFS - LATEX FREE

HEINE ORIGINAL CUFFS 1 TUBE

High quality single tube cuffs with velcro fastener. Perfectly balanced for comfort and ease of use. Compatible with all Heine sphygmomanometers. Water resistant, washable material. Latex free.

| GIMA code | HEINE code | HEINE CUFFS 1 TUBE | Colour | Size cm | Circum. cm | Compatibility |
|-----------|--------------|--------------------|--------|---------|------------|------------------|
| 49801 | M-000-09-702 | Child - 2 | grey | 10x34 | 13-20 | All Heine 1 tube |
| 49802 | M-000-09-704 | Adult small - 3 | grey | 14x53 | 20-29 | |
| 49803 | M-000-09-619 | Adult - 4 | grey | 14x58 | 29-41 | |
| 49804 | M-000-09-706 | Adult large - 5 | grey | 17x70 | 35-47 | |



ERKA ORIGINAL GREEN CUFFS SUPERB RAPID - LATEX FREE

ERKA GREEN CUFFS SUPERB RAPID 1, 2 OR TWIN TUBES

Superb cuffs are easy to disinfect and clean, equipped with a waterproof surface. Cuffs are comfortable and consist of a cover and a removable inflation bladder. Super soft adhesive and fleece strap on the cuff. Made of durable patented German PU material. Optimal in blood pressure measurement because of pressure distribution and the highly precise values. PVC, DEHP and latex free. Available in different sizes with 1, 2 or twin tubes.

| GIMA code | ERKA code | ERKA GREEN CUFFS SUPERB RAPID 1 TUBE | Colour | Size cm | Circum. cm | Compatibility |
|---|-----------|--------------------------------------|--------|----------|------------|---------------------------|
| 32895 | 025-28493 | Infant | grey | 21.5x6.6 | 10-15 | Switch and Kobold |
| 32896 | 025-27493 | Baby | grey | 30.3x9.1 | 14-21.5 | |
| 32891 | 025-26493 | Child | grey | 40x11.6 | 20.5-28 | |
| 32892 | 025-20493 | Adult | grey | 47x14.3 | 27-35 | |
| 32893 | 025-24493 | Adult large | grey | 56x17.5 | 34-43 | |
| 32894 | 025-29493 | Thigh | grey | 75x21.5 | 42-54 | |
| ERKA GREEN CUFFS SUPERB RAPID 2 TUBES | | | | | | |
| 32898 | 025-26481 | Child | green | 40x11.6 | 20.5-28 | Vario |
| 32899 | 025-20481 | Adult | green | 47x14.3 | 27-35 | |
| 32900 | 025-24481 | Adult large | green | 56x17.5 | 34-43 | |
| ERKA GREEN CUFFS SUPERB RAPID TWIN TUBE | | | | | | |
| 32877 | 025-26482 | Child | green | 40x11.6 | 20.5-28 | Perfect Aneroid and Profi |
| 32878 | 025-20482 | Adult | green | 47x14.3 | 27-35 | |
| 32879 | 025-24482 | Adult large | green | 56x17.5 | 34-43 | |



PRESS^{SAN} HYGIENIC PROTECTION OF SPHYGMOMANOMETER

• 32858 PRESS^{SAN} HYGIENIC PROTECTION OF SPHYGMO - box of 10

A practical hygienic accessory for the sphygmomanometer to protect patient's skin whilst measuring the blood pressure. Made of a white non-woven fabric.

Ensures patient's comfort and hygiene.

Can fit any kind of sphygmomanometer or pressure holter.

MADE IN ITALY

32858

