

## KOBOLD SPHYGMOMANOMETER BY ERKA

- 32691 BLACK KOBOLD with grey cuff - 223.20493
  - 32695 GREY KOBOLD SET with 3 ped. "Superb" cuffs - 227.71.393
  - 32696 GREY KOBOLD SET with 6 "Superb" cuffs - 227.71.195
- Elegant and compact, single tube sphygmomanometer for frequent use in clinic and consulting room. Precision manometer, anodised aluminium Ø 48 mm. Sensitive adjustable air reduction valve, high quality chrome finish. Scale: 300 mmHg, 2 mm scale division. New technology with air intake via the manometer capsule in the specially tempered copper-beryllium membrane. Codes 32693, 32695, 32696 are latex free.

### ERKA "SUPERB" GREEN CUFFS

All ERKA Green cuff products are PVC, demp and latex-free, each is individually tested prior to delivery and guaranteed for 50,000 inflation cycles. They are equipped with clearly marked size and positioning features for the most accurate possible measurements. The soft, easy to disinfect, biocompatible materials ensure maximum patient comfort. Conform to ANSI/AAMI SP10 standard

- 32891 PEDIATRIC "SUPERB" CUFF 1 TUBE 20.5 - 28 cm - grey
- 32892 ADULT "SUPERB" CUFF 1 TUBE 27-35 cm - grey - spare



32691

**ERKA.**

MADE IN GERMANY



32696

CE 0124

Includes 6 cuffs for below indicated circumferences of upper arm

Thigh 42-54 cm    Adult large 34-43 cm    Adult 27-35 cm    Child 20.5-28 cm    Baby 13-20 cm    Infant 10-15 cm



CE 0124

Child 20.5-28 cm    Baby 14-21.5 cm    Infant 10-15 cm

32695

## PERFECT ANEROID

## ERKA SWITCH

## ERKA SWITCH SIMPLEX



32688

**ERKA.**

MADE IN GERMANY

- 32688 ERKA PERFECT - 201.204.82 dial Ø 48 mm - 2 tubes green
- Sphygmomanometer with screw valve for sensitive adjustments, velcro cuff and case. - manometer housing brass chromed plated - membrane-protection 2-tube technology - measuring range 0-300 mmHg ±1% 3 year warranty.
- 32877 PEDIATRIC "SUPERB" CUFF TWIN TUBE 20.5 - 28 cm - grey
  - 32878 ADULT "SUPERB" CUFF TWIN TUBE 27-35 cm - grey - spare



32698

**ERKA.**

MADE IN GERMANY

- 32698 ERKA SWITCH 2.0 COMFORT - dial Ø 56 mm - 1 tube grey - 293.204.93
- It offers outstanding ergonomics and simple adaptation for left and right-handle users. It is manufactured from shock-resistant polycarbonate and has passed the shock test according to AAMI SP9. The sight glass is chemically stable and unbreakable. The brass valve screw works smoothly and is adjustable to exact settings. MR compatibility. Measuring range 0-300 mmHg ±1%. Supplied with pouch. 4 year warranty hand held, 2 year manometer.
- 32891 PEDIATRIC "SUPERB" CUFF 1 TUBE 20.5 - 28 cm - grey
  - 32892 ADULT "SUPERB" CUFF 1 TUBE 27-35 cm - grey - spare



32697

**ERKA.**

MADE IN GERMANY

- 32697 ERKA SWITCH SIMPLEX - dial Ø 56 mm - twin tubes grey - 294.408.94
- It is manufactured from shock-resistant polycarbonate and has passed the shock test according to AAMI SP9. The sight glass is chemically stable and unbreakable. The brass valve screw works smoothly and is adjustable to exact settings. Measuring range 0-300 mmHg ±1%. Black velcro cuff. Supplied with pouch. Magnetic Resonance (MR) compatibility. 2 year warranty.