


GLOVE DISPENSER - STAINLESS STEEL



- **25431 GLOVE DISPENSER - single - stainless steel**
 - **25432 GLOVE DISPENSER - triple - stainless steel**
- Hold 1 or 3 standard boxes of 100 gloves.
Made of AISI 304 stainless steel. Autoclavable at 121 to 134 °C.
Can be easily hung to the wall and into ambulances.
Screw included.

GLOVE DISPENSER - PLASTIC



- **25430 GLOVE DISPENSER - plastic**
- Holds one box of 100 gloves.
Made of unbreakable washable thermoplastic material.
Fixed base and removable cover once the box of gloves needs to be replaced. Can be easily hung to the wall and into ambulance in vertical or horizontal position.
Compatible with box of gloves up to 25.5x7xh 12 cm

DRACULA

- **25700 DRACULA**
- Universal holder for gloves and tissue boxes, with little sharp teeth which automatically adjust to tightly hold any size of box.
One-hand operation.



For minimum 500 pcs we can put your logo here

COTTON GLOVES

100% cotton gloves. Extra long and coloured cuff for easy identification of the size.

GIMA code	Size	Cuff colour	Length cm		Weight g		Minimum order
			min	max	min	max	
25864	6.5	red	20	22	18	20	box of 10 pairs
25865	7	blue	21	23	19	21	box of 10 pairs
25866	7.5	green	23	25	21	23	box of 10 pairs
25867	8	white	24	26	23	25	box of 10 pairs
25868	8.5	orange	25	27	25	27	box of 10 pairs




LATEX FINGER COTS

- **25941 LATEX FINGER COTS - medium**
 - **25942 LATEX FINGER COTS - large**
- Latex finger cots to cover 1 finger.
Non sterile, single use, for finger protection, with talc.
minimum order: 1.000 pcs. (10 bags of 100 pcs.)




NON-WOVEN UNIFORM AND PANTS



- **26076 UNIFORM COAT+PANT - S - box of 50**
- **26077 UNIFORM COAT+PANT - M - box of 50**
- **26078 UNIFORM COAT+PANT - L - box of 50**
- **26079 UNIFORM COAT+PANT - XL - box of 50**



COLONOSCOPY PANTS



- **21030 COLONOSCOPY PANTS - box of 50**
- Colonscopy pants, one size, non sterile.
Elastic waist, rear opening with velcro.
Made of water-repellent non-woven SMS 35g/m², dark blue, non transparent.

DISPOSABLE KITS FOR PATIENT/VISITOR DRESSING - NON STERILE

- **26017 KIMONO - blue - box of 100**
- Opaque non-woven SMS kimono 30 g/m² with short-sleeves. Front opening and closing at the waist with laces.
One size: length 120 cm, width 150 cm.



- **26018 PATIENT DRESSING KIT - box of 25**
- Kit includes kimono, round cap with elastic band and shoe-covers.



- **26019 VISITOR DRESSING KIT - box of 25**
- Kit includes gown for visitor, round cap with elastic band, three layers mask with ear loops, polyethylene shoe-covers.