

AUDIOMETRIC BOOTHS - CONFORM TO EU 2017/745 - Class I



- **33648 NEW PRO 25 - 96x96x197 cm**
 - **33650 PRO 25 - 96x96x197 cm**
- Portable, lightweight, sound proof booth, easy to install in few minutes. Interior finishing made of acoustic pyramid foam in anthracite grey (33650) or made of sound-absorbing polyurethane, coated with acoustic lining washable and sanitizable (33648). External finishing made of highly resistant, fire proof washable panel. Module for cable connection included.
- **33651 TABLE - 79x50 cm**

- **33649 PRO 28F - 96x96x197 cm**
- New compact booth with an outstanding level of soundproofing and excellent acoustic thanks to innovative materials. Internal finishing made of fire proof, washable, sanitizable and very resistant fabric. Internal ventilation system. External finishing made of highly resistant, fireproof, washable panel. Module for cable connection included.
- **33651 TABLE 79x50 cm**

- **33654 PRO 30 - 107x107x244** cm**
 - **33653 PRO 30 - 132x132x244** cm**
 - **33655 PRO 30 - 211x107x244** cm**
 - **33656 PRO 30 - 211x211x244** cm**
 - **33657 PRO 30 - size on request**
- available any size from 1x1 to 3x3
- Interior finishing made of auto extinguishing acoustic pyramid foam in anthracite grey. Range with different sizes, with high noise reduction and ventilation system. External finishing made of highly resistant, fireproof, washable panel and aluminium frame. Module for cable connection included. Sanitizing air-exchanging system on request.
- **33651 TABLE 79x50 cm**
 - **33652 RAMPS**
- for wheelchair access

MADE IN ITALY

ACOUSTIC INSULATION								
Hertz	125	250	500	1000	2000	4000	8000	
PRO 25 NEW	dB 19.1	18.7	30.8	31.5	27.9	35.2	43.3	
PRO 25	dB 18.0	16.0	26.7	30.8	28.8	30.9	39.0	
PRO 28	dB 21.7	25.6	32.6	36.2	34.0	41.1	43.0	
PRO 30	dB 23.4	24.1	27.8	34.5	31.3	36.4	37.7	

STANDARD ACCESSORIES

- module for silent audiometer cable connection
- internal illumination
- ventilation system (PRO 28, PRO 30 only)
- ISPM-15 wooden box for transport
- manual: GB, IT

TECHNICAL SPECIFICATIONS							
	33648	33650	33649	33654	33653	33655	33656
External size (cm):	96x96xh 197	96x96xh 197	96x96xh 197	107x107xh 244**	132x132xh 244**	211x107xh 244**	211x211xh 244**
Internal size (cm):	85x85xh 180	81x81xh 180	81x81xh 173.5	92x92xh 201	117x117xh 201	196x92xh 201	196x196xh 201
Door size cm:	55xh 168	55xh 168	55xh 168	75xh 200	75xh 200	75xh 200	75xh 200
Door opening:	Left/Right	Left/Right	Left/Right	Right*			
Window size cm:	52xh 52	52xh 52	52xh 52	56xh 88			
Window material:	2 - crystal insulating safety glass						
Wall thickness cm:	5.5	7.5	7.5	7.5	7.5	7.5	7.5
Weight (kg):	198	198	233	285	355	435	645
Ventilation:	no	no	yes	yes	yes	yes	yes
Tested noise reduction dB at 1,000 Hz:	31.5	30.8	36.2	34.5			
Internal illumination power V-Hz:	230 - 50/60			86-265 - 50/60			
Norms:	Class I Medical Devices European Directive EU 2017/745						
	*Left opening available on request				**Height includes air ventilation system (27 cm)		

UNIVERSAL NOSE MASKS AND FACE MASKS FOR CPAP - SINGLE PATIENT

NOSE MASKS AND FACE MASKS FOR CPAP

Airtight, light and soft masks adjust automatically according to different face shapes. Honeycomb vents designed to be gentle and quiet. Quick-release connector helps separate the tubing from the mask easily. Headgear is made of soft and skin-friendly fabric to ensure maximum comfort. Masks are suitable for Yuwell, Resmed, Philips, Fisher & Paykel and most of CPAP machines in the market.

Therapy pressure range: 4-30 cm H₂O. Tube connection: 22 mm

GIMA code	BREATHING NOSE MASKS
53710	Breathing nose mask - small
53705	Breathing nose mask - medium
53712	Breathing nose mask - large
BREATHING FACE MASKS	
53713	Breathing face mask - small
53714	Breathing face mask - medium
53715	Breathing face mask - large

Breathing nose masks

Breathing face masks

