

DISPOSABLE AND RE-USABLE TURBINES FOR MIR SPIROMETERS

• 33507 DISPOSABLE TURBINE WITH INTEGRATED MOUTHPIECE - box of 60

Single patient use. Low cost and high accuracy. Individually packed to guarantee maximum hygiene. Easy to replace, no calibration needed.

Compatible with GIMA codes: 33509, 33511-14, 33518-22, 33527-29, 33530, 33532, 33534.

• 33526 RE-USABLE TURBINE - bayonet

Easy to replace, no calibration needed. Compatible with GIMA codes: 33510-14, 33518, 33520-22, 33528-32, 33534.

33526




33507



MOUTHPIECES FOR MOST SPIROMETER BRANDS

Paper mouthpieces (disposable) and plastic mouthpieces (re-usable) for most brands of spirometers and peak flow. Paper mouthpieces are packed in individual transparent pouch.



MADE IN EUROPE

33432

GIMA PAPER MOUTHPIECES					
GIMA code	Sizes - cm Ø int x ext x h	Suitable for	Compatibility		Minimum order
			GIMA codes	Other brands	
33410	2.8x3x6.5	adult	33431, 33433, 33471-2, 33510, 33512-14, 33518-22, 33528-32, 33534	Micromedical Vitalograph, Mir Ferraris (Koko legend)	box of 500
33897				Cardiette (Pneumos 300)	dispenser of 100
33411	2.0x2.2x6	paediatric		Micromedical Vitalograph, Mir	box of 250
33412	Plastic adapter		33411	Mir	-
33506	2.6x2.8x6.5	adult	33484-7	Sibel, BTL	box of 100
33890	2.5x2.7x7	adult	-	Fukuda Sangjo, Thor	dispenser of 100
33892	2.8x3.05x7	adult	-	Cosmed (Pony model)	dispenser of 100
33898	2.8x3.05x7	adult	-	Cosmed (Pony model)	box of 500
33895	3.03x3.3	adult	-	Sensormedics, Microcard	box of 500
33899	2.42x2.6	adult	-	Sorin Lifewatch, Fukuda Denshi	box of 500
GIMA PLASTIC MOUTHPIECES					
33432*	1.7x2.1x5 2.7x3x5	adult/ped	33431, 33433	Ferraris	-
33691		adult	33682-4, 6, 8	Thor	box of 10

*Made in USA

SPIROMETRY PAPERS

• 33493 MIR SPIROLAB PAPER - roll 112mm x 25m - box of 10

Suitable also for SIBEL DATOSPIR 120 B, C, D, DATOSPIR Touch and GIMASPIR 120 models B, C, D.

• 33891 COSMED PONY PAPER - roll 57 mm x 11 m - box of 25 rolls

• 33894 COSMED PONY FX PAPER - roll 112 mm x 13m - box of 25 rolls

Compatible also with Fukuda Spirosift.

• 35674 THERMAL PAPER - roll 57 mm x 25 m - box of 20

Suitable for SIBEL DATOSPIR 120 A and GIMASPIR 120-A.

• 33404 THERMAL PAPER - roll 111 mm x 8m - box of 25 rolls

Suitable for Microlab 3500

BACTERIAL FILTER

• 33525 BACTERIAL FILTER - box of 50 for MIR, Medisoft spirometers

• 33557 SPIROGUARD BACTERIAL FILTER - box of 50

Compatibility of Spiroguard

Erich Jaeger GmbH	Masterlab and Masterscreen
Ferraris Medical - P.K. Morgan	Autolink-Benchmark Transflo and USA Transfer System
Sensormedics	2200 - 2400 - 2450
Vitalograph	Wedge-Bellow-Alpha and Escort Spirometers

BACTERIAL FILTER - MOUTHPIECE FOR MOST SPIROMETER BRANDS



CE 0476

BACTERIAL FILTER - MOUTHPIECE

Disposable anatomical filter with integrated mouthpiece to prevent bacterial and viral contamination. Easy to use, direct connection to spirometer on one side, to patient on the other side. Suitable for most of spirometers and devices. Class II A. Individually packed in transparent bag.

MADE IN ITALY

BACTERIAL FILTERS - MOUTHPIECES - box of 100					
GIMA code	Size - cm Ø int x ext	Compatibility with			
		GIMA codes	Other brands		
33884	2.75x3	33341, 33433, 33471-2, 33509-14, 33518-22, 33528-32, 33534	Micromedical, Vitalograph, Mir, Ferraris (koko legend), Cardiette (Pneumos 300)		
33885	2.7x3.05	-	Cosmed		
33886	2.95x3.3	-	Customed, Jaeger		
33887	3.05x3.4	-	Sensormedics, Thor, Morgan, Chest, Microgard		
33888	3.4x3.6	-	Medical Graphics		
33889	2.5x2.7	-	Fukuda Sangjo		
33902	2.22x2.8	-	BTL, Sibel		

DISPOSABLE NOSE CLIPS FOR SPIROMETRY TESTS



CE 0476

33475

33476

33477

Rubber pads

Foam pads

Rubber large pads

• 33475 NOSE CLIPS - rubber - 70x40 mm - bag of 100

• 33476 NOSE CLIPS - foam - 70x40 mm - bag of 100

• 33477 NOSE CLIPS - rubber large - 70x45 mm - bag of 100

Gima rubber or foam-padded disposable nose clips help prevent patients from inhaling or exhaling through the nose during spirometry. The tension in the nose clip body keeps it firmly in place while the rubber or foam pads are designed for the best patient comfort. Made of thermoplastic, latex free. Rubber nose clip is available in 2 sizes: standard or large.

MADE IN ITALY



33476