

NON WOVEN SURGICAL SPECIALISTIC GOWNS - STERILE

SURGICAL SPECIALISTIC GOWNS

Specialistic, water repellent, breathable gowns for a barrier effect to fluids and bacteria

- cotton cuffs 8 cm
- ultrasonic welding
- neck velcro closure
- 4 waist straps for internal and external closure

Pre folded

For easy programmed dressing

Welded reinforcement (only 21015, 21016)

Water-proof, breathable reinforcement on front and on sleeves ultrasonic welded.

All gowns include two dry-paper towels. Conforms to UNI EN ISO 13795.




21020-21023

HIGH QUALITY
50 g/m²

MADE IN ITALY

GIMA code	SURGICAL GOWNS STERILE	Minimum order
SURGICAL GOWNS 35 g		
21010	Size M - length 115 cm	box 30
21011	Size L - length 130 cm	box 30
21012	Size XL - length 135 cm	box 30
SURGICAL GOWNS 50 g		
21020	Size M - length 120 cm	box 30
21021	Size L - length 130 cm	box 30
21022	Size XL - length 140 cm	box 30
21023	Size XXL - length 150 cm	box 30
REINFORCED GOWNS 35 g		
21015	Size L - length 130 cm	box 30
21016	Size XL - length 135 cm	box 30



SMMS DISPOSABLE SURGICAL GOWNS - NON STERILE OR STERILE

SURGICAL GOWNS

Disposable surgical gowns that guarantee comfort and protection.

- fluid repellent
- abrasion resistant
- ultrasonic sewing
- book folding
- knitted cuff/velcro at neck
- anti microorganism and anti body fluids



21046-21049

HIGH QUALITY
50 g/m²


- antistatic treatment
- great air permeability
- high strength and reliability
- reinforced pad 38 g (PP + PE) on the chest (50x70 cm 21052-21054 or 50x80 cm 21055) and sleeve (30 cm 21052-21054, or 40 cm 21055)
- set-in sleeve

Material: SMMS (spunbond + 2 meltblown + spunbond non wovens).

Weight per piece: 120 ± 2g
Sterile gowns are individually packed.

Sterilized by oxide Ethylene.
Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.

GIMA code	SURGICAL GOWNS				Minimum order
	Size cm	Length cm	Width cm	Length of sleeve cm	
SURGICAL GOWNS 40g - NON STERILE					
21040	M	120	140	58	box of 100*
21041	L	120	150	59	box of 100*
21042	XL	130	150	59	box of 100*
21043	XXL	140	160	61	box of 100*
Packing size 66x35x43 cm. *10 pouches of 10 pieces					
SURGICAL GOWNS 50g - STERILE					
21046	M	120	140	58	box of 50
21047	L	120	150	59	box of 50
21048	XL	130	150	59	box of 50
21049	XXL	140	160	61	box of 50
Packing size 66x37x38 cm.					
REINFORCED SURGICAL GOWNS 50g - STERILE					
21052	M	120	140	58	box of 50
21053	L	120	150	59	box of 50
21054	XL	130	150	59	box of 50
21055	XXL	140	160	61	box of 50
Packing size 66x37x40 cm.					



INSULATION GOWNS - WATERPROOF

WATERPROOF DISPOSABLE GOWNS PP 25 g

It can be used in clinics, patient rooms, laboratory, check point, community...

Not suitable for ICU and operating room.


- collar with two ties
- cuffs with elastic
- long sleeve
- waist with ties

Packing size: 52x33x30 cm

Gross weight: approx. 6.8 kg

Multilingual instructions:

GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.
Sold in a box of 100 (10 pouches of 10 pieces).



GIMA code	POLYTHENE GOWNS Size	Minimum order
21131	Medium 115x137 cm	box of 100
21132	Large 120x150 cm	box of 100

