

## DISPOSABLE INSULATION COVERALLS - TYPE 5B-6B AND TYPE 4B-5B-6B



### BASIC COVERALLS

Breathable, non-woven microporous coveralls with hood, elastic cuffs, elastic waist, zipper front.

#### Protection

Solid particulate resistant.

#### Applications

General outdoor duty, common insulation in out-patient service, inpatient ward, checkout room.

Each coverall is individually packed in a plastic bag. Multilingual instructions (GB, FR, IT, ES, PT, DE, GR, SE).

Meet PPE (EU) 2016/425 CAT III



**TYPE 5B-6B**  
EN 14126 - EN 13034



### TAPED SEAM COVERALLS

Breathable, non-woven microporous coveralls with hood, with taped seam, elastic cuffs, elastic waist, zipper front.

#### Protection

Dust proof, dry particulate/light spray/virus and blood permeate resistant.

#### Applications

Apply to General field, medical institutions outpatient, ward, laboratory and other common insulation use.

Each coverall is individually packed in a plastic bag. Multilingual instructions (GB, FR, IT, ES, PT, DE, GR, SE).

Meet PPE (EU) 2016/425 CAT III



**TYPE 4B-5B-6B**  
EN 14605 - EN 14126  
EN 13034



GIMA code	Size	Chest (cm)	Height (cm)
21070	Size S	84-92	164-170
21071	Size M	92-100	170-176
21072	Size L	100-108	176-182
21073	Size XL	108-116	182-188
21074	Size XXL	116-124	188-194

GIMA code	Size	Chest (cm)	Height (cm)
21080	Size S	84-92	164-170
21081	Size M	92-100	170-176
21082	Size L	100-108	176-182
21083	Size XL	108-116	182-188
21084	Size XXL	116-124	188-194

TYPE-4B, spray-tight seams, offer protection against saturation of liquid chemicals  
TYPE-5B, against solid dust, provide protection to the full body against airborne solid particulates  
TYPE-6B, against minor splashing, offer limited protection against a light spray of liquid

## MICROPOROUS BOOT COVERS - TYPE 5B-6B AND TYPE 4B-5B-6B



**TYPE 5B-6B**  
EN 14126

**TYPE 4B-5B-6B**  
EN 14126

**21090**

Size 45x36x26 cm

PVC antislip sole

**21091**

Size 45x40x26 cm

**21092**

Size 45x36x26 cm

- 21090 BOOT COVERS - TYPE 5B-6B - box of 500 (25 bags of 20 pcs)
  - 21091 BOOT COVERS - TYPE 5B-6B - with antislip sole - box of 40
  - 21092 TAPED SEAM BOOT COVERS - TYPE 4B-5B-6B - box of 20
- Microporous disposable boot covers for general purposes as: cleanrooms, hospitals, labs, food industry, medical and pharmaceutical. Optimized body fit and lightweight (60 g) to improve comfort and safety. The PP material provides high level of comfort and protection against solid particles such as dust.

#### Features:

- PP + PE, microporous
- breathable and dust proof
- liquid penetration resistance
- barrier to fine particulates
- air and water vapor permeable to reduce risk of heat stress and fiber contamination
- with ties, elastic open.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE.



## SHOE COVERS

- 21094 SHOE COVER - box of 2000 (20 bags of 100)
  - 21095 ANTI-SKID SHOE COVER - box of 1000 (10 bags of 100)
- made of lightweight PP
  - with elastic band
  - anti-skid model with elastic and printing at the bottom (21095)
  - size: 17x40 cm (21094), 15x40 cm (21095) - light, only 35 g



**21094**

**ANTI-SKID**

**21095**

## CAPS WITH ELASTIC OR TIES

- 21097 CAP WITH ELASTIC - box of 1000 (10 bags of 100)
  - 21098 CAP WITH TIES - box of 1000 (10 bags of 100)
- operating cap made of SMS
  - size: 62x14,5 cm
  - light, 35 g (21097) or 40 g (21098)
  - breathable and waterproof



**21097**

Cap with elastic

**21098**

Cap with ties