

GIMA ECG THERMAL PAPERS FOR DIFFERENT ECG BRANDS

GIMA code	ECG MANUFACTURER AND MODEL	Type	Width mm	Length m	Sheet	Marker	Grid'	Min. order	GIMA code	ECG MANUFACTURER AND MODEL	Type	Width mm	Length m	Sheet	Marker	Grid'	Min. order
33012	MAC 400, MAC 600	pack	80	90	280	yes	or/red	25	32999	Page Writer TC/50 - TC/70	pack	210	300	200	hole	or/red	1
33013	MAC 800	pack	110	140	200	yes	or/red	25	32997	Touch, TRIM I / II / III, Page Writer TC/20, TC/30	pack	210	300	100	hole	or/red	1
33058	MEDIANA YM-303	roll	80	30	-	-	orange	10	32953	PROGETTI EPG 612	roll	112	23	-	no	or/red	10
29528	MEDICAL ECONET SMART - 3		152	25	-	-	green	10	32960	EPG 12 view	roll	210	30	-	-	orange	5
32950	Cardio M - Cardio M Plus	roll	215	25	-	-	orange	5	32964	EPG Project	roll	80	25	-	-	orange	10
29522	SMART - 1	roll	215	23	-	-	green	5	32996	EPG Basik-pocket	roll	58	25	-	no	or/red	20
33058	MEDIGATE MeCA 303i	roll	80	30	-	-	orange	10	33037	SCHILLER AT 1, AT 4	pack	90	90	390	yes	green	25
33014	MINDRAY Beneheart R3	roll	80	20	-	no	or/red	10	32982	AT 2, CS 200	pack	210	280	215	yes	or/red	1
33016	Beneheart R3	pack	80	70	200	yes	or/red	25	33038	AT 10	pack	210	140	250	yes	green	10
32966	BeneHeart R12	pack	210	295	150	yes	or/red	1	32983	AT 101	pack	80	70	300	yes	or/red	20
33018	Mindray D60	roll	112	20	-	no	orange	10	32984	AT 102	pack	210	280	200	yes	or/red	1
33017	MORTARA Eli 50	pack	50	70	200	yes	or/red	25	33039	AT110, AT 10 Plus, AT 2 Plus	pack	210	140	150	yes	green	5
33010	Eli 100	roll	108	23	-	no	or/red	10	32987	Cardiopad, MS 2010	pack	114	90	200	yes	green	20
33019	Eli 210, 250, 280, 350, 380	pack	210	300	250	yes	or/red	1	33040	MS-2015	pack	210	140	180	yes	or/red	5
33020	Eli 150	pack	108	140	200	yes	or/red	25	33041	Cardioart FT/1	pack	114	150	64	yes	orange	20
33021	Eli 230	roll	210	20	-	no	or/red	5	33042	Cardioart AT/1	pack	80	70	300	yes	orange	20
32988	NIHON KOHDEN ECG-6651 - Def. Series TEC 7100 / 7200/ 7300	roll	50	30	-	-	orange	20	33043	Cardioart AT/102-G2	pack	210	280	100	yes	orange	1
32994	Cardiofax 5403 - ECG-8240		145	30	-	-	orange	10	32990	SIEMENS Cardiostat 1	pack	50	75	400	yes	or/red	20
33025	Cardiofax 8110, 9620, 1150K	roll	63	30	-	no	or/red	20	32991	Cardiostat 11	pack	63	75	400	yes	or/red	25
33026	Cardio-fax GEM/1250/1950, ECG-9010/9020/9022	pack	210	140	215	yes	white	10	32988	TRISMED Cardipia 200 / 201N	roll	50	30	-	-	orange	20
32962	Cardio-fax GEM/1250/1950, ECG-9010/9020/9022	pack	210	140	215	yes	or/red	10	32989	Cardipia 200, 203N	roll	63	30	-	no	or/red	20
33027	GEM 9010, 9020, 9022, 1250, 1950	pack	110	140	143	yes	white	20	33045	Cardipia 400, 406RS	roll	110	30	-	no	or/red	10
33028	1500, 1550, 9110, 9130, 9132, 9320, 2250	pack	210	295	170	yes	white	1	33047	Cardipia 800, 812RS	roll	215	30	-	no	or/red	5
32998	PHILIPS Page Writer TC/10	roll	110	25	-	no	or/red	10	33647	ZOLL R-Series, X-Series	roll	80	25	-	-	white	10
									32964	X-Series	roll	80	25	-	-	orange	10
									32962	ZONCARE iMac 120 / iMac 1500 / iMac 1800	pack	210	140	215	-	orange	10
									32964	ZQ-1203G / iMAC 300 / iMAC 30	roll	80	25	-	-	orange	10

FIAB ECG ELECTRODES WITH OFF-CENTERED SNAP CONNECTION

DISPOSABLE ELECTRODES WITH OFF-CENTERED SNAP CONNECTOR						
GIMA code	FIAB code	Material	Size	Gel	Patient	Minimum order
33308	F2017P	PE-Foam	42 mm	Liquid	Pediatric	box of 1,000 (20 bags of 50)
33309	F2017	PE-Foam	50 mm	Liquid	Adult	box of 1,000 (20 bags of 50)
33310	F2017H	PE-Foam	60 mm	Liquid	Adult large	box of 750 (25 bags of 30)


All electrodes have silver chloride sensor



3M™ RED DOT™ ELECTRODES

GIMA code	3M code	3M™ RED DOT™ ELECTRODES	Minimum order
Gel electrodes - short term diagnostics			
33378	2330	3M™ Red Dot™ - 3.2 x 2.2 cm	bag of 100
Foam monitoring electrodes			
33293	2570	3M™ Red Dot™ - 4 x 3.5 cm	bag of 50
Non woven monitoring electrodes			
33299	2248-50	3M™ Red Dot™ - Pediatric - Ø 4.5 cm	box of 50
33297	2239	3M™ Red Dot™ - Ø 6 cm	box of 1,000*

*20 bags of 50 pieces



ECG GEL AND SPRAY

- 33266 ECG GEL - 250 ml - box of 40 tubes
- 33267 ECG GEL - 5 l - box of 2 bags**
- 33268 ECG SPRAY GEL - 250 ml - box of 20 tubes

Suitable for non pre-gelled electrodes, ECG latches and pendant switches in any ECG and exertional applications. Allows ECG waves clear transfer to equipment.

It is hypoallergenic, does not irritate skin, does not contain oil and oily substances, does not contain formaldehyde and salt, has no toxic effects and is odourless.

- 33283 PUMP FOR BAG 5 l

**Includes 1 empty bottle 250 ml and 1 nozzle for easy refill

