




ECG CABLES FOR BM GIMA AND BIONET MONITORS

GIMA code	ECG CABLES FOR BM3 (33718, 33719, 33720, 33780, 33781)	ACCESSORIES FOR MONITOR BM3, BM5, BM7	
		NEW (after 6/2006)	OLD (before 6/2006)
33711	3-lead patient cable kit (cable+snap connectors 33748)	 33713 Vet cable kit 33760 Vet clip connect 33762 Cable 33748 Snap connect	 33752 Vet clip connect 33750 Snap connect
33713	3-lead veterinary cable kit (cable+3 clips connectors)		
33748	3-lead snap connectors for 33711 - spare - new		
33749	5-lead snap connectors for 33712 - spare - new		
33750	3-lead snap type connectors for 33726 - spare - old		
33751	5-lead snap type connectors for 33728 - spare - old	 33761 Vet clip connect 33763 Cable 33749 Snap connect	 33753 Vet clip connect 33751 Snap connector
33752	Vet 3-lead clip type connectors for 33726-7 - spare - old		
33753	Vet 5-lead clip type connectors for 33728 - spare - old		
33760	Vet 3-lead clip type connectors for 33713 - spare - new		
33761	Vet 5-lead clip type connectors for 33712 - spare - new		
33762	3-lead patient cable		
33763	5-lead patient cable		
ECG CABLES FOR BM5-BM7 (33776, 33777, 33778, 33779)			
33764	3-lead patient cable kit (cable + snap connector)		
33765	10-lead patient cable kit (cable + snap connector) only 33778		
33766	3-lead veterinary cable kit (cable + clip connector)		



BM1, BM3, BM5, BM7 SpO₂ SENSORS

GIMA code	SpO ₂ SENSORS AND EXTENSION CABLE FOR ALL BM MONITORS	SpO ₂ REUSABLE PROBES		SpO ₂ EXTENSION CABLE
		7 pin	33730 33731	33716
33730	Adult reusable finger sensor* (> 45 kg) - 7 pin			
33731	Pediatric reusable finger sensor* (15-45 kg) -			
33732	Vet reusable sensor* - clip *Needs extension cable 33716			
33716	Extension cable**, 2 m (7 pin big SpO ₂ connector) since 2007 **Necessary for all probes			

NELLCOR COMPATIBLE DISPOSABLE SpO₂ SENSORS

						All sensors are supplied in sterile pack with 70 cm cable.
33757		33754	33755	33756		
GIMA code	NELLCOR COMPATIBLE DISPOSABLE SpO ₂ SENSORS (not compatible with Nellcor Oximax)					TECHNICAL SPECIFICATION SpO₂ SENSORS Range of SpO ₂ measurement: 70%~100% Accuracy: 70%~100% ±2%; below 70%: unspecified Range of pulse measurement: 30~250 bpm Accuracy: ±2 bpm of 2% (select larger) Optical Sensor: red light (wavelength is 650~670 nm, 6.65 mW) infrared (wavelength is 880~910 nm, 6.75 mW)
33755	Disposable neonate/adult sensor - under 3 kg (neonate foot) - over 40 kg (adult finger)					
33756	Disposable infant sensor - 3-20 kg					
33757	Disposable pediatric sensor - 10-50 kg					
33754	Disposable adult sensor - over 30 kg					

BM1, BM3, BM5, BM7 MONITOR ACCESSORIES

GIMA code		ACCESSORIES FOR BM1, BM3, BM5, BM7					
33738	NIBP	Extension cable, 3 m (necessary for all cuffs)	 Kit of 3 neonatal cuffs* 33737	 33735	 33736		
33735		Adult cuff 27.5-36.5 cm need ext. cable 33738					
33736		Child cuff - need ext. cable 33738					
32848		Obese cuff - need ext. cable 33738					
33737		Kit of 3 neonatal cuffs* - needs cable 33738					
32810	Temp	Mini cuff size 16x5 cm - needs 33738+33739	 32810	 32811			
32811		Mini cuff size 23x5 cm - needs 33738+33739					
33739		Metal connector (necessary for mini cuff)					
33771		Temperature sensor 3.0 m, skin type***					
33772		Temperature sensor 3.0 m, rectal type					
33743	Paper	Thermal paper 57.5 mm, 38 Ø	 33738	 33771	 33772		
33715	Battery	Ni-MH 12V battery - Since January 2006** - BM3, BM3 Plus					
33707		Li-ion battery 10.8 V for BM3 Pro, BM5 Pro					
33717		Li-ion battery - Since Sept. 2007** - new BM3, Plus, BM5					
33710		Li-ion battery 11.1 V for BM1					
33709		Li-ion battery 11.1 V for BM7					
33747	Adapter	AC adapter for BM3, BM5, BM7	 33715	 33717			
*Cuff range 3.3-5.6 cm, 4.2-7.1 cm, 5-10.5 cm							
**Date is approximative. Please select battery according to picture and specifications							
***Compatible with Drager MM15114, MM15822, GE 21733601, Philips HP21059B, Nikon Kohden B3424, Datascan 0226-02-0002							

*Cuff range 3.3-5.6 cm, 4.2-7.1 cm, 5-10.5 cm

**Date is approximative. Please select battery according to picture and specifications

***Compatible with Drager MM15114, MM15822, GE 21733601, Philips HP21059B, Nihon Kohden P242A, Datascope 0206-02-0003