

HEINE CLASSIC+® F.O. XHL OR LED^{HQ} LARYNGOSCOPES

HEINE CLASSIC+® F.O. XHL OR LED^{HQ} LARYNGOSCOPES

The "Classic+" range of F.O. blades with HEINE HiLite® Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative HEINE LED^{HQ} light in the handle or the traditional XHL Xenon Halogen Bulb, the exclusive HEINE HiLite® Fibre Optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs. **5 years warranty or 4,000 standard autoclave cycles (134°C/5 min.)** more than 1,000 lux on HEINE HiLite® fibre bundle blades. **Long lasting bulbs** HEINE XHL® Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40% more light! XHL 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. HEINE LED^{HQ} handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding bright illumination. **Long lasting 3.5 V rechargeable batteries** HEINE 3.5 V systems (LED^{HQ} or XHL) use a long-lasting 3.5 V lithium rechargeable battery. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).

• **34374 NT 4 CHARGER UNIT FOR 2 HANDLES - X-002.99.494**




34374 Two independent charging compartments:
 - fast charge capability:
 - 2 hours for Li-ion handle
 - automatic charging
 - charger indicator

• **34449 NT300 CHARGER UNIT FOR 2 HANDLES - X-002.99.495**

Size: 14x6x6.6 cm



GIMA code	HEINE code	HEINE FIBRE OPTIC BLADES	Overall length	Distal length
34411	F-000.22.119	MILLER 00	76 mm	9.5 mm
34412	F-000.22.120	MILLER 0	80 mm	10.5 mm
34413	F-000.22.121	MILLER 1	100 mm	11.5 mm
34414	F-000.22.122	MILLER 2	152 mm	13 mm
34415	F-000.22.110	PAED 0	82 mm	12 mm
34416	F-000.22.111	PAED 1	95 mm	13 mm
34417	F-000.22.100	MAC 0	82 mm	9 mm
34418	F-000.22.101	MAC 1	93 mm	11 mm
34419	F-000.22.102	MAC 2	115 mm	13 mm
34420	F-000.22.103	MAC 3	135 mm	15 mm
34421	F-000.22.104	MAC 4	155 mm	15 mm
34422	F-000.22.105	MAC 5	176 mm	15 mm
34423	F-000.22.313	MAC 3 - Flex Tip+	115 mm	13.2 mm
34368	F-000.22.314	MAC 4 - Flex Tip+	137 mm	13.2 mm
34361	F-000.22.763	XP MAC 3 - disposable - box of 25	134 mm	14 mm
34362	F-000.22.764	XP MAC 4 - disposable - box of 25	154 mm	14 mm



HEINE STANDARD LARYNGOSCOPE HANDLES AND ACCESSORIES

34430	F-001.22.800	Slim F.O. handle 2.5 V - pediatric
34400	F-001.22.812	Short F.O. handle 2.5 V - pediatric and obese
34401	F-001.22.860	Standard F.O. handle 2.5 V
34405	F-008.22.860	Standard F.O. LED handle 2.5 V
34370	F-008.22.820	EasyClean® LED handle - 100% waterproof
31769	X-001.88.035	Spare bulb 2.5 V for 34400, 34401, 34430
31768	X-002.88.044	Spare bulb 3.5 V for 34369

HEINE RECHARGEABLE LARYNGOSCOPE HANDLES AND BATTERIES

34371	F-008.22.804	Pediatric LED handle 2.5 V Li-ion with recharg.* battery
34372	F-007.22.894	Standard F.O. handle 3.5 V Li-ion with recharg.* battery
34373	F-008.22.894	Standard LED handle 3.5 V Li-ion with recharg.* battery
34402	F-008.22.891	Standard LED handle 3.5 V Li-ion with recharg.** battery
34378	X-007.99.104	Re-chargeable 2.5 V battery pediatric for 34371 - spare
34448	X-007.99.383	Re-chargeable 3.5 V battery for 34372-3, 34402, 34367 - spare

*Need NT4 charger unit code 34374 **Need NT300 charger unit code 34449

HEINE CLASSIC+® FIBRE OPTIC LARYNGOSCOPE SETS**

with 2.5V handle		
34403	F-119.10.860	3 blades set with Paed 1, Mac 2, Mac 3, XHL handle
34404	F-120.10.860	4 blades set with Paed 1, Mac 2, Mac 3, Mac 4, XHL handle
34446	F-120.18.860	4 blades set with Paed 1, Mac 2, Mac 3, Mac 4, LED ^{HQ} handle
with 3.5V rechargeable handle		
34367	F-119.94.865	3 blades set with Paed 1, Mac 2, Mac 3, LED HQ rechargeable handle

**Set includes spare bulb and F.O. handle 2.5 V (34401) or F.O. LED handle 2.5 V

STANDARD HANDLE 3.5V LED Li-ion®
 Rechargeable 3.5V F.O. LED handle with charge status display in bottom cap. State of charge is indicated for 10 seconds after switching on or after removing from the charger. Two hours charging time with NT4 charger; up to 13 hours ON time. 12,000 LUX at 40 mm.

EASYSOFT® LED HANDLE
 LED Handle for low-temperature gas plasma sterilisation (STERRAD®/STERIS®). 100% waterproof (IP68) for easy wiping and immersion disinfection. No disassembly required. Easy-to-wipe "wave design" Autoclavable LED technology: no light removal during disinfection, just batteries. 10 hours operating time with 2xC batteries (not included).

LED HANDLES:
 - 5,000 Kelvin
 - 50,000 hours bulb life

34373

34370

