

# GIMA GREEN F.O. XENON/HALOGEN AND LED LARYNGOSCOPES

## GIMA GREEN LARYNGOSCOPES

- maximum light transmission thanks to durable and high quality fibres
- specially designed blades for easy intubation
- matt finished to reduce glare
- autoclavable up to 134°C for 5 minutes. Blades autoclavable up to 2,000 times
- compatible to plasma sterilization
- all blades are compatible to BS EN ISO 7376 Green standard handles
- all blades are made of Austenitic Stainless steel thereby making them rust proof
- long life blades. GIMA offers 3 years Warranty
- a large selection of power sources to use: 'AA' and 'C' type batteries
  - or 2.5 V and 3.5 V rechargeable batteries
  - 2.5 V handles can be used with standard or rechargeable batteries
- handles are available with 2.5 V and 3.5 V Xenon/Halogen and LED bulbs
- all blades are non magnetic

### Mc INTOSH BLADES



### MILLER BLADES

### MAXLITE

58056

58052

### GIMA GREEN



#### Gima Green and Maxlite: Miller and Mc-Intosh blades

- Gima Green 4.0 mm Ø F.O. bundle with more than 5,500 individual fibre strands
- Miller Maxlite 3.0 mm Ø bundle with more than 5,000 individual fibre strands
- McIntosh Maxlite 4.5 mm Ø bundle with more than 6,500 individual fibre strands

**High quality fibres ensuring excellent light transmission Gima Green and Miller Maxlite:** 7,000 lux with 2.5 V Xenon/Halogen lamp, 12,500 lux with 3.5 V Xenon/Halogen lamp, 18,000 lux with 2.5 LED lamp, 30,000 lux with 3.5 V LED lamp

McIntosh Maxlite: 8,000/14,000/20,000/40,000 lux

GIMA code	"GIMA MAXLITE" FIBRE OPTIC BLADES	A	B
58051	MILLER N° 0 blade - Newborn	80	55
58052	MILLER N° 1 blade - Infant	100	75
58056	Mc INTOSH N° 0 blade - Newborn	82	57
58057	Mc INTOSH N° 1 blade - Infant	93	68
58058	Mc INTOSH N° 2 blade - Child	115	90
58059	Mc INTOSH N° 3 blade - Adult	135	110
58060	Mc INTOSH N° 4 blade - Adult large	155	130
58061	Mc INTOSH N° 5 blade - Adult x-large	176	151
GIMA code	"GIMA GREEN" FIBRE OPTIC BLADES*	A	B
34460	MILLER N° 00 OXY blade - Premature	65	40
34461	MILLER N° 0 blade - Newborn	80	55
34462	MILLER N° 1 blade - Infant	100	75
34463	MILLER N° 2 blade - Child	152	127
34464	MILLER N° 3 blade - Adult	192	167
34469	Mc INTOSH N° 0 blade - Newborn	80	55
34470	Mc INTOSH N° 1 blade - Infant	95	70
34471	Mc INTOSH N° 2 blade - Child	115	90
34472	Mc INTOSH N° 3 blade - Adult	134	109
34473	Mc INTOSH N° 4 blade - Adult large	155	130
34474	Mc INTOSH N° 5 blade - Adult x-large	175	150

\*equipped with detachable fibre optic light guide

#### F.O. LIGHT GUIDE FOR GIMA GREEN BLADES

- 34483 F.O. Light guide for "GIMA Green" Mc Intosh N° 3
- 34484 F.O. Light guide for "GIMA Green" Mc Intosh N° 4
- 34486 F.O. Light guide for "GIMA Green" other - specify code

• **31542 GIMA GREEN CHARGING STATION - 100/240 V - 50/60 Hz**  
It can charge either 2.5 V / 3.5 V NiCa, NimH or Li-ion rechargeable batteries independently of their status of charge. Unit is designed to prevent over charging of batteries. Simultaneously charging of pediatric and adult handle.

Multilingual manual:  
GB, FR, IT, ES

31542




Charging station includes pediatric adaptor,  
wall support, UK, US and EU plug

377

### HALOGEN/XENON BATTERY HANDLE

### LED BATTERY HANDLE



#### Xenon/Halogen and LED handles 2.5 V and 3.5 V

- closed battery compartment prevent liquid penetrating inside the battery
- handles are supplied with lamp. Battery not included
- solid metal handle of chromium plated metal

#### GIMA GREEN RECHARGEABLE BATTERY HANDLES

GIMA code	F.O. XENON/HALOGEN BATTERY HANDLES - 2.5 V
34480	Pediatric rechargeable metal handle (2xAA) - 2.5 V
34481	Adult rechargeable metal handle (2xC) - 2.5 V
34466	Pediatric rechargeable Li-ion battery 2.5 V 1,000 mAh
34467	Adult rechargeable Li-ion battery 2.5 V 2,200 mAh
34485	Halogen bulb 2.5 V (minimum order 6 pcs)
"GIMA GREEN" RECHARGEABLE HANDLES AND BATTERIES 3.5 V - OLD up to lot 8 (ped handle) up to lot 7 (adult handle)	
34489	Pediatric rechargeable NiCa battery 3.5 V/700 mAh for 34487
34491	Halogen Bulb 3.5 V

#### F.O. XENON/HALOGEN BATTERY HANDLE 3.5 V

34456	Pediatric rechargeable metal handle (3.5 V)
34457	Adult rechargeable metal handle (3.5 V)
31540	Pediatric rechargeable Li-ion battery** 3.5 V for 34456, 34487***
31541	Adult rechargeable Li-ion battery** 3.5 V for 34457, 34482, 34488****
**for handles sold after may 2013 ***since lot 10 ****since lot 12	
34491	Halogen Bulb 3.5 V

#### F.O. LED BATTERY HANDLES 2.5 V AND 3.5 V\*

34476	Pediatric LED rechargeable handle 2.5 V
34477	Adult LED rechargeable handle 2.5 V
34482	Adult rechargeable metal handle - LED light 3.5 V
*LED bulb is assembled with handle in a fixed system and is not replaceable	

#### COMPATIBILITY HANDLE - BATTERY - BULB

GIMA code	Handle	Battery	V	Standard	Rechargeable	Bulb
	Xenon	LED		2xAA	2xC	Li-ion*
34480	Ped.	-	2.5	yes	no	34466
34481	Adult	-	2.5	no	yes	34467
34456	Ped.	-	3.5	no	no	31540
34457	Adult	-	3.5	no	no	31541
34476	-	Ped.	2.5	yes	no	34466
34477	-	Adult	2.5	no	yes	34467
34482	-	Adult	3.5	no	no	31541
34485						
34486						
34487						
34488						
34489						
34490						
34491						

\*Li-ion batteries have twice the capacity of an equivalent size NimH batteries.