

GIMA RESUSCITATION KITS

- 34181 SPEED 1 RESUSCITATION KIT - UNI cylinder
- 34182 SPEED 2 RESUSCITATION KIT - NF cylinder
- 34183 SPEED 3 RESUSCITATION KIT - without cylinder
- 34180 SPEED EMPTY CASE - with sponge
- 34222 FACEMASK N° 2 - child small

First aid kits for reanimation, very useful for medical and dental rooms, sport clubs, factories, ambulances... Case made of anti-shock plastic. Size: 43 x 37 x 11 cm.

| STANDARD ACCESSORIES | SPEED 1 | SPEED 2 | SPEED 3 |
|---|---------|---------|---------|
| Silicone resuscitation bag - autoclavable | • | • | • |
| Silicone Face mask n° 4 | • | • | • |
| Oxygen cylinder 0.5 l (empty). UNI standard | • | • | • |
| Oxygen cylinder 0.5 l (empty). NF standard | • | • | • |
| High pressure manometer | • | • | • |
| Guedel airways (small + medium + large) | • | • | • |
| Helicoidal mouthwedge + Tongue holding forceps | • | • | • |
| Atoxic tubing for connection to oxygen cylinder | • | • | • |
| Plastic case | • | • | • |



34181

GIMA 3 FIRST AID CASE

- 34151 GIMA 3 - FIRST AID CASE + OXYGEN CYLINDER

Advised for up to 50 people

Size: 32x44.5x15 cm

Meets Italian law D.M. 388 (15.7.03) art. 1 + emergency kit.

Emergency kit includes:

- autoclavable resuscitation with security valve and silicone mask n. 4
- oxygen bottle 0.5 l (empty) with reductor and oxygen distributor - UNI standard
- 3 guedel airways, atoxic tube for connection to oxygen, mouth opener and wall support.

See complete case contents page 384

- 34148 GIMA 3 - EMPTY PLASTIC CASE



34151

34148



MEDICAL SPORT BAG



- 34135 MEDICAL SPORT BAG
- Special kit designed for first aid in sporting clubs (tennis clubs, swimming pools, gymnasium,...).

Soft nylon bag has 2 side pockets and zip closure.

Size: 470x240x190 mm
Weight: 6.5 kg

STANDARD ACCESSORIES

- | | |
|---|---|
| 1 isothermic blanket Gold/Silver 160x210 cm | 1 pack of 20 adhesive plasters 2x7 cm |
| 1 Disinfectant liquid 250 ml | 8 adhesive plasters 10x6 cm |
| 2 adhesive plasters 5 m x 2,5 cm | 1 sterile TNT bandage 60x40 cm |
| 2 pairs of sterile gloves | 1 sterile TNT bandage 60x80 cm |
| 2 VENTO mouth to mouth resuscitators | 6 sterile TNT pads 10x10 cm |
| 6 sterile gauze pads 18x40 cm | 5 instant ice packs |
| 2 instant ice sprays 150 ml | 1 hydrogen peroxide solution |
| 3 triangular TNT bandages 96x96x136 cm | 1 elastic net bandage |
| 3 Elastic bandages 6 cm | 1 haemostatic lacet |
| 3 Elastic bandages 8 cm | 1 sterile tweezers |
| 1 cotton wool pack | 5 ammonia towelettes |
| 1 pair of LISTER scissors 14,5 cm | 1 sterile saline solution bottle 250 ml |
| 2 compressive pad bandages 80x100 mm | 1 gel hands disinfectant 100 ml |
| 2 compressive pad bandages 100x120 mm | guidance leaflet 5 languages |

FIRST AID FAMILY CASE



Size: 235x190x43 mm

34136

Case in polypropylene

Guidance in:
GB, FR, IT, ES,
DE, TR

MADE IN ITALY

- 34136 FAMILY FIRST AID CASE
 - 34450 EMPTY CASE
- Small first aid case suitable for use anywhere and anytime: (at home or office, in the car or boat...).

STANDARD ACCESSORIES

- | | |
|--|-----------------------------------|
| 2 sterile gauze pads 18x40 cm | 3 antiseptic wipes |
| 1 pair of scissors. 1 plastic tweezers | 1 adhesive dressing 10 x 6 cm |
| 1 sterile bandage 60x40 cm | 1 gauze bandage 5 cm |
| 1 cotton wool pack 5 g | 1 adhesive plaster roll 2 cm x 1m |
| 1 pack of 10 adhesive plasters | 1 pair of protection gloves |
| 1 tissue pack 10 pcs | Guidance Leaflet 6 languages |

FIRST AID CASE AND BAG



34184

Advised for offices or workshops of 2 people
Meets Italian law 388-Art. 2

- 34184 SMALL KIT - plastic case
- Size: 290x215x90 mm - 1.4 kg
See content of kit at page 334
- 34188 REFILL PACK for 34184-34185-34194
 - 34186 EMPTY PLASTIC CASE
- 290x215x90 mm
- 34185 MIZAR BAG - nylon
- Size: 250x180x90 mm
Same contents as 34184 see page 384

34185

