

DOCTOR'S KIT

STANDARD ACCESSORIES

	34063	34066	34067
Utility bag - blue	•		
Multiuse bag - blue/grey		•	•
Aeon A200 non contact infrared thermometer (GB, FR, IT, ES)	•	•	•
Oxy-3 finger oximeter (GB, FR, IT, ES, PT)	•		
Oxy-4 finger oximeter - blue		•	•
Kobe palm aneroid sphygmomanometer	•	•	•
Obese cuff + 1 tube bladder		•	•
Trad stethoscope - Y black	•		
Classic dual head stetho - Y black		•	•
Sigma-C LED otoscope - black - pouch		•	•
GIMA glucose monitor kit GB, FR, ES, PT		•	
GIMA glucose monitor kit IT, DE, GR, Arabic			•
GIMA LED penlight	•	•	•
Dejerine GIMA 2000 - neurological hammer	•	•	•
Fast tourniquet	•	•	•



• 34063 DOCTOR'S KIT 1



• 34066 DOCTOR'S KIT 4A
• 34067 DOCTOR'S KIT 4B

S/S UTILITY AND BANDAGES SCISSORS

UTILITY AND BANDAGES SCISSORS

All purpose, first aid scissors, easily cut bandages, leather and plaster casts. Autoclavable at 136°C. Plastic coloured finger-contour handles for maximum comfort. Instructions: GB, FR, IT, ES, PT, DE.

• 20590 STEEL CHAIN for scissors and forceps



UTILITY SCISSORS - box of 10				
12 cm - 4.7"	16.5 cm - 6.5"	19 cm - 7.5"	Colour	
20560	20570	20580	Black	
-	20571	20581	Blue	
-	20572	20582	Red	
-	20573	20583	Purple	
-	20574	20584	Pink	

SOFT BAG KIT DIN



• 34137 SOFT BAG KIT DIN 13164-2014
Contents conform to DIN norms 13164.
Size: 260x172x82 mm. Weight: 0.7 kg.

STANDARD ACCESSORIES

- 1 adhesive plaster roll 5 m x 2.5 cm
- 4 disposable vinyl gloves
- 2 gauze bandages 6 cm x 4 m
- 3 gauze bandages 8 cm x 4 m
- 1 pair of lister scissor 14.5 cm
- 1 isothermic blanket gold/silver 160x210
- 1 sterile compressive pad 60x80 mm
- 2 sterile compressive pads 80x100 mm
- 1 sterile compressive pads 100x120 mm
- 1 sterile TNT burn sheets 60x40 cm
- 1 sterile bandage TNT 60x80 cm
- 3 sterile gauzes 10x10 cm
- 2 alcohol 70% wipes
- 1 pack 14 assorted patches
- 2 triangular bandages TNT 96x96x136 cm


Instructions (GB, FR, IT, ES, DE, TR)



EMERGENCY SCISSORS

• 34127 EMERGENCY SCISSORS 23.5x10.8 cm

Autoclavable Stainless steel AISI 420.



MADE IN ITALY

Supplied with transport bag



34127

Weight 540g

SPLINTER REMOVAL KIT



34190

• 34190 SPLINTER REMOVAL KIT

Includes everything to remove splinters and for first aid medications. In a polypropylene box size 235x190x43 mm:

- 1 Eye wash solution 100 ml with an eye glass
- 1 magnifying glass, 2 sterile gauze pads 18x40 cm
- 1 Professional removal pencil, 1 mirror.
- 1 stainless steel tweezers 11.5 cm

• 34192 PROFESSIONAL REMOVAL PENCIL

Small metallic splinters eye tool remove. On one side a magnetic point to remove any metallic splinter avoiding any risk for the eye. On the other side an antiabrasive nylon loop for non metallic splinters.

Small metallic splinters, eye tool remove

Anti-abrasive nylon loop

34192

"NEODIMIO" magnetic point

The body is made of anodised aluminium and the neoprene seals guarantee a better pencil hold. Supplied in a box and CE marked.

EYE WASH SYSTEMS

• 34889 DROP - eye wash station

Thermo-formed wall kit with all necessary items for efficient first emergency eye wash:

- 2 sterile saline solutions of 500 ml
- 2 integrated eye cups and one mirror
- screw and plug to be wall-mounted.

• 34887 EYE WASH SOLUTION 250 ml

- sterile - spare - box of 18

• 34888 EYE WASH SOLUTION 500 ml -

sterile - spare
Disposable, 0.9% sodium chloride sterile solution for emergency eye wash. Simple to use. The bottle is ready to use and fitted with a cap.



34887

MADE IN ITALY



34889