

SILICONE NURSE WATCHES



NURSE WATCHES

Silicone nurse watches, provided with a safety pin.

Suitable to measure patient pulse, beat and sterilization procedure timing.

Conforming to procedures not allowing nurses to wear wrist watches.

Watch case 35 mm, watch length 85 mm.

Shelf life of battery: 12-15 months from first use.

Manual: GB, FR, IT, ES, PT, DE, GR, PL, Arabic.



20403

20411

GIMA code	SILICONE WATCHES Round - single colour	GIMA code	SILICONE WATCHES Round - single colour with aluminium box
20400	White	20410	White
20401	Blue	20411	Blue
20402	Red		Square - single colour with aluminium box
20403	Purple		
20404	Pink	20455	Blue
20405	Black		
	Round - multicolour		Square - multicolour
20420	Red & white	20460	Violet and sky blue
20421	Patchwork	20461	Hearts
20422	Zebra	20462	Flowers
20423	Bubbles	20463	Leopard
20424	Leaves	20464	Roses
20425	Caterpillar	20465	Foot print
	Silicone straps - spare		Digital**
20418	Kit of 6 silicone straps*	20480	Blue

*Spare silicone straps for round watches

**Manual: GB, FR, IT, ES

20480



NURSE STAINLESS STEEL SCISSORS AND FORCEPS



20513

20500

20512

Nurse bandage scissors - 14 cm



20542

20551

20552

Straight artery forceps - 16 cm



20521

20529

20520

Straight scissors blunt/sharp - 14 cm

Model	BANDAGE* 14 cm	SCISSORS 18 cm	STRAIGHT* SCISSORS - 14 cm	ARTERY* FORCEPS - 16 cm
Standard	34129**	34130	26725***	26717
Standard with belt clip	20500	-	20520	-
Green ring	20501	-	20521	20541
Lilac ring	20502	-	20522	20542
Foot print	20509	-	20529	20549
Leopard	20511	-	20531	20551
Bubbles	20512	20492	20532	20552
Hearts	20513	-	20533	20553
Caterpillar	20514	20494	20534	-

**15 cm

***14.5 cm

*These powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C

STAINLESS STEEL UTILITY AND BANDAGE SCISSORS



20581

20582

20573

20574

20560

Plastic coloured finger-contour handles for maximum comfort.

All purpose, first aid scissors, easily cut bandages, leather and plaster casts. Autoclavable at 136°C. Instructions: GB, FR, IT, ES, PT, DE.

UTILITY SCISSORS - box of 10			
12 cm - 4.7"	16.5 cm - 6.5"	19 cm - 7.5"	Colour
20560	20570	20580	Black
-	20571	20581	Blue
-	20572	20582	Red
-	20573	20583	Purple
-	20574	20584	Pink
20590	Steel chain for scissors and forceps		



20590