

AIR CUSHION SPLINTS



• **34106 FIRST AID AIR CUSHION SPLINTS** - set of 4 pieces
Indispensable instrument for first aid in case of accidents, cutting and crushing injuries. Prevent surface bleeding without stopping the circulation and secure painless transportation of the injured part. The bandage is a transparent cover of PVC easy to clean and reusable, resistant against water and shifting weather conditions. X-rays can penetrate and medical diagnosis can be made without removing the bandage. The long tube (14.5 cm) with the valve make the inflation of the bandage easy, even by the injured part himself.

Set of 4 different pieces: full leg, foot/ankle, full arm, hand/wrist splints in a 4 compartment bag, with inflating pump.

ALUMINIUM SPLINTS soft lining

Soft lining for comfort.

Aluminium thickness 1.0 mm.


• **34120 SPLINTS** - 15x495 mm

• **34121 SPLINTS** - 20x495 mm

• **34122 SPLINTS** - 25x495 mm



BOSTON SPLINTS



• 34123 BOSTON SPLINTS

Adaptable light-weight splints, x-ray translucent, coated with a non-toxic material which can be easily washed and disinfected. Size cm 92x11. Multilingual instructions: GB, FR, IT, ES, DE.

FINGER SPLINTS

FINGER SPLINTS






Available in 5 different sizes in clear colour.

Made of high quality rigid plastic.

With holes for air circulation and cooling.

Multilingual instructions:

GB, FR, IT, ES, PT, DE, NL, PL, DK, SE, FI, CZ, SK, HU, HR, SI, RO, BG, GR, LT, LV, EE, Arabic.

FINGER SPLINTS				
				
42102	42103	42104	42105	42106
GIMA code	FINGER SPLINTS			Minimum order
42102	Finger splint - size 3			box of 20
42103	Finger splint - size 4			box of 20
42104	Finger splint - size 5			box of 20
42105	Finger splint - size 5.5			box of 20
42106	Finger splint - size 6			box of 20

GIMASPLINTS IN 5 SIZES: LEG, ARM, FOREARM, WRIST AND ELBOW



GIMASPLINTS these supports immobilize your limbs in an efficient and practical way. Each splint can be shaped to fit the patient's personal need. It ensures a high degree of rigid support: thanks to inner part made of neoprene lined with 100% laminated, absorbent-proof nylon assembled to a rigid metal splint. A velcro band allows a quick fixing. Available in 5 sizes or in a bag set of 5 pieces.

GIMA code	GIMASPLINTS
34680	Splint for leg - 57 cm
34681	Splint for arm - 60 cm
34682	Splint for forearm - 39 cm
34683	Splint for wrist - 33 cm
34684	Splint for elbow/ankle - 54 cm
34685	Set of 5 splints in a blue bag
34054	Bag for 5 splints - blue



VACUUM SPLINT KIT

• 34689 VACUUM SPLINT KIT

The set includes vacuum splints for arm (70x50 cm), leg (100x70 cm) forearm (50x30 cm) and a nylon bag. Produced with high quality materials, these splints allow perfect immobilization of the lower and upper limbs by removing the air inside. Equipped with an outside valve, to avoid troubles to the patient, and three velcro straps. Made of vinyl covered nylon, x-ray translucent.

• **34095 ALUMINIUM PUMP** - for vacuum splints and mattress

Capacity: 0.75 l, rotating handle. Supplied with hose and fitting.



Carrying bag included

