


PLASTIC OR ALUMINIUM SCOOPING STRETCHERS

Supplied with multilingual manual: GB, FR, IT, ES, PT, DE



Both stretchers comply with EN1865 (patient handling equipment used in road ambulance)

- 34103 SCOOPING STRETCHER - aluminium
 - 34104 SCOOPING STRETCHER - plastic
- Practical scooping stretchers allow to put the stretcher under the patient without moving him/her much. The shape of the stretchers has been designed to immobilize the patient properly during transport.
- effective hooking and unhooking device for stretcher uncoupling
 - adjustable stretcher length with automatic locking device
 - bag with 3 quick unhooking safety belts for patient immobilization
 - once folded it can be easily stored inside the ambulance
 - built in light and sturdy aluminium (34103) or plastic (34104)
 - X-ray translucent (34104 only)

TECHNICAL SPECIFICATIONS				
	34103		34104	
Open size cm:	min 171x43	max 205x43	min 164x44.5	max 200x44.5
Folded size cm:	121x43x10		118x44.5x8.5	
Weight:	8.5 kg		7.5 kg	
Load:	159 kg		159 kg	


HEAD IMMOBILIZERS



34108

Waterproof plastic coating makes it easy to clean and prevents bacterial growth

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, SE, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic.



34109



Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, GR, BG, RO.



34079

Radiolucent material, suitable for x-ray chamber

HEAD IMMOBILIZERS

Adult and child models of head immobilizers, different in shape, material and colour. All are composed of two head blocks and two straps to prevent any head movement of the injured person during transport. Both sides have an ear hole. They can be used with scooping and PE stretchers.

GIMA code	Type	Colour	Unfolded size (cm)	Weight (kg)	Material
34108	FERMO 3 - child	blue	26x34x14	0.5	high density waterproof plastic
34109	FERMO 3 - adult	red	40x26x18	1	high density waterproof plastic
34079	CONTOUR - adult	yellow	-	1.5	moulded polyurethane and PVC

IMMOBILISATION BELTS FOR STRETCHERS AND SPINAL BOARDS



33840

33842

33843

33844

34081

32794



33847 set of 4



33849

76 cm
90 cm 90 cm 90 cm

Adjustable spider strap for child and adult total length 173 cm

GIMA code	Type	Colour	Size (cm)		Belt material	Buckle	Compatibility
			Width	Length			
33840	Safety belt B2	black	4.8	180	polypropylene	plastic buckle	PE stretchers
33842	Safety belt B4	black	4.8	180	polyester	plastic buckle	PE stretchers
33843	Safety belt B8	black	4.8	130+40*	polyester	plastic buckle, end with hooks	PE stretchers
33844	Safety belt B9	black	4.8	180	polyester	metal buckle	PE stretchers
33847	Set of 4 safety belts	mix	4.8	150+30*	polyester	metal buckle, end with hook	PE stretchers, boards
33849	Adjust. Ragno/Spider	mix	-	-	polypropylene	-	spinal boards
34081	Set 3 belts	orange	-	-	-	quick release plastic buckle	stretchers + boards
32794	Nylon bag - empty	red	-	-	-	-	contains up to 3 belts

* Two parts belt