

## LANCETS AND LANCET DEVICE

### • 23919 AUTO LANCET DEVICE

Reusable auto lancet plastic device, adjustable tip, allows precise configuration of lancet penetration, with 10 adjustable depths. Lancet ejector for safe disposal of the lancet. Instructions: GB, FR, IT, ES.



23919

- 23917 GIMA LANCETS 28G - sterile - box of 25 pcs
  - 23916 GIMA LANCETS 28G - sterile - box of 100 pcs
  - 23918 GIMA LANCETS 30G - sterile - box of 100 pcs
- Fine needles usable with code 23919 and most of lancet devices in the market to draw a blood sample for glucose testing.



Multilingual box: GB, FR, IT, DE, ES, PT, GR, PL, CZ, Arabic

## ABBOTT GLUCOSE STRIPS


| ABBOTT FREESTYLE OPTIUM BLOOD GLUCOSE TEST STRIPS |                                                                                                   |               |
|---------------------------------------------------|---------------------------------------------------------------------------------------------------|---------------|
| GIMA code                                         | Compatibility                                                                                     | Minimum order |
| 23942                                             | FreeStyle Optium Neo, FreeStyle Libre, FreeStyle Optium, Optium Xceed, Optium, Optium Easy, Boots | box of 25     |
| 23943                                             |                                                                                                   | box of 50     |
| ABBOTT FREESTYLE LITE BLOOD GLUCOSE TEST STRIPS   |                                                                                                   |               |
| 23944                                             | FreeStyle Lite, FreeStyle Freedom Lite, FreeStyle InsuLinx, FreeStyle Navigator II                | box of 25     |
| 23945                                             |                                                                                                   | box of 50     |



## GIMALANCET SQUARE AND ROUND SAFETY LANCETS - STERILE

Easy to use, comfortable, safety lancets to meet various capillary blood requirements.

- **single-use:** the build-in-self-destructing design can prevent reuse of the device
- **no accidental injury risk:** needle is hidden before and after use to prevent accidental injuries
- **sterilization guaranteed:** the lancet is Gamma sterilized and the needle tip is protected with a sterility tab until its use allowing a long shelf life
- **no blood contamination**



Easy to use pre-loaded and pressure activated device

High quality tri-bevel needle high-speed penetration with single spring

| GIMA ROUND SAFETY LANCETS - NEEDLES  |            |                      |       |            |              |                    |          |               |
|--------------------------------------|------------|----------------------|-------|------------|--------------|--------------------|----------|---------------|
| GIMA code                            | Colour     | Depth of penetration | Gauge | Blood Flow | Finger Stick | Glucose Hematorrit | Bleeding | Minimum order |
| 24451                                | Yellow     | 2.2 mm               | 21G   | 40-60 µl   | •            | •                  | High     | box of 100    |
| 24453                                | Green      | 2.2 mm               | 23G   | 40-60 µl   | •            | •                  | High     | box of 100    |
| 24454                                | Light blue | 1.8 mm               | 26G   | 5-10 µl    | •            | •                  | Medium   | box of 100    |
| 24455                                | Purple     | 1.8 mm               | 28G   | 5-10 µl    | •            | •                  | Low      | box of 100    |
| GIMA SQUARE SAFETY LANCETS - NEEDLES |            |                      |       |            |              |                    |          |               |
| 24461                                | Orange     | 2.2 mm               | 21G   | 40-50 µl   | •            | •                  | High     | box of 100    |
| 24463                                | White      | 1.8 mm               | 23G   | 20-40 µl   | •            | •                  | Medium   | box of 100    |
| 24465                                | Purple     | 1.8 mm               | 28G   | 5-10 µl    | •            | •                  | Low      | box of 100    |

Multilingual manual on [www.medipurpose.com](http://www.medipurpose.com)



Multilingual box: GB, FR, IT, ES, DE, PT, NL, SE, DK, RU, GR.

24515

## SURGILANCE™ SAFETY AUTOMATIC LANCETS - STERILE

### SURGILANCE™ SAFETY AUTOMATIC LANCETS - NEEDLES

Multi purpose capillary blood sampling devices for all testing requirements.

The one-step safety lancet is easier to use, no arming is required, and safer, once lancet is used it is rendered inoperable.

High-speed delivery and penetration method minimizes patient pain and operator error.

Latex-free, hypo-allergenic, disposable, sterile.

| SURGILANCE™ SAFETY LANCETS - NEEDLES |        |                      |                     |                  |            |              |                    |               |
|--------------------------------------|--------|----------------------|---------------------|------------------|------------|--------------|--------------------|---------------|
| GIMA code                            | Colour | Depth of penetration | Entire needle Gauge | Needle tip Gauge | Blood Flow | Finger Stick | Glucose Hematorrit | Minimum order |
| 24513                                | Yellow | 1.0 mm               | 21G                 | 26G              | 5-10 µl    | •            | •                  | box of 100    |
| 24514                                | Grey   | 1.8 mm               | 21G                 | 23G              | 10-20 µl   | •            | •                  | box of 100    |
| 24515                                | Orange | 2.2 mm               | 21G                 | 22G              | 20-40 µl   | •            | •                  | box of 100    |
| 24517                                | Pink   | 2.8 mm               | 21G                 | 21G              | 40-60 µl   | •            | •                  | box of 100    |