

STERILE ABSORBENT PADS



HIGH ABSORPTION

Multilingual box:
GB, FR, IT, ES,
PT, DE, NL, SE

Ideal for situations where a great capacity of absorption is needed, such as bleedings and chronic wounds. Made of cellulose pulp, white and blue polypropilene and tissue paper.

GIMA code	STERILE ABSORBENT PADS	Minimum order
59402	10x10 cm	box of 60
59403	10x20 cm	box of 30
59404	20x20 cm	box of 30
59405	20x40 cm	box of 25

DENTAL COTTON ROLLS

• 35000 COTTON ROLLS 1x3.8 cm
100% cotton interior, soft, pliable, non-linting, non-irritating rolls. Highly absorbent up to 6 times their weight. Will not adhere to mucous membranes, maintains shape better when wet, non sterile. Minimum order: carton of 10 boxes of 1,000 pcs.



35000



100% COTTON WOOL



100% COTTON WOOL

Hydrophilic 100% cotton for sanitary and cosmetic use, hygiene, non-injured skin disinfection. Non-invasive Class I Medical device. Non sterile. Multilingual pack: GB, FR, IT, DE, ES, FI, GR, HR, PL, PT, RO, SE.

GIMA code	COTTON WOOL	Type	Closing laces	Minimum order
34831	Cotton wool 100 g	Z folded	yes	box of 100 packs
34834	Cotton wool 250 g	Z folded	yes	box of 40 packs
34835	Cotton wool 500 g	Z folded	-	box of 20 packs
34838	Cotton wool 1000 g	roll	-	box of 10 packs
34839	Cotton wool 1000 g	Z folded	-	box of 10 packs

SURGIFIX® ELASTICATED TUBULAR NETTING

Surgifix® is a unique elastic tubular netting, light, soft and breathable. It is the ideal solution for gauze or dressing retention. Its unique structure makes it suitable for any part of the body, thus assuring perfect dressing adhesion under ideal compression grade and offering unrivalled patient's comfort.

Surgifix can be cut at the requested length. Two ranges: Poly (polyester 38%, polyamide 35%, latex 27%) and Zerolatex (polyamide 86%, elastomer 14%).



Multilingual box (GB, FR, IT, ES, DE) and instructions (GB, FR, IT, ES, PT, DE, HU, PL, BG, GR, RO, SE, SI, NO, NL, LT, LV, CZ, DK, FI, RU, TR, SK, EE, JP, Arabic).

MADE IN ITALY



POLY			SURGIFIX® ELASTICATED TUBULAR NETTING Length under tension: 25 m	ZEROLATEX		
Code	Range	Width* mm		Code	Range	Width* mm
34701	1	18	Net for fingers	-	-	-
34702	2	20	Net for wrist, hand	-	-	-
34703	3	22	Net for hand, elbow, arm	58903	C	38
34704	4	28	Net for elbow, arm, foot	58904	C/D	48
34705	5	30	Net for foot, knee, leg	58905	E	62
34706	6	56	Net for head, shoulder, pelvis, leg	58906	F	78
34707	7	61	Net for body, pelvis, thigh	58907	G	88
34708	8	65	Net for shoulder, body, back	-	-	-
34709	9	85	Net for shoulder, body, back	-	-	-
34700	10	114	Net for shoulder, body, back	-	-	-

*Different width between Poly and Zerolatex is due to different elasticity of materials