

## ZINCA BAND E



- 34715 ZINCA BAND E ELASTIC BANDAGE 5 m x 10 cm
  - 34717 ZINCA BAND E ELASTIC BANDAGE 10 m x 10 cm
- Elastic bandage impregnated with zinc oxide paste. Support for treatment of leg ulcers, varices, thrombosis and post-care for fracture, dislocations, bone and joint surgery. Each band is packed in sealed aluminium foil bag and multilingual box: GB, IT. Medical device.

### TEXTILE COMPOSITION

Viscose 69% - Polyamide 31% - Soaked with zinc oxide, gelatine, glycerin, propyl-hydroxybenzoate, methyl-hydroxybenzoate, water. Elasticity: 75%

## PREVIZINC A and E



Previzinc "A" and "E" are bandages impregnated with zinc oxide paste indicated in the treatment of thrombophlebitis, oedema, chronic venous insufficiency and in the healing phase and follow-up care of fractures, sprains and strains.

GIMA code	PREVIZINC "A" INELASTIC BANDAGES	Minimum order
34714	Previzinc "A" 6 m x 8 cm	box of 10
34716	Previzinc "A" 5 m x 10 cm	box of 10
PREVIZINC "E" ELASTIC BANDAGES		
34718	Previzinc "E" 6 m x 8 cm	box of 10
34719	Previzinc "E" 5 m x 10 cm	box of 10
34724	Previzinc "E" 7 m x 10 cm	box of 10

### TEXTILE COMPOSITION

Previzinc A: 100% cotton zinc oxide paste

Previzinc E: 54% polyamide, 46% viscose zinc oxide paste

## 3M™ RESTON™ FOAM PAD

Self-adhering open cell foam designed to protect skin from splints, casts and prostheses. Not intended to be used as a pressure reducing device. It can also be used as a padding for beds and chairs. Can be cut, layered and shaped.



GIMA code	3M code	3M RESTON™ FOAM PAD	Minimum order
34920	1560	RESTON™ 30x20x1.25 cm	box of 10
34921	1561	RESTON™ 30x20x2.5 cm	box of 5
34922	1563	RESTON™ ROLL 10cmx5mx5mm	1 roll



## PLATRIX® HIGH QUALITY PLASTER OF PARIS BANDAGES



Multilingual box:  
GB, FR, IT, ES, DE, NL

### PLATRIX® PLASTER OF PARIS

High quality plaster of Paris bandages. Economic thanks to the high strength and the minimal plaster loss.



GIMA code	BSN code	PLATRIX® PLASTER OF PARIS	Minimum order
34686	71996 - 05	Platrix® 2 m x 15 cm	box of 40
34688	71996 - 14	Platrix® 3 m x 20 cm	box of 20
34691	71996 - 28	Platrix® 4 m x 20 cm	box of 20

## OSCILLATING PLASTER SAW

### • 59099 HEBU GOLD II PLASTER SAW

The new oscillating professional saw can quickly, safely, quietly and accurately remove all plaster of Paris and synthetic dressings.

High-performance motor with electronic speed regulation and a gradual start-up, gently on materials. Motor is cooled using the circulation principle.

Cool air is suctioned in at the end of the motor, to protect power aggregate from resin dust and residues. Made of high quality materials.

- 59101 SAW BLADE Ø 50 mm
- 59102 SAW BLADE Ø 65 mm



59099




MADE IN GERMANY

UP TO 21,000 OSCILLATIONS

### STANDARD ACCESSORIES

Saw blade Ø 50 mm  
Saw blade Ø 65 mm  
One pair of wrenches  
User manual: GB, FR, IT, ES, DE



59101

59102

### TECHNICAL SPECIFICATIONS

Voltage: 220-240 V, 50-60 Hz  
Power: 250 W  
Strokes: 12,000 - 21,000 oscillations/min  
Length without cord: 340 mm  
Weight with/without cable: 1.3/1.7 Kg  
Noise: 68-82 dB  
Cord length: 4.5 m  
IP code: IP20