

HYPOR FILM ADHESIVE WATERPROOF DRESSING FOR FIXATION

Sterile adhesive dressing on transparent polyurethane film apt to fix catheters, gauzes, bandages and other dressings. The films acts as water and bacterial barrier. Permeable to air and water vapor, lets the skin breathe and avoids the risk of maceration. Individually packed. Extremely elastic substrate and highly conformable to sensitive areas. Reduces the risk of fluid leaking and infection. Multilingual box: GB, FR, IT, ES, DE. Medical Device Class I.

GIMA code	HYPOR FILM DRESSING FOR FIXATION - STERILE	Minimum order
58840	Hypor film 6 x 8 cm	box of 50 pcs
58841	Hypor film 10 x 12.5 cm	box of 50 pcs
58842	Hypor film 10 x 20 cm	box of 25 pcs
58843	Hypor film 10 x 25 cm	box of 25 pcs

CE 0546



Rounded edges to prevent the dressing from peeling off

Hypoallergenic acrylic adhesive, dermatologically tested

Protective liners easy removable with gloves

MADE IN ITALY

ADHESIVE WOUND DRESSING I.V. FIXATION DRESSING



CE 1370

Soft transpirant non woven fabric coated with hypoallergenic adhesive polyacrylate. Central absorbent, non adherent pad in folded non woven fabric. X-ray transparent

GIMA code	STERILE - NON WOVEN ADHESIVE WOUND DRESSING HYPOALLERGENIC	Minimum order
34860	Sterile adhesive dressing 5x7 cm	box of 100 pcs
34861	Sterile adhesive dressing 8x10 cm	box of 50 pcs
34862	Sterile adhesive dressing 8x15 cm	box of 50 pcs
34849	Sterile adhesive dressing 10x25 cm	box of 25 pcs

• 34845 STERILE ADHESIVE DRESSING for cannula fixation - box of 50 pcs
Designed for wide applications. Hypoallergenic acrylic adhesive, dermatologically tested. Size 8x5.8 cm. Wound pad size 2x1.5 cm. Multilingual box: GB, FR, IT, ES, DE. Medical Device Class I.

MADE IN ITALY



CE 0546

34845

TAILORED BANDAGES - SENSITIVE SKIN

Hypoallergenic extensible non woven adhesive dressing to be cut, with highly absorbent central wound pad. Recommended for sensitive skin. Multilingual box: GB, FR, IT, ES, DE. Medical Device Class I.

MADE IN ITALY

GIMA code	TAILORED BANDAGES	Wound pad	Minimum order
34855	1 m x 6 cm	1 m x 2.4 cm	1 roll
34856	5 m x 6 cm	5 m x 2.4 cm	1 roll



34856

SURGICAL ADHESIVE DRESSINGS BWT ADHESIVE DRESSING



CE 0546



Multilingual box: GB, FR, IT, ES, DE.

Sterile hypoallergenic adhesive non woven dressing with central non-wound-sticking absorbant pad. Medical Device Class I.

GIMA code	NON WOVEN DRESSING PHARMACY LINE - STERILE	Wound pad	Minimum order
34818	Surgical dressing 5x7.2 cm	2.8x3.8 cm	box of 5 pcs
58822	Surgical dressing 10x8 cm	5x4 cm	box of 5 pcs
34819	Surgical dressing 10x10 cm	5x5 cm	box of 5 pcs
34820	Surgical dressing 8x15 cm	3.8x10.5 cm	box of 5 pcs
NON WOVEN DRESSING HOSPITAL LINE - STERILE			
58830	Surgical dressing 5x7.2 cm	2.8x3.8 cm	box of 50 pcs
34853	Surgical dressing 10x8 cm	5x4 cm	box of 50 pcs
58832	Surgical dressing 10x10 cm	5x5 cm	box of 50 pcs
34822	Surgical dressing 8x15 cm	3.8x10.5 cm	box of 50 pcs
58834	Surgical dressing 10x20 cm	5x15 cm	box of 50 pcs
58835	Surgical dressing 10x25 cm	5x20 cm	box of 25 pcs



WATERPROOF

BWT material: transpiring, waterproof, thin

Round corners avoid dressing from coming off easily

CE 0546

Polyurethane film acts as an antibacterial barrier

MADE IN ITALY

Sterile hypoallergenic adhesive dressing with central non-wound-sticking absorbant pad. Sterile impermeable protection of surgical wounds, superficial wounds, cuts and abrasions. Multilingual box: GB, FR, IT, ES, DE. Medical Device Class I.

GIMA code	DRESSING PHARMACY LINE - STERILE	Wound pad	Minimum order
58800	BWT dressing 5x7.2 cm	2.8x3.8 cm	box of 5 pcs
58801	BWT dressing 10x6 cm	6.5x3.5 cm	box of 5 pcs
58802	BWT dressing 10x8 cm	5x4 cm	box of 5 pcs
58803	BWT dressing 10x10 cm	5x5 cm	box of 5 pcs
58804	BWT dressing 10x12.5 cm	5x7.5 cm	box of 5 pcs
58805	BWT dressing 8x15 cm	3.8x10.5 cm	box of 5 pcs
DRESSING HOSPITAL LINE - STERILE			
34866	BWT dressing 5x7.2 cm	2.8x3.8 cm	box of 50 pcs
34867	BWT dressing 10x8 cm	5x4 cm	box of 50 pcs
34868	BWT dressing 8x15 cm	3.8x10.5 cm	box of 50 pcs