

## 3M™ INADINE™

### PVP-I NON ADHERENT DRESSING

Non adherent dressing, made of low adherent knitted viscose fabric impregnated with a polyethylene glycol (PEG) base, containing 10% Povidone Iodine. Designed to protect wounds even when they are infected, providing long-lasting antimicrobial effect against bacteria, mycobacteria, fungi, protozoa and viruses, including resistant species such as MRSA. Inadine dressing minimizes adherence to the wound bed, therefore reducing pain for patient and the risk of damage to the granulation tissue at dressing removal. Colour evolution from orange to white indicates loss of antiseptic efficacy and dressing needs to be changed. Weight: 3.25g. Indicated for ulcerative wounds, prevention of infection in minor burns, traumatic skin loss injuries and heavily infected wounds.



GIMA code	3M™ code	3M™ INADINE™	Minimum order
59575	P01481	Inadine™ 5x5 cm	box of 25
59576	P01491	Inadine™ 9.5x9.5 cm	box of 10

## 3M™ KERRAMAX™ CARE

### SUPER-ABSORBENT DRESSING

The unique wicking layer distributes fluid both horizontally and vertically throughout the dressing, utilizing the full absorption capacity. Exu-Safe™ Technology, built into the dressing core, is an advanced method of absorbing high level of exudate from wounds, including bacteria and MMPs, creating an optimal environment for wound healing. It protects delicate wound tissue and surrounding, helping to reduce the risk of maceration, even under compression. Its heat-sealed border prevents exudate leakage. It can be left in place for up to 7 days. Absorption capacity: 0.95 g/cm². Indicated for: pressure injuries, venus leg ulcers and diabetic foot ulcers.



GIMA code	3M™ code	3M™ KERRAMAX™ CARE	Minimum order
59580	PRD500-050	Kerramax™ Care 10x10 cm	box of 10
59581	PRD500-100	Kerramax™ Care 13.5x15.5 cm	box of 10
59582	PRD500-120	Kerramax™ Care 10x22 cm	box of 10
59583	PRD500-240	Kerramax™ Care 20x22 cm	box of 10
59584	PRD500-380-B10	Kerramax™ Care 20x30 cm	box of 10

## OCULAR ADHESIVE OCCLUSIVE BANDAGES

100% polyester extensible nonwoven occlusive fabric protection. Suitable for safe and comfortable dressing, the central swab has a protective action and the hypoallergenic acrylic adhesive is suitable for sensitive eye skin area. Individually packed. Medical Device Class I. Multilingual box: GB, FR, IT, ES, DE.



GIMA code	Model	Colour	Size	Wound pad	Minimum order
58750	Regular Kids	patterned	7.2x5.7cm	4.9x3.2 cm	box of 21*
58751	Junior	skin	6.0x4.8 cm	3.8x2.6 cm	box of 21
58752	Regular	skin	7.2x5.7cm	4.9x3.2 cm	box of 21
58753	Large	skin	8.0x5.7cm	4.9x3.2 cm	box of 21

\*Mixed patterns: 7 pcs x 3 patterns

## EYE DRESSING

### • 34848 OCULAR ADHESIVE EYE DRESSING

Sterile, extremely absorbent ocular dressing on hypoallergenic permeable adhesive backing. White soft polyester non woven suitable for a safe and comfortable eye medication. Medical Device Class I. Multilingual box: GB, FR, IT, ES, DE.



GIMA code	Size	Wound pad	Minimum order
34848	6.5x9.5 cm	3.8x6.6 cm	box of 10

## 3M™ ADAPTIC™ DIGIT/TOE DRESSINGS

Non-adherent dressing that combines a low tack silicone coated material with a tubular bandage specifically designed for the treatment of finger and toe wounds. Conformable and comfortable dressing, which allows for flexibility and maximum movement of the fingers and toes. Easy to apply, helps minimize trauma associated with dressing changes. It can be left on for several days.



GIMA code	3M™ code	3M™ ADAPTIC™ DRESSING - NON ADHERENT - STERILE	Minimum order
34994	MAD013	ADAPTIC™ DIGIT - for finger medium	box of 10
34995	MAD023	ADAPTIC™ DIGIT - for finger large	box of 10
34996	MAD062	ADAPTIC™ TOE - for toe	box of 10

## SINGLEFIX® FINGER DRESSING

- 34800 SINGLEFIX® FINGER DRESSING - small - box of 100
- 34801 SINGLEFIX® FINGER DRESSING - medium - box of 100
- 34805 SINGLEFIX® FINGER DRESSING - adult - box of 100

Finger dressing consisting of a cotton sheath and soft inner absorbent padding. Convenient wrist knot fixing. Easy to wear, it protects the injured finger, offering protection and comfort.

Multilingual instructions: GB, FR, IT, ES, PT, DE, HU, PL, BG, GR, RO, SE, SI, NO, NL, LT, LV, CZ, DK, FI, RU, TR, SK, EE, JP, Arabic.

