

3M™ SOLUTIONS FOR CHRONIC WOUNDS

TEGADERM™ HYDROCOLLOID

Sterile hydrocolloid dressing for low to moderate draining wounds. Also available in sacral shape with highly breathable waterproof border.



34964

TEGADERM™ HYDROGEL

Preservative-free hydrogel wound filler, easy to apply. Formulated to help providing a moist wound healing environment. Single use tube.



34969

TEGADERM™ ALGINATE

High integrity, gel forming alginate dressing for moderate to heavily exuding wounds. Compact and comfortable.



34972

GIMA code	3M code	3M™ TEGADERM™	Minimum order
34969	90007	HYDROCOLLOID™ 16x17 cm sacral	box of 6
34964	90112	TEGADERM™ ALGINATE - 10x10 cm	box of 10
34972	91110	TEGADERM™ HYDROGEL 15 g	box of 10



34954

34957



34955

3M

3M™ TEGADERM™ HIGH PERFORMANCE FOAM

Sterile, multi-layer foam dressing that combines high absorbency with high breathability to reduce the risk of wound maceration.

Well suited for changing levels of exudate to maintain moisture balance for optimal wound healing.

Highly breathable and waterproof, also available in heel/elbow shape with adhesive border and simplified application system

GIMA code	3M code	3M™ TEGADERM™ FOAM - non adhesive or adhesive dressing	Minimum order
34954	90601	TEGADERM™ FOAM - 10x10 cm	box of 10
34955	90604	TEGADERM™ FOAM - 8.8x8.8 cm for tracheostomy	box of 10
34957	90619	TEGADERM™ FOAM adhesive for heel and elbow - 14x14 cm	box of 5

HYPOR STERILE POLYURETHANE FOAM DRESSING

HYPOR STERILE POLYURETHANE FOAM DRESSING

Sterile hydrophilic polyurethane foam is indicated for the treatment of non-infected chronic or acute wounds, moderate to heavily exuding.

The dressing may be used during the healing process on pressure ulcers, leg and foot ulcers, diabetic ulcers and surgical wounds.

- high vertical absorption and retention capacity - minimize the risk of maceration
- semipermeable barrier - lets the skin breathe
- water and bacteria barrier - reduces the risk of infection
- soft, elastic and extremely conformable



Hypor Foam can stay in place for 3-4 days and up to a maximum of 7 days.

Individually packed. Medical Device Class IIB

GIMA code	HYPOR FOAM HEEL - NON ADHESIVE	Minimum order
58850	HYPOR FOAM heel	box of 5
	HYPOR FOAM PAD - NON ADHESIVE	
58852	HYPOR FOAM pad 5x5 cm	box of 10
58853	HYPOR FOAM pad 10x10 cm	box of 10
58854	HYPOR FOAM pad 15x15 cm	box of 10



58850



58852



58856

Anatomical shape. Quick and simple to use. Provides cushioning and pressure relief.

MADE IN ITALY

Provides cushioning protection

GIMA code	HYPOR FOAM DRESS - ADHESIVE	Size	Wound pad	Minimum order
58856	HYPOR FOAM DRESS	7.5x7.5 cm	4.5x4.5 cm	box of 10
58857	HYPOR FOAM DRESS	12.5x12.5 cm	7.5x7.5 cm	box of 10

NaCl 0.9% B-BRAUN STERILE IRRIGATION SOLUTIONS



GIMA code	B-BRAUN code	B-BRAUN	Minimum order
36550	3570310	B-BRAUN ECOLAV® - 30 ml	box of 100
36551	3570350	B-BRAUN ECOLAV® - 100 ml	box of 20
36552	3570380	B-BRAUN ECOLAV® - 250 ml	box of 20
36553	3570410	B-BRAUN ECOLAV® - 500 ml	box of 10
36554	3570110	B-BRAUN ECOTAINER® - 250 ml	box of 12
36555	3570130	B-BRAUN ECOTAINER® - 500 ml	box of 10
36558	3513576	B-BRAUN OMNIFLUSH® - 10 ml	box of 100

ECOLAV®

The spray nozzle of Ecolav® allows a precise and direct irrigation of the target area by simply applying a slight pressure on the container, helping optimize the time for irrigation processes. Multilingual label: GB*, FR, IT, ES, DE, PT, NL. *Excluded 36551

ECOTAINER®

Suggested for quick irrigation with high volume of solution. Multilingual label: GB**, FR, IT, ES, DE, PT, NL, DK**, FI**, NO**, SE**. **Only 36554

COMPOSITION - ECOLAV and ECOTAINER

1,000 ml solution contain: sodium chloride 9 g, water, injection electrolytes: mmol/l Na⁺ 154 Cl⁻ 154

B-BRAUN OMNIFLUSH®

Pre-filled flush syringe 10 ml with NaCl 0.9% - EM-3513576

The Omniflush® syringe is designed specifically to minimize unintended blood reflux into the catheter lumen thus reduces the risk of catheter blockage. With Omniflush® the fluid amount used for flushing can be reduced to the required minimum.

For flushing of all intravenous accesses.

Multilingual manual and label:

GB, FR, IT, ES, DE, PT, PL, CZ, BG, DK, EE, FI, GR, HU, LT, LV, NL, NO RO, SK, SI, SE, TR.



36558



36558