


SURGICAL SCRUB BRUSHES/SPONGES

MADE IN ITALY



- **36570 SCRUB KIT - sterile - dispenser of 30 pcs**
Sterile hand-cleaning hospital kit. It includes a scrub, a sponge and a plastic nail cleaner. Dry, without detergent. Sterilization method: ETO. Multilingual label (blister and dispenser): GB, FR, IT, ES, DE, PT, PL, NL, GR, RU, Arabic.



- **36571 IODO-POVIDONE SCRUB KIT dispenser of 30 pcs - IT**
Italian Health Ministry registration PMC 19531
- **36572 IODO-POVIDONE SCRUB KIT dispenser of 30 pcs - Multilingual***
*Label: GB, FR, ES, DE, PT, NL, GR (blister and dispenser), PL, RU, Arabic (blister only), IT (dispenser only)
same as 36570 with solution of PVP Iodine 7.5 g in 100 ml, not sterile.



- **36573 CHLOREXIDINA SCRUB KIT dispenser of 30 pcs - IT**
Italian Health Ministry registration PMC 19575
- **36574 CHLOREXIDINA SCRUB KIT dispenser of 30 pcs - Multilingual***
*Label: GB, FR, ES, DE, PT, NL, GR (blister and dispenser), PL, RU, Arabic (blister only), IT (dispenser only)
same as 36570 with solution of Chlorhexidine Gluconate 20% 4 g in 100 ml, not sterile.

All scrub kits are individually blistered in 30 pieces dispenser



36570



36571



36572



36573



36574

PERFECTION BRUSHES

STERILIZABLE IN AUTOCLAVE



- **25715 PERFECTION BRUSHES**
Nylon brushes for cleaning before and after operation. These brushes are autoclavable (at 100°C for not more than 5 minutes) and can be stored in any standard brush dispenser in hospitals. Resistant and reliable. Handle in PP. Size: 109x41x34 mm. Supplied in boxes of 12 units.

SURGERY DISINFECTANT

- **36632 NEO STERIXIDINA SOAP - 1 l**
 - **36690 DISPENSER PUMP**
Ready-to-use disinfectant liquid soap for hands and skin with antibacterial activity, based on Chlorhexidine Digluconate 4%.
- It can be used for surgeon and medical staff's hands and also patient disinfection, before surgical operations. It reduces the risk of cross-infection. It removes microorganisms and keeps bacteria flora at lowest level. Label in GB, IT.

4% CHLORHEXIDINE

MADE IN ITALY

CHEMICAL COMPOSITION

100 g contain:
Chlorhexidine digluconate 4 g - Coformulants
- Depurated water as needed for 100g
Italian Health Ministry registration 18980

36632
+ 36690

CHLORHEXIDINE DISINFECTION LINE

36634

- **36637 GERMOXID SPRAY DISINFECTION - 100 ml - box of 12**
- **36635 GERMOXID DISINFECTION - 250 ml - box of 12**
- **36636 GERMOXID DISINFECTION - 1,000 ml - box of 12**
- **36634 GERMOXID WIPE - 12 boxes of 10 wipes**
- **36638 GERMOXID WIPE - box of 400 wipes**

Ready-to-use disinfectant for skin that contains Chlorhexidine with long-lasting bactericidal action. Chlorhexidine has a proven broad-spectrum effectiveness against Gram⁺ and Gram⁻ bacteria. Thanks to its structure, Chlorhexidine can be easily absorbed by the skin. Safe for medical operators thanks to its high epidermis compatibility, it neither irritates the skin nor stains or burns. Available in 3 versions: liquid (250 or 1,000 ml), spray (for quick and efficient action) or wet wipes with chlorhexidine.

MADE IN ITALY

CHEMICAL COMPOSITION

100,0 g of solution contain: Chlorhexidine digluconate 1,0 g;
soluble bergamot essential oil and depurated water q.s. to 100,0 g.
Contain: limonene, Linalool.
Italian Health Ministry registration 11571



36638

14x20 cm
wipe
suitable
for hands
and skin
disinfection
(bag 6x8 cm)



36636

36635

36637

Label in
GB, IT