

DIGITAL PET SCALES FOR LARGE AND SMALL ANIMALS

SOEHNLE


• 80220 SOEHNLE 6858 VETERINARY SCALE

Digital vet scale with large and flat stainless steel weighing platform. Designed for small and medium-sized animals. Quick and easy operation. Easy to read LCD display with 25 mm high digits and aluminium housing. Auto-off function. Hold function and taring function. Supplied with rubber mat and multilingual manual: GB, FR, IT, DE.

TECHNICAL SPECIFICATIONS

Capacity:	200 kg	Sensitivity:	50 g
Display:	25 mm backlit LCD	Weighing units:	kg, lb, oz, lb.oz
Size:	platform 900x600xh 40 mm - indicator 215x100xh 70 mm		
Power supply:	DC 5V power supply (included) or 4x 1.5 V AA batteries (not included)		



• 27258 DIGITAL PET SCALE - medium

MADE IN EUROPE

• 27259 DIGITAL PET SCALE - large

Precise full electronic multifunctional scale, designed for measuring your pet's weight available in two sizes. Quality product with very competitive price. Switchable kg, lb, st. Minimal measuring weight: 3 kg. Automatic switch-off. Multilingual manual: GB, FR, IT, ES, PT, DE, PL, CZ, SE, FI, RO, GR, HR, BG.

TECHNICAL SPECIFICATIONS

Capacity/Sensitivity:	100 kg / 100 g
Display:	LCD, 3"
Dimensions:	355x635 mm (27258) - 620x830 mm (27259)
Lithium battery:	included, CR2032, 3 V

SOEHNLE

MADE IN GERMANY

Taring functions
switchable to kg/lb

• 27266 SOEHNLE 8320 FOLDABLE SCALE

• 27271 CARRYING BAG

Foldable scale with magnetic lock for easy transport.

Manual in DE, GB, FR, ES, IT, GR, CZ.

TECHNICAL SPECIFICATIONS

Capacity:	15 kg	Sensitivity:	10 g	Weight:	1.9 kg
Dimensions:	395x315x105 mm (closed) - 740x315x53 mm (open)				
Alkaline battery	2x1.5V AAA - not included				
Warranty:	2 years				

Tray 530x290 mm

Overload and
low battery
indication

Lb/kg/st
switchable
by user

LCD display 25 mm

Automatic
switch-off

Tare and
automatic hold
functions

Multilingual manual:
GB, FR, IT, ES, SE, FI,
CZ, RO, GR, HR.

• 27261 DIGITAL SMALL PET SCALE

TECHNICAL SPECIFICATIONS

Capacity/Sensitivity:	20 kg / 5 g
Minimum measurable weight:	20 g
Tray size:	530x290 mm
Lithium battery:	CR2032, 3 V (not included)

DIAGNOSTIC AND OPERATION VETERINARY TABLES



Pack size: 121x15x90 cm
Size: 120x58xh 85 cm
Weight: 23 kg
Load: 60 kg

MADE IN ITALY

• 80300 VETERINARY EXAMINATION TABLE

Veterinary operating table with enamelled tubular steel frame, with AISI 304 15/10 stainless steel top with hole and a liquid tray holder underneath. Lateral hook to immobilize animals.



Size: 118x56xh 85/102
Weight: 48 kg
Max load: 90 kg

MADE IN ITALY

• 80301 VETERINARY OPERATION TABLE

Veterinary operating table with electric lift (82/102 cm), baked enamelled steel tubular frame and AISI 304 stainless steel top (120x56 cm) with hole and a liquid bucket holder underneath. Manual Trendelenburg 15° and Antitrendelenburg up to 15°. Lateral hook to immobilize animals.