

MD SCOPE VET VIDEO OTOSCOPE

• 32164 VET MD SCOPE VIDEO OTOSCOPE WITH 3 PROBES

Instant View

A clear and sharp image can be acquired and displayed instantly on a built-in screen. Just press a single button and wait for 3 seconds.

Real Time Interaction

The medical image can be captured and displayed simultaneously on TV monitor or downloaded to Mac/PC system.

Electronic medical record

A medical image is automatically generated to support most of the stand-alone health information system. Probes and accessories at page 254.

• 32175 CRADLE

TECHNICAL SPECIFICATIONS

2.4" colour LCD display
Video out port
Size/weight: 8x20 cm (2.5"x8") - 218 g
Battery lifespan: 4 hours



STANDARD ACCESSORIES

- 3 camera probes 75/150/300 mm
- image dongle
- cable set
- 2 AA alkaline batteries

32175



HEINE VET OTOSCOPE

• 80500 BETA OTOSCOPE SET 2.5 V - G-112.10.118

Set complete with slit illumination head, 1 set long closed metal specula (57x4,65x6,90x7 mm), battery handle, 1 spare bulb, hard case.

• 31772 HEINE 077 OTOSCOPE BULB 2.5 V - spare

• 31781 HEINE 069 OPHTHALMOSCOPE BULB 2.5 V - spare



80500



VET OTOSCOPES, OPHTHALMO AND DIAGNOSTIC SETS



31442

Open otoscope head is supplied with codes 31422, 31424, 31425

31429



31424

• 31424 XENON HALOGEN DIAGNOSTIC SET

Includes: xenon halogen open otoscope head+handle, xenon-halogen ophthalmoscope head, kit of 5 reusable plastic cannulas (31451), metal tongue depressor holder illuminator, nasal specula.

Rechargeable battery handle 2.5 V

All sets have 2.5 V handle for use with 2xC batteries or rechargeable batteries (34467)



31422

• 31422 XENON HALOGEN OTOSCOPE SET

Includes: xenon-halogen open otoscope head+handle, kit of 5 reusable plastic cannulas (31451), kit of 3 reusable plastic cannulas (31447)



31425

• 31425 XENON HALOGEN DIAGNOSTIC SET

Includes: xenon-halogen open otoscope head+handle, xenon-halogen ophthalmoscope head, kit of 3 reusable metal cannulas + 2 metal specula (31453)

• 31447 PARKER VETERINARY OTOSCOPE

Includes: kit of 3 reusable plastic cannulas (31447)



31447



31451

Suitable for small animals



31453

For open otoscope head



31426

Slit illumination head



Easy adjustable with finger clip

31454

GIMA code	VETERINARY SETS AND ACCESSORIES
31451	Kit of 5 reusable plastic cannulas (Ø 2.5, 3.5, 4.5, 5.5, 9 mm)
31447	Kit of 3 reusable plastic cannulas (Ø 4, 5, 7 mm)
31453	Kit of 3 reusable metal cannulas (Ø 4, 5, 6 mm) + 2 slotted metal specula (Ø 7.5, 9.4 mm)
31454	Chromed brass nasal specula 79 mm
31478	Xenon halogen otoscope bulb 2.5 V
31479	Xenon halogen ophthalmoscope bulb 2.5 V
31449	Standard otoscope bulb vacuum 2.5 V for 31442
31542	Charging station
34467	Rechargeable Li-ion battery 2.5 V / 3000 mAh
31429	Swivel lens 2X for Vet oto sets 31422, 31424, 31425
31426	Slit illumination head xenon/halogen for 31422, 31424, 31425