

VET TOURNIQUET

• 80150 QUICK RELEASE VET TOURNIQUET - grey

Quick release VET tourniquet for small animals.

Extra slim, ideal for cats and small dogs.

- easy to put on
- practical handling
- smooth and elastic cord
- plastic cord clamp

Multilingual instructions: GB, FR, IT, ES, PT, DE, PL, RO, GR, SE



TECHNICAL SPECIFICATIONS

Cord clamp size:	43x33x8 mm
Cord size:	200 mm
Weight:	10 g

COMBISCREEN URINE STRIPS

• 24100 VETERINARY URINE STRIPS -

11 parameters - bottle of 100 strips

For examination of

glucose, ascorbic acid, ketones, protein, ph-value, blood, nitrite, leukocytes, specific gravity, bilirubin, urobilinogen. Each box includes a manual for interpretation of results of the following animals: horse, pig, cattle, sheep, goat, dog, cat, bunny and guinea pig.

Multilingual box and instructions: GB, FR, IT, ES, PT, NL, DE.

ASCORBIC ACID PROTECTED



MADE IN GERMANY

Performs highest quality standard

ABSORBENT PADS FOR PETS

• 80900 ABSORBENTS PAD FOR PETS 60x60 cm - box of 10

• 80901 ABSORBENTS PAD FOR PETS 60x90 cm - box of 10

Pads designed for pets with superabsorbent polymer, contained in the inner pulp layer, to retain fluids and help alleviate the onset of bad smells. Practical adhesive tapes positioned on the back side of the pad, to ensure strong fixation to surfaces avoiding any pad shifting.

Suitable for physiological needs or involuntary losses of organic liquids due to incontinence or other pathologies.

- anti-odour control
- 100% PEFC certified cellulose
- non-toxic & waterproof back sheet
- hygienic, practical and resistant

Components: fluff pulp PEFC, tissue PEFC, top sheet made of non-woven polypropylene, back sheet made of PE FILM

TECHNICAL SPECIFICATIONS

	80900	80901
Weight:	54.3 g	84.2 g
Absorption capacity according to ISO11948-1:	1100 g	1700 g

MADE IN ITALY



80901

80900

KORA VET MAGNETOTHERAPY

• 80770 KORA VET MAGNETOTHERAPY BASIC KIT

• 80771 KORA MAT MATTRESS 64x40 cm

Kora is a veterinary device designed for pulsed electromagnetic field (PEMF) therapy.

Ideal for reduction of pain and inflammation in pets.

Suitable for the treatment of neuropathy, ulcers and diseases of the musculoskeletal system of cats and dogs of any size.

Easy to use programmes guarantee a safe and effective therapy for pets. We suggest the use of a single Kora Mat mattress for small and medium-sized animals and use of two mats for large or medium-sized animals, combinable through velcro and connectable to a single device.



80770

Kora device



Cover can be removed for cleaning



STANDARD ACCESSORIES
Carrying bag
User manual: GB, IT, ES, PT
Magnet for verifying therapy operation



Combination of 2 mattresses through velcro for large or medium-sized animals

80770

+ 80771

TECHNICAL SPECIFICATIONS

Channels:	2
Display:	colour display
Magnetic field intensity:	up to 100 G
Frequency:	up to 75 Hz
Programmes:	18
Power supply:	15VDC-2A
Size:	180x110x50 mm

18 PROGRAMMES

Analgesic	Sciatic nerve
Arthritis	Neuropathy
Arthrosis	Prostheses
Muscular atrophy	Osteonecrosis
Discopathy	Osteoporosis
Fractures	Treat. 1 Hz
Healing	Treat. 18 Hz
Ligament lesion	Treat. 50 Hz
Tendon lesion	Treat. 75 Hz