

STERILE MEDICATION, SUTURE, CATHETERIZATION KITS




• 26930 DRESSING KIT 1 - sterile




• 26929 DRESSING KIT 2 - sterile



• 26933 CATHETERIZATION KIT - sterile



• 26936 SUTURE REMOVAL KIT 1 - sterile



• 26945 SUTURE REMOVAL KIT 2 - sterile



• 26932 SUTURE REMOVAL KIT 3 - sterile

KIT CONTENT	Dressing		Suture		Suture removal			Catheters
	26930	26929	26931	26949	26932	26936	26945	26933
Field cloth	•		•		•			•
Cushing surgery forceps			•					
Needle Holder			•					
Kocher forceps			•	•				
Anatomy forceps	•		•	•	•		•	
Plastic forceps			•	•		•	•	
Haemostatic flow clamp		•	•					
Ring shaped forceps	•	•	•		•			•
Styptic clamp (Kelly forceps)			•					
Haemostatic clamp			•	•				
Straight scissors			•					
Adson anatomy forceps 12 cm		•	•	•				
Light green plastic scissors		•	•	•		•		
Littauer scissors					•			
Gauze swabs 7.5x7.5 cm 8 ply	•	•	•	•	•			
Gauze swabs 10x10 cm 8 ply			•	•	•			
4 or 5 Gauze plugs	•		•	•	•			
3 Gauze plugs		•						•
Gauze swabs 10x20 cm - 16 plies		•						
Gauze swab 18x40 cm						•		
Latex gloves - medium								
Vinyl gloves			•		•			•
Disposable big bag			•		•			•
Medication plaster 15x12 cm					•			
WNW (TNT)+PE drape 32x50 cm		•						
WNW (TNT)+PE drape 50x60 cm			•	•				
WNW (TNT) compresses								•
WNW (TNT)+PE drape 50x60 cm with hole			•	•				•
Water syringe								•
Single dose lubricant								•
2 disinfectant wipes							•	•
2 adhesive WNW (Woven Non Woven) strips				•				
Waste bag		•						
Scalpel with handle							•	
2 cotton applicators							•	
PE bag minigrip				•				



• 26931 SUTURE KIT 1 - sterile



• 26949 SUTURE KIT 2 - sterile

Labellinstructions in GB, IT for 26930-1-2-3
Labellinstructions in GB, FR, IT, ES for 26929, 26936,
26945, 26949

GENERAL SURGERY PROCEDURE PACK



26985

• 26985 GENERAL SURGERY PACK - sterile

Ready to use pack including all necessary for a general surgery procedure.

Peel open PE + tyvek envelope.

Sterile EO.

KIT CONTENT

1 table cover 160x200 cm	4 hand towels
2 surgical gowns XL (1 standard, 1 reinforced)	1 bowl 120 ml
2 adhesive drapes 75x90 cm	1 kidney dish
1 foot drape (ADH side) 150x180 cm	1 scalpel n°15
1 head drape (ADH side) 150x240 cm	5 gauze swabs 10x20 cm
2 syringes 10 ml	4 gauze balls 30 mm
2 needles 21Gx1.5"	2 foam sticks 15 cm
1 scalpel n°11	1 basic tray