

# SURGICAL NON ABSORBABLE SUTURES - STERILE MDR Certified

**TRULON - polyamide nylon non absorbable monofilament suture**

**Equivalent Ethicon brand: ETHILON / NUROLON**

**Technical details**

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

**Applications**

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRULON - box of 12								
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length	
22800	2/0	3.0	Blue	3/8	Reverse Cutting	26 mm	45 cm	
22801	2/0	3.0	Black	3/8	Reverse Cutting	45 mm	70 cm	
22803	3/0	2.0	Blue	3/8	Conventional Cutting	20 mm	45 cm	
22804	3/0	2.0	Blue	3/8	Reverse Cutting	24 mm	45 cm	
22805	3/0	2.0	Blue	3/8	Conventional Cutting	25 mm	45 cm	
22807	3/0	2.0	Black	3/8	Reverse Cutting	19 mm	45 cm	
22808	3/0	2.0	Black	3/8	Reverse Cutting	24 mm	45 cm	
22810	4/0	1.5	Blue	3/8	Reverse Cutting	19 mm	45 cm	
22811	4/0	1.5	Black	3/8	Conventional Cutting	10 mm	38 cm	
22812	4/0	1.5	Black	3/8	Reverse Cutting	19 mm	45 cm	
22814	5/0	1.0	Blue	3/8	Reverse Cutting	16 mm	70 cm	
22815	5/0	1.0	Blue	3/8	Reverse Cutting	19 mm	70 cm	
22816	5/0	1.0	Black	3/8	Reverse Cutting	19 mm	45 cm	
22817	5/0	1.0	Black	3/8	Reverse Cutting	19 mm	75 cm	
22820	6/0	0.7	Blue	3/8	Reverse Cutting	12 mm	45 cm	
22821	6/0	0.7	Blue	3/8	Reverse Cutting	16 mm	45 cm	
22822	6/0	0.7	Black	3/8	Reverse Cutting	16 mm	45 cm	

**TRULENE - monofilament non absorbable suture - polypropylene**

**Equivalent Ethicon brand: PROLENE**

**Technical details**

- monofilament polypropylene
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

**Applications**

- hernia repair
- permanent fixtures
- general surgery
- micro surgery
- orthopedic surgery
- cardiovascular surgery
- plastic surgery
- can also be used in ophthalmic and neurological procedures

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



**TRUSILK - braided black silk non absorbable multifilament suture**

**Equivalent Ethicon brand: MERSILK / PERMA-HAND SILK**

**Technical details**

- made from 100% protein fiber spun
- coating of wax/silicone
- sterilized by ethylene oxide
- shelf life of 5 years

**Applications**

- general surgery
- ophthalmic surgery
- cardiovascular surgery
- gastrointestinal surgery
- orthopedic surgery

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUSILK - box of 12								
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length	
22830	0	3.5	Black	1/2	Round Body	30 mm	76 cm	
22831	2/0	2.0	Black	1/2	Taper Cutting	19 mm	75 cm	
22832	2/0	3.0	Black	3/8	Reverse Cutting	26 mm	45 cm	
22833	2/0	3.0	Black	1/2	Reverse Cutting	22 mm	45 cm	
22834	2/0	3.0	Black	1/2	Taper Cutting	25 mm	75 cm	
22835	2/0	3.0	Black	1/2	Round Body	30 mm	76 cm	
22836	3/0	2.0	Black	1/2	Conventional Cutting	22 mm	45 cm	
22837	3/0	2.0	Black	1/2	Round Body	22 mm	45 cm	
22838	3/0	2.0	Black	3/8	Reverse Cutting	19 mm	70 cm	
22839	3/0	2.0	Black	3/8	Reverse Cutting	24 mm	45 cm	
22840	3/0	2.0	Black	1/2	Taper Cutting	19 mm	75 cm	
22842	3/0	2.0	Black	1/2	Taper Cutting	25 mm	75 cm	
22844	4/0	1.5	Black	1/2	Reverse Cutting	17 mm	45 cm	
22845	4/0	1.5	Black	1/2	Taper Cutting	19 mm	75 cm	
22846	4/0	1.5	Black	1/2	Taper Cutting	25 mm	75 cm	
22847	4/0	1.5	Black	Straight	Round Body	19 mm	45 cm	
22848	4/0	1.5	Black	3/8	Reverse Cutting	13 mm	45 cm	
22849	4/0	1.5	Black	3/8	Reverse Cutting	19 mm	45 cm	
22851	5/0	1.0	Black	3/8	Reverse Cutting	19 mm	70 cm	
22852	6/0	0.7	Black	3/8	Reverse Cutting	16 mm	45 cm	

TRULENE - box of 12								
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length	
22870	1	4.0	Blue	1/2	Round Body	30 mm	70 cm	
22871	0	3.5	Blue	1/2	Round Body	30 mm	70 cm	
22872	2/0	3.0	Blue	3/8	Reverse Cutting	24 mm	70 cm	
22874	3/0	2.0	Blue	3/8	Reverse Cutting	24 mm	70 cm	
22875	3/0	2.0	Blue	Straight	Conventional Cutting	60 mm	70 cm	
22877	4/0	1.5	Blue	3/8	Reverse Cutting	19 mm	45 cm	
22879	4/0	1.5	Blue	3/8	Reverse Cutting	19 mm	75 cm	
22881	5/0	1.0	Blue	3/8	Reverse Cutting	16 mm	70 cm	
22882	5/0	1.0	Blue	3/8	Reverse Cutting	19 mm	70 cm	
22884	6/0	0.7	Blue	3/8	Reverse Cutting	13 mm	45 cm	
22885	6/0	0.7	Blue	3/8	Conventional Cutting	16 mm	45 cm	

# SURGICAL ABSORBABLE SUTURES - STERILE

**TRUBARB - absorbable surgical needled-suture (synthetic) knotless tissue-closure device**

**Equivalent Ethicon brand: STRATAFIX**

**Technical details**

- monofilament polydioxanone
- secure and effective triangular end-stopper
- single angle cut uni-directional helical barb design
- world class 300 series QNPL needle
- sterilized by ethylene oxide
- shelf life of 3 years

**Applications**

- gynaecological and urology procedures
- general surgery procedures
- oncology procedures
- plastic and cosmetic surgeries

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.



TRUBARB - box of 6								
GIMA code	Gauge	Metric	Colour	Circle	Needle profile	Needle length	Suture length	
22890	0	3.5	Violet	1/2	Taper point	37 mm	30 cm	
22891	0	3.5	Violet	1/2	Taper point	37 mm	45 cm	
22892	2/0	3	Violet	1/2	Taper point	37 mm	30 cm	
22893	2/0	3	Violet	1/2	Taper point	37 mm	45 cm	
22894	2/0	3	Violet	1/2	Taper point	26 mm	30 cm	
22895	3/0	2	Violet	1/2	Taper point	26 mm	45 cm	
22896	3/0	2	Violet	3/8	Reverse cutting	19 mm	45 cm	