

ETHICON SUTURES - ABSORBABLE AND NON ABSORBABLE




ETHICON

NON ABSORBABLE STERILE SURGICAL SUTURES

GIMA code **ETHICON code** **Gauge** **Metric** **Colour** **Circle** **Needle length** **Needle** **Suture length** **Minimum order**

Non absorbable surgical sutures are indicated for soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures

ETHILON - monofilament suture



22296	1669H	3/0	2.0	Black	3/8	19 mm	PS-2	45 cm	box of 36
22297	EH7145H	4/0	1.5	Black	3/8	19 mm	FS-2	45 cm	box of 36
22298	1162H	4/0	1.5	Black	3/8	16 mm	FS-3	45 cm	box of 36
22299	2808G	8/0	0.4	Black	3/8	6.5 mm	BV130-5	13 cm	box of 12
22300	660H	6/0	0.7	Black	3/8	16 mm	FS-3	45 cm	box of 36
22301	661H	5/0	1.0	Black	3/8	19 mm	FS-2	45 cm	box of 36
22302	662SLH	4/0	1.5	Black	3/8	19 mm	FS-2S	45 cm	box of 36
22303	663H	3/0	2.0	Black	3/8	24 mm	FS-1	45 cm	box of 36
22304	EH7506H	3/0	2.0	Black	3/8	19 mm	FS-2	45 cm	box of 36
22305	FH5171H	5/0	1.0	Black	3/8	19 mm	FS-2	75 cm	box of 36
22306	668H	5/0	1.0	Black	3/8	13 mm	C-2	45 cm	box of 36
22307	698H	5/0	1.0	Black	3/8	13 mm	P-3	45 cm	box of 36
22308	669H	3/0	2.0	Black	3/8	24 mm	FS-1	75 cm	box of 36

PERMA-HAND - silk sutures - braided



22310	18501G	4/0	1.5	Black	1/2	17 mm	J-1	45 cm	box of 12
22311	632G	3/0	2.0	Black	1/2	22 mm	X-1	45 cm	box of 12
22312	737G	2/0	3.0	Black	1/2	22 mm	X-1	45 cm	box of 12
22313	18507G	4/0	1.5	Black	Straight	19 mm	ST-4	45 cm	box of 12
22314	641G	4/0	1.5	Black	3/8	13 mm	P-3	45 cm	box of 12
22315	18509G	3/0	2.0	Black	1/2	22 mm	SH-1 Plus	45 cm	box of 12
22316	EH7149G	4/0	1.5	Black	3/8	19 mm	FS-2	45 cm	box of 12
22317	K872	3/0	2.0	Black	1/2	17 mm	RB-1	75 cm	box of 36
22318	683H	4/0	1.5	Black	3/8	19 mm	FS-2S	45 cm	box of 36
22319	684H	3/0	2.0	Black	3/8	24 mm	FS-1	45 cm	box of 36

PROLENE - blue pigmented monofilament sutures



22320	EH7150H	4/0	1.5	Blue	3/8	19 mm	FS-2	75 cm	box of 36
22321	EH7151H	4/0	1.5	Blue	3/8	19 mm	FS-2	45 cm	box of 36
22322	W8873T	6/0	0.7	Blue	3/8	16 mm	PC3 Prime	45 cm	box of 24
22323	8661H	5/0	1.0	Blue	3/8	19 mm	FS-2	45 cm	box of 36
22324	8665H	3/0	2.0	Blue	3/8	19 mm	FS-3	45 cm	box of 36
22325	8682H	4/0	1.5	Blue	3/8	19 mm	PS-2	45 cm	box of 36
22326	8683H	4/0	1.5	Blue	3/8	19 mm	FS-2	45 cm	box of 36
22327	8687H	3/0	2.0	Blue	3/8	19 mm	PS-2	45 cm	box of 36
22328	8698H	5/0	1.0	Blue	3/8	13 mm	P-3	45 cm	box of 36
22329	8699H	4/0	1.5	Blue	3/8	13 mm	P-3	45 cm	box of 36

ABSORBABLE STERILE SURGICAL SUTURES

GIMA code **ETHICON code** **Gauge** **Metric** **Colour** **Circle** **Needle length** **Needle** **Suture length** **Minimum order**

MONOCRYL - complete absorption in 91-119 days



Monocryl is indicated for general soft tissue approximation and/or ligation, but not for cardio, neuro, micro and ophthalmic surgery

22350	Y493H	5/0	1.0	Undyed	3/8	13 mm	P-3 prime	45 cm	box of 36
22351	MCP3205G	4/0	1.5	Undyed	3/8	16 mm	PS-3	70 cm	box of 12


VICRYL PLUS BRAIDED - complete absorption in 56-70 days



Vicryl Plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular

22360	VCP392ZH	4/0	1.5	Violet	3/8	19 mm	FS-2	45 cm	box of 36
22361	VCP305H	3/0	2.0	Violet	1/2	17 mm	RB-1 Plus	70 cm	box of 36
22362	VCP452H	3/0	2.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22363	VCP453H	2/0	3.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22364	VCP311H	3/0	2.0	Violet	3/8	22 mm	SH-1	70 cm	box of 36
22365	VCP317H	2/0	3.0	Violet	1/2	26 mm	SH	70 cm	box of 36
22366	VCP494H	4/0	1.5	Undyed	3/8	13 mm	P-3	45 cm	box of 36


VICRYL RAPID BRAIDED - complete absorption in 42 days



Vicryl Rapid Braided indicated only for superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.


22370	V76H	3/0	2.0	Undyed	1/2	17 mm	J-1	45 cm	box of 36
22371	V2920G	4/0	1.5	Undyed	3/8	19 mm	FS-2	45 cm	box of 12
22372	V2930G	3/0	2.0	Undyed	3/8	19 mm	FS-2	45 cm	box of 12
22373	V2150G	3/0	2.0	Undyed	1/2	17 mm	RB-1 Plus	70 cm	box of 12
22374	V2190G	3/0	2.0	Undyed	1/2	22 mm	SH-1 Plus	70 cm	box of 12
22375	VR2298	3/0	2.0	Undyed	3/8	19 mm	FS-2	75 cm	box of 36
22376	V2910H	5/0	1.0	Undyed	3/8C	19 mm	FS-2	45 cm	box of 36

VICRYL



22380	V453H	2/0	3.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22381	V452H	3/0	2.0	Violet	3/8	24 mm	FS-1	70 cm	box of 36
22382	V304H	4/0	1.5	Violet	1/2	17 mm	RB-1	70 cm	box of 36
22383	V305H	3/0	2.0	Violet	1/2	17 mm	RB-1	70 cm	box of 36
22384	V311H	3/0	2.0	Violet	1/2	22 mm	SH-1	70 cm	box of 36
22385	18531G	4/0	1.5	Violet	Compaure	25 mm	Da-visi black	45 cm	box of 12
22386	V994H	4/0	1.5	Violet	1/2C	17 mm	V5	70 cm	box of 36
22387	V422H	3/0	2.0	Undyed	3/8C	24 mm	FS-1	70 cm	box of 36

PDS II



22390	Z340H	0	3.5	Violet	1/2	36 mm	CT-1	70 cm	box of 36
22391	Z317H	2/0	3.0	Violet	1/2	26 mm	SH	70 cm	box of 36
22392	Z292ZE	4/0	1.5	Undyed	3/8	19 mm	FS-2	45 cm	box of 24
22393	Z494E	4/0	1.5	Undyed	3/8	13 mm	P-3	45 cm	box of 24
22394	Z496E	4/0	1.5	Undyed	3/8	19 mm	PS-2	45 cm	box of 24