


DISPOSABLE SKIN STAPLER



Multilingual manual and box:
GB, FR, IT, ES, PT, DE.



• **25891 DISPOSABLE STERILE SKIN STAPLER** - box of 5
Easy to use device with 35 staples made of High Quality 317 L Stainless Steel for medical use.

High sewing up and stapling speed help shorten the time of sewing and avoid cross infection of contagious diseases. Reduced tissue reaction, thinner scar after operation and faster wound healing. Painless and convenient staples removal after the incision has healed up.

For skin closure in operating room, delivery, emergency, outpatient surgery, clinics and physician's office.

3M™ PRECISE™ SKIN STAPLERS


- **25896 3M PRECISE™ SKIN STAPLER** 5 staples - sterile - box of 12
- **25897 3M PRECISE™ SKIN STAPLER** 15 staples - sterile - box of 12

The disposable Precise skin staplers offer numerous benefits over sutures


- significantly quicker, simpler to use and safer than sutures
- better cosmetic effect
- six times faster than suturing
- no risk of accidental stick injury
- less painful than sutures
- easy to remove.

It comes pre-loaded with 5 or 15 staples (6.5x4.5 mm).

Suitable for simple, recent, open, unbruised wounds. Staples are particularly suitable for scalp wounds as they do not catch on surrounding hair and are more easily identifiable at the time of removal.



SKIN STAPLE REMOVERS



25893

• **25893 DISPOSABLE SKIN STAPLE REMOVER** 11 cm - sterile - box of 50
Suitable for use with most brands of skin staples. It combines functionality with cost effectiveness.

26722

26737



- **26737 S/S SPENCER SCISSORS** - 13 cm
- **26722 S/S MICHEL STITCH REMOVER FORCEPS** - 12 cm



25892
CE 0123

• **25892 DISPOSABLE REMOVER FORCEPS** - sterile - box of 20 pcs.
Practical disposable forceps for removing regular and wider size staples. Plastic handle and stainless steel tips.



25895

• **25895 SR-3 DISPOSABLE SKIN STAPLE REMOVER** - sterile - box of 10 pcs.
The scissor-like action of the 3M SR-3 provides ergonomic grip and removes staples in the same direction they were implanted, making removal simple. Plastic handle and stainless steel tips.

BONEWAX - STERILE

• **23000 TRUWAX STERILE SURGICAL BONEWAX** 2.5 g - box of 12

Equivalent Ethicon brand: **BONE WAX**

Truwax is a sterile bone wax used to control bleeding from bone surfaces during surgical procedures by providing hemostasis from bone injuries and cuts. It works like a mechanical barrier, through a tamponade effect, it blocks pores and channels of cut or damaged bone. Made from sterile mixture of beeswax paraffin wax and isopropyl palmitate. Materials ensure excellent malleability. Multilingual manual and box:

GB, FR, IT, ES, PT, DE, GR.

MDR Certified



23000

TECHNICAL SPECIFICATIONS

Composition:	Bees wax I.P., White Hard Paraffin Wax, I.P., Iso Propyl Palmitate, U.S.P.
Colour - Odour:	opaque - waxy odour
Sterilization:	gamma radiation
Shelf life:	5 years
Packing:	2.5 g per unit pack



CLINICAL OXIDIZED REGENERATED CELLULOSE

CLINICAL OXIDIZED REGENERATED CELLULOSE

Equivalent Ethicon brand: **SURGICEL**

Features

- absorbable hemostat
- made from oxidized regenerated cellulose (plant based product)
- malleable and flexible

Applications

- Neurosurgery: operative intervention into the cerebrum, ventricular system, spinal surgery
- Cardiac Surgery: valve replacement surgery, sterna medullary bleeding, approaches to the skull base, cerebrovascular surgery
- Gynecological surgeries: uterine wall perforation
- ENT: posterior epistaxis
- General surgeries: iatrogenic perforations

TECHNICAL SPECIFICATIONS

Effective Hemostasis:	in 3-5 minutes
Sterilization:	gamma radiation
Shelf life:	3 years



23005



23008

Multilingual manual and box:
GB, FR, IT, ES, PT, DE, GR.

MDR Certified

GIMA code	CLINICEL STANDARD - box of 6	cm	inches
23003	Clinicel standard	5.1x7.6	2x3
23004	Clinicel standard	5.1x35.5	2x14
23005	Clinicel standard	10x20	4x8
GIMA code	CLINICEL FIBRIL TYPE - box of 6	cm	inches
23007	Clinicel fibril type	2.5x5.1	1x2
23008	Clinicel fibril type	5.1x10	2x4