

# SOEHNLE PROFESSIONAL DIGITAL SCALES

- 27275 SOEHNLE 7830
- 27276 SOEHNLE 7831 with height rod and BMI display
- 27277 SOEHNLE 7730  $\text{C} \text{E} \text{M}$  Class III

The new range of Soehnle professional personal scales with stand is a special highlight with its gently curved glass indicator. While stepping on the scale, the display is activated automatically and the weight is shown with red backlighting numbers. It does not require any operation and is perfectly suitable for simple weighing in medical institutions, wellness areas and pharmacies.

**Scale with integrated height rod Soehnle 7831 (27276) with stand and height rod switches automatically from kg to cm if the height rod is raised.** Thus, patients can be simultaneously weighed and measured.

Features:

- indication of weight, height and BMI (27276 only)
- integrated height rod from 120 cm to 202 cm (27276 only)
- simple to use, tare key on the back of the display
- glass scale indicator (glass thickness 8 mm)
- bright LED display: 5 digits, height 15 mm

Made in Germany.

Manual in: DE, GB, FR, ES, PT (only 27277), IT, NL, DK, FI, SE, NO, GR, TR, CZ, SK, PL, RU.

### TECHNICAL SPECIFICATIONS

	27275	27276	27277
Capacity:	250 kg	250 kg	250 kg
Sensitivity:	100 g	100 g	50/100 g
Dimension mm:	365x490x890	365x490x1,280	365x490x890
Weight:	14.1 kg	17.2 kg	14.1 kg
Power Supply:	230 V - 50-60 Hz		
	(available on request 110 V)		

No keys on the glass plate



# SECA MEDICAL SCALES - CLASS III




- 27284 SECA 756 - mechanical with height rod (60-200 cm) - class III
- Tradition combined with modern design. With a very sensitive precision mechanism to show weight through the aid of a well visible red pointer on the circular dial.
- User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

### TECHNICAL SPECIFICATIONS

Capacity: 160 kg  
Sensitivity: 500 g  
Weight: 18.5 kg  
Dimensions: 628x360x1,000 mm  
Platform: 360x350 mm with wheels

Keys to symbols:  $\text{M}$  = Metrological Approval for use in medical field

$\text{C} \text{E} \text{C}$  = Calibration according to class III sensitivity required by 90/384/EEC



- 27296 SECA 711 - mechanical with height rod (60-200 cm) - class III
  - 27297 SECA 700 - mechanical with height rod (60-200 cm)
- Mechanical column scale with weight adjustment at eye level. Weight can be measured in a conventional way: by aligning two chrome plated indicators.
- User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.
- 27298 SECA 220 HEIGHT ROD - spare

### TECHNICAL SPECIFICATIONS

Capacity: 220 kg  
Sensitivity: 100 g (27296) - 50 g (27297)  
Weight: 21.7 kg  
Dimensions: 520x1,556x520 mm  
Platform: 350x360 mm



- 27292 SECA 799 - electronic with BMI - class III
  - 27298 SECA 220 HEIGHT ROD
- Versatile, low maintenance and site independent with economic battery operation, ideal for daily use in hospital and doctors practices. Column scale is easy to move around with its transport castors. BMI function permits a reliable evaluation of the nutritional conditions of the patient. Equipped with several functions: TARE, HOLD, weighing range switch-over, automatic switch off. User manual in GB, FR, IT, ES, DE, PT, PL, NL, GR, DK, SE, FI, NO.

### TECHNICAL SPECIFICATIONS

Capacity: 200 kg  
Sensitivity: 100 g < 150 kg > 200 g  
Weight: 6.3 kg  
Dimensions: 294x831x417 mm  
Platform: 272x280 mm

$\text{C} \text{E} \text{C}$  0109 = Approved scales according to 93/42 EEC