

GIMA SHARPS CONTAINERS

Full range of Sharps containers for the disposal of sharp or cutting waste. Equipped with temporary and definitive closure, certified hermetic seal, do not release any harmful emission during incineration.

The containers hold the UN mark for the transportation of dangerous materials.

- practical opening system
- temporary and permanent closure
- needles removal
- insulin pen removal device
- maximum filling level
- shock resistant
- puncture resistant
- easy to use

- easy to assemble
- excellent stability
- square or round shaped
- sturdy handle for easier transport
- customisable label for quantities
- autoclavable at 134°C for 18'
- no emission of any harmful gas during the incineration process



| GIMA code | CS LINE SHARPS CONTAINERS | Minimum order | GIMA code | PBS SHARPS CONTAINERS | Minimum order |
|--------------------------------|---|---------------|---------------------|------------------------------------|---------------|
| 25782 | CS 2 l - 14.8x14.8xh 19.5 cm | box of 60 | 25790 | PBS 0.6 l - Ø 10xh 14.5 cm | box of 100 |
| 25783 | CS 3 l - 16.7x16.7xh 23.2 cm | box of 50 | 25791 | PBS 1.5 l - Ø 14xh 16.5 cm | box of 100 |
| 25784 | CS 4 l - 16.7x16.7xh 28.2 cm | box of 40 | 25792 | PBS 2 l - Ø 14xh 21 cm | box of 75 |
| 25785 | CS 5 l - 21x21xh 20 cm | box of 30 | 25793 | PBS 5 l - Ø 21.4xh 19.5 cm | box of 55 |
| CS PLUS LINE SHARPS CONTAINERS | | | 25796 | PBS 6 l - Ø 21.4xh 24 cm | box of 55 |
| 25787 | CS Plus 10 l - 22x22xh 32.5 cm | box of 33 | 25794 | PBS 7 l - Ø 25xh 19.5 cm | box of 40 |
| 25788 | CS Plus 22 l - 37.5x27.2xh 34.5 cm | box of 10 | 25795 | PBS 12 l - Ø 25xh 31.3 cm | box of 25 |
| 25799 | Fixing system* - horizontal or vertical | | TRAY FOR CONTAINERS | | |
| *Suitable for wall or trolley | | | 37725 | Compartment tray* 32x24.5xh 4.4 cm | box of 24 |

CERTIFICATIONS



FDA K071517

ISO 23907-1:2019

*Suitable for 25790-2, 25782

37725 + 25790

WASTE CONTAINERS - 30-50-60 L

WASTE CONTAINERS - 3 sizes

Containers for disposal of special hospital waste in 3 sizes, with 2 different lids.

- extremely solid and shock resistant at high and low temperature
- hermetic seal by peripheral adhesive on the lid (single or double)
- lateral and central handles tested for twice the max weight
- stacking up to a maximum of 3 m with full load
- conformity: ISO 23907-1:2019

MADE IN ITALY



| GIMA code | SC LINE WASTE CONTAINERS | | | |
|-----------|--------------------------|--------|-------------|---------------|
| | Tide | Lid | Max load kg | Size |
| 25774 | Container 30 l | single | 15 | 40x30xh 39 cm |
| 25775 | Container 30 l | double | 15 | 40x30xh 39 cm |
| 25776 | Container 50 l | single | 20 | 40x30xh 56 cm |
| 25777 | Container 50 l | double | 20 | 40x30xh 56 cm |
| 25778 | Container 60 l | single | 25 | 40x30xh 65 cm |
| 25779 | Container 60 l | double | 25 | 40x30xh 65 cm |

Single lid with adhesive


Double lid with adhesive

CONTAINERS UP TO 72 CM

LONG SHARPS CONTAINERS

Disposable containers for the elimination of long objects up to 72 cm, such as biopsy forceps and sharps. Large opening, for elimination of bulky objects, temporary and definitive closure, ergonomic handle. Autoclavable at 134°C for 18'.

MADE IN ITALY



| GIMA code | LONG SHARPS CONTAINERS |
|-----------|---|
| 25797 | Magnum L 6.3 l - 17.2x12xh 61.7 cm |
| 25798 | Magnum XL 8.3 l - 17.2x12xh 84.5 cm |
| 25799 | Fixing system* - horizontal or vertical |

*Suitable for wall or trolley