


TOURNIQUETS IN LATEX, RUBBER AND LATEX-FREE



GIMA TOURNIQUETS
 Efficient and simple solutions for professional medical use.
Latex: Ø 5 mm, for manual use
Velcro: for easy fastening
Velcro "D" buckle: self applicable made of a nylon band threads through plastic "D" buckle
Lock system: self applicable with a quick release system.
 Multilingual box (excluding 25733): GB, FR, IT, DE, ES, PT, GR, Arabic.




GIMA code	TOURNIQUETS	Material	Colour	Length cm	Width cm	Fastening systems	Origin
25733	TUBOLARE Tourniquet - box of 200	Latex	Yellow	33	-	Manual	Europe
25721	GIMA Tourniquet	Rubber band, terylene cover	Yellow	40	2.5	Velcro	Asia
25725	STRIP-VELCRO Tourniquet	Velcro	Green	45	2.5	Velcro+ring: use with one hand	Asia
25727	FAST Tourniquet	Rubber band, terylene cover	Red	40.5	3	Single lock system: quick release	Asia
25728	FAST Tourniquet	Rubber band, terylene cover	Green	40.5	3	Single lock system: quick release	Asia
25726	FAST Tourniquet - latex free	Spandex band, terylene cover	Blue	40.5	3	Single lock system: quick release	Asia
37500	FAST Tourniquet - daisy	Rubber band, terylene cover	Blue floral	40.5	3	Single lock system: quick release	Asia
37501	FAST Tourniquet - peace	Rubber band, terylene cover	Multi-colored	40.5	3	Single lock system: quick release	Asia
37502	FAST Tourniquet - space	Rubber band, terylene cover	Blue space	40.5	3	Single lock system: quick release	Asia

JETPULL 2 TOURNIQUETS WITH GRADUAL ADJUSTMENT



JETPULL 2 WITH GRADUAL ADJUSTMENT
 Easy to use, allows single-hand operation. Sophisticated 2 button system allows an easy opening of the fastening with 2 choices: smooth or fast, even when under tension. Fastening device is made using high class medical plastics (PVC free). Ribbon is in special material to avoid any skin bruises. Washable at 60/90° C and autoclavable.



GIMA code	TOURNIQUETS	Material	Colour	Length cm	Width cm	Fastening systems	Origin
25722	JETPULL 2	Polyester/rubber	Blue	56	2.4	Double lock system: quick and slow release	Germany
25723	JETPULL 2 - latex free	Polyester/elastan	Green	56	2.4	Double lock system: quick and slow release	Germany
25724	JETPULL BABY	Polyester/rubber	Child design	56	2.4	Double lock system: quick and slow release	Germany

TAPE TOURNIQUET

• 37515 TAPE TOURNIQUET - latex free - disposable - box of 25

Latex free disposable tape tourniquet made of TPE with dot surface. Each box contains 25 pre-cut tourniquets 45x2.5 cm (18"x1") with 0.636 mm thickness. Multilingual box: GB, FR, IT, ES, PT, DE, NL, PL, RO, GR, EE, SE, HU, Arabic.



PRE-CUT TOURNIQUETS

• 37510 PRE-CUT TOURNIQUET 45x2.5 cm - blue - box of 250
 • 37511 PRE-CUT TOURNIQUET 45x2.5 cm - orange - box of 250
 • 37512 PRE-CUT TOURNIQUET 45x2.5 cm - violet - box of 250

Latex free disposable tape tourniquet made of TPE with smooth surface. Thickness 0.636 mm
 Multilingual box: GB, FR, IT, ES, PT, DE, PL, NL, RO, GR, EE, SE, HU, Arabic.

