

## FLORIDA LEATHER/PU BAGS

External dimensions:  
38x20xh 20 cm



- 27121 FLORIDA leather - burgundy
- 27122 FLORIDA leather - black
- 27123 FLORIDA leather - brown
- 27124 FLORIDA leather - chestnut



Rigid bottom, inside zip fastener central pocket, lateral inside pocket, burnished lock. Outside made of leather, inside made of PU (polyurethane). Made in Asia - Italian design.

## KANSAS BAGS - 5 COLOURS

- 27095 KANSAS skay - black
- 27096 KANSAS skay - chestnut
- 27097 KANSAS skay - blue
- 27098 KANSAS skay - ochre
- 27099 KANSAS skay - brown

Made in Asia - Italian design

Rigid bottom, outer velcro pocket, inside partition, inside zip pocket, lock with key.



27099

27097

27096

27098

External dimensions:  
42x20xh 24 cm

External pocket with velcro closure

27097

## HEAT REGULATING BAG

- 27218 ENERGY HEAT REGULATING BAG

Heat regulating bag for the warming of IV liquids. Used as an active source of heat is indicated for the prevention and treatment of hypothermia. Ideal for emergency medicine, hospitals and field camps. Its compact design allows for use both hung on the stretcher or in the ambulance as a bottle holder or inside bigger containers. The infusion liquids remain at a constant temperature and the IV liquids can be transported out of the ambulance for a short period. The main compartment is heated by an electrically fed resistive pad. Supplied with front pocket for the electric power supply cable, transparent identity pocket, cable with plug, adjustable shoulder strap and manual GB, IT.



27218

SPENCER  
MADE IN ITALY

### TECHNICAL SPECIFICATIONS

Padding material: outside PVC, inside nylon 210  
Voltage supplies: 12 V DC, absorption 3,75 A  
Car lighter: DIN ISO 4165 Protection: IP3x  
Size: 24x14xh 33 cm Weight: 1,2 kg  
Isolation: Class II

## VACCINE CARRIER

- 27211 PROFESSIONAL VACCINE CARRIER

Compact, extra strong with special insulation guaranteed by high density PU foam, suitable for transportation of vaccines. Comes complete with 8 empty ice-packs 400 ml suitable for a cold life of 40 hours at 8°C with external temperature of 43°C without opening. Vaccine storage capacity: 2.6 l. Body and lid with gaskets for better insulation.



27211

External size: 28.5x25.5xh 32 cm.  
Internal size: 22x18xh 24 cm.  
MADE IN ITALY

## THERMAL BAGS

- 27209 THERMAL BAG - yellow - nylon
  - 27212 THERMAL BAG - blue - nylon
- Crafted, durable and insulated thermal bag with detachable shoulder strap. Supplied with pockets and 3 temperature packs, which maintain cold and hot temperature for long period. External size: 34x25x24 cm - Internal size: 30x23x21 cm
- 27206 THERMAL BAG TEMPERATURE PACK - pack of 2
  - 27216 DUAL SOFT ICE PACK 100+100 g
  - 27215 SOFT ICE PACK 200g

MADE IN ITALY



27209

27206



27215



27216

