

## LOGIC-1 RUCKSACK

27170



- 27170 LOGIC-1 RUCKSACK - polyester - medium - 40x20xh 47 cm
- Supplied with 5 coloured pouches with transparent window and velcro strips, 2 frontal pockets and 2 lateral pockets. Internal velcro strips and elastics, to retain instruments. Provided with yellow reflective strips.



27170

## LOGIC-2 RUCKSACKS

POLYESTER  
PVC COATED

27171



PVC TARPAULIN

27172



- 27171 LOGIC-2 RUCKSACK - polyester - large - 40x25xh 47 cm
  - 27172 LOGIC-2 RUCKSACK - PVC - large - 40x25xh 47 cm
- Rucksack externally made of strong water resistant 600D polyester PVC coated (27171) or PVC tarpaulin (27172). Supplied with 6 coloured pouches with transparent window and velcro strips, 2 frontal pockets and 2 lateral pockets. Internal velcro strips and elastics, to retain instruments. Central division with zip and elastics. Provided with yellow reflective strips and shoulder straps.



27171



27172



## AMBULANCE RUCKSACK

27174



Large 10 mm zips

27174

- 27174 AMBULANCE RUCKSACK - polyester
- Modular extra large rucksack provided with 8 internal pouches, 4 coloured with handle (model E5) and 4 coloured with transparent window (model E2). Velcro internal strips for pouches fixing. Provided with double yellow reflective strip. Size: 40x30xh 57 cm



Rubber bottom



27174

## SILOS-2 RUCKSACKS

POLYESTER

27168



Velcro removable plate

POLYESTER WITH  
OUTER PVC COATING

27169



- 27168 SILOS-2 RUCKSACK - polyester - blue - 38x24xh 50 cm
  - 27169 SILOS-2 RUCKSACK - PVC - red - 38x24xh 50 cm
- Professional rucksacks made of strong water resistant 600D polyester PVC coated and with outer PVC (27169 only)
- removable compartment with 4 transparent pockets and many elastic strips
  - central storage space with 3 pockets, velcro strips and elastic strips
  - two external pockets with elastic strips and velcro strips: one 35x16 cm, one 35x32 cm with 4 internal transparent pockets
  - rubber bottom with 4 feet
  - two removable external plates for printing name
  - comfortable back shoulder straps and plastic handle on both sides.



27169