

BMI TAPES



- **27341 BMI - Tape Measure**
Pocket retractable tape measure 150 cm.
BMI (Body Mass Index) calculator on front panel.
- **27342 BMI - Tape Measure**
Economical version of 27341 with same features.

BMI CALCULATOR

Can be personalized with your logo for a minimum of 1,000 pieces.



- **27336 BMI CALCULATOR - box of 10**
Used in hospital and surgeries gives an indication of general health through the ratio of weight to height.
Made of 2 rotating discs: 108 mm diameter (base) and 76 mm (top).

CIRCUMFERENCE MEASURING TAPE

- **27345 SECA 203 CIRCUMFERENCE MEASURING TAPE with extra Waist-To-Hip-Ratio (WHR) calculator**
Retractable and ergonomic measuring tape for accurate determination of body circumference. Measuring range: 15-205 cm, division 1 mm. High-quality metal scale to obtain the Waist-To-Hip-Ratio (WHR). The WHR is an indirect method of determining the amount of abdominal fat tissue. The ratio allows conclusions about the distribution of fat stores and thus contributes to a diagnosis. Medical device, Class IM. Multilingual manual: GB, FR, IT, ES, DE, PT, PL.




BODY TAPE

- **27343 BODY TAPE MEASURE**
Sturdy vinyl tape to track body measurements and progress. Measures in inches and centimeters any body part: arm, thigh, calf, chest, waist, hip... Total measurement length: 60 inches / 152 cm. Push-button retraction and locking



ensure easy measurement and superior accuracy.


Multilingual box and manual: GB, FR, IT, ES, PT, GR.



PROFESSIONAL BODY - HEAD TAPE



- **27347 PROFESSIONAL BODY TAPE**
- **27350 PROFESSIONAL BODY TAPE with pediatric gauge box**
Professional body tape measure for adult and pediatric use. Tape is coloured on both sides to help to identify metabolic risk in abdominal measure for male and female.
Length of tape: 8-150 cm.
Thanks to its rounded shape, it can be used also as a head measure tape. Packed in single box.
Body tape with pediatric gauge box (27350)
Specific tape for pediatric use. Supplied in a box with a ruler to compare the measured head circumference with the standard percentiles of children's growth from 0 to 36 months.
Suitable for both sexes.



Multilingual manual: GB, FR, IT, DE, ES, GR, PT, Arabic.

AISI 430 STAINLESS STEEL GONIOMETERS


All goniometers have manual in GB, FR, IT, ES




- **24900 FINGER GONIOMETER 6" - 15 cm**
Finger goniometer measures range of motion of finger joints. Linear scale in inches and centimeters on arms. 180° scales marked in 5° increments.



- **24901 180° DIGIT GONIOMETER 6" - 15 cm**
Allows taking readings using only one hand. Measures 110° flexion through 40° hyperextension in 5° increments.



- **24903 360° GONIOMETER 14" - 35 cm**
This goniometer has two opposing 180° scale and one 360° scale all marked in 1° increments. Sturdy adjustable knob, controlling tension of arms.



- **24902 180° ROBINSON GONIOMETER 7" - 17.5 cm**
Pocket size 180° Robinson goniometer for doctor and therapist. Linear scale marks 5° increments.



- **24904 180° EXTRA LONG ARMS GONIOMETER 14" - 35 cm**
Extra long arms, useful for measuring all joints. Arm has a linear scale marked in 1° increments.