

MOBILE STADIOMETER

• 27317 SECA 213 MOBILE STADIOMETER

Mobile stadiometer for measuring height, simple and easy to transport with integrated handle. The measuring rod can be dismantled into several pieces and set up easily and quickly. The spacer keeps the rod straight and stable without any fittings. Flexible position due to large floor plate which ensures stability. It guarantees precise and clearly visible reading of results. Suitable for adults and children and for measuring patients lying down.

Multilingual manual in:
GB, FR, IT, ES, PT, DE,
PL, NL, GR, DK, SE,
NO, FI.

Integrated handle and
secure catches on all
individual parts for
transport ease.

TECHNICAL SPECIFICATIONS

Measuring range: 20-205 cm
Graduation: 1 mm
Size (open): 337x590xh 2,130 mm
Size (folded): 337x177x624 mm
Weight: 2.4 kg



Large platform
for stability.

27317

SECA INFANTOMETER



The unique V-shape automatically places
the baby in the right measuring position

seca



Size: 110x16.9x40.2 cm.

Weight: 3.8 kg.



The head and foot
positioners lock
into place to capture
results. Indication
can be read after
measuring

• 27316 SECA 416 - INFANTOMETER

Infantometer for measuring babies up to two years old. Sturdy construction and stability makes it quick and easy to use. Rounded sides guide the baby easily and safely into the right position. Permanently mounted with rollers on two guide rails, the footpiece slides smoothly along the scale. Footpiece is lockable so the baby can be removed and the result noted later. Measuring range: 33-100 cm. Graduation: 1 mm.

BABY MEASURE


• 27326 BABY MEASURE "CALIBRO"

Made of oxidized aluminium, metric graduation 0-900 mm, division 5 mm. Certified by Italian "weights and measure institute". Weight: 300 g.



27326

ELECTRONIC HEIGHT ROD - ADULT AND PEDIATRIC



27333

SOEHNLE
PROFESSIONAL

LCD digit
height: 9 mm

Size:
345x92x28 mm

Power supply:
2x1.5 V AA batteries

Weight:
300 g

• 27333 SOEHNLE ELECTRONIC HEIGHT ROD

Measuring range: 50-240 cm (20-94.5 inches), division: 1 cm. Manual in: DE, GB, FR, ES, IT, GR, CZ.

Switchable
to inches

Includes adapter for attachment to
baby scales 27306, 27260

Power supply:
2x1.5V AAA batteries

SOEHNLE

Size:
630x290x70 mm

• 27264 SOEHNLE BABY HEIGHT ROD

Digital and mechanical rod 35 to 80 cm, division 1 cm. LCD display, digit height 16 mm. Manual in 23 languages.

27264

HEIGHT MEASURE

• 27335 TAPE HEIGHT MEASURE

Clear scale 0-200 cm.
Easy wall fixation.
Supplied with screw for wall mounting.
Multilingual box and manual:
GB, FR, IT, ES, PT, DE, PL, RO, HU, FI, GR,
Arabic.



27335



CRANIOMETRO - HEAD MEASURE

• 27327 CRANIOMETRO - HEAD MEASURE

Made of special, soft, non stretchable plastic. Measures child head girth up to 60 cm. Instructions: GB, FR, IT, ES.



27327

MADE IN ITALY

BONE MEASURE

• 27325 BONE MEASURE


Slide-knot caliper in anticorodal to measure knee, shoulder, wrist, elbow... Measuring capacity: 150 mm. Burner length: 100 mm. Tested by Italian "weights and measures institute".

MADE IN ITALY



27325

ORCHIDOMETER



27337

A string of 12 numbered
wooden beads of increasing size
from about 1 to 25 mm used to
determine size of testicles.

• 27337 ORCHIDOMETER - wood

Used by pediatricians for screening of children's testicles. Multilingual box and instructions: GB, FR, IT, ES, PT, DE, GR, NL, SE, PL, FI, CZ.