

HOSPITAL DIGITAL BABY SCALE



- **27267 SECA 376 HOSPITAL SCALE - digital class III**
- high-quality Seca scale technology with fine 5-gram graduation
- TARE, HOLD, auto-HOLD and BMIF for extra convenience
- extra large and ergonomic tray for safe weighing
- adjustable damping
- large display for easy reading

TECHNICAL SPECIFICATIONS

Capacity: 20 kg
Sensitivity: 5 g up to 7.5 kg, 10 g from 7.5 kg up
Dimensions - weight: 620x190x358 mm - 3.7 kg
Platform: 615 x 130 x 275 mm
Battery operated

SECA ELECTRONIC BABY SCALES



- **27282 SECA 354 - electronic**
 - **27283 SECA 384 - electronic - CE M Class III**
- Two scales in one: with a comfortable tray to weigh in a reclining position, easily detachable to be converted into a floor scale.
Breast-milk-intake function determines how much milk the baby has suckled.
Other functions: tare, auto-hold, kg/lbs switch over automatic switch off.
User manual: GB, FR, IT, ES, DE (27282-3), PT, PL, NL, GR, DK, SE, FI, NO (27283 only).

TECHNICAL SPECIFICATIONS

Capacity: 20 kg
Dimensions: 552x156x332 mm
Batteries operated
Sensitivity: 10 g < 10 kg > 20g
Weight: 2.3 kg (27282), 2.8 kg (27283)
Made in Asia

GIMA MECHANICAL BABY SCALE




- **27300 GIMA MECHANICAL BABY SCALE**
- High quality sturdy mechanical baby scale for accurate and sensitive weighing.
Supplied with curved weighing tray, aluminium levers and pins, stable casing for safe weighing.
User manual in GB, FR, IT, ES, PT, DE, RO, GR.

TECHNICAL SPECIFICATIONS

Capacity: 20 kg
Dimensions: 570x330x285 mm
Sensitivity: 10 g
Weight: 6.5 kg

SECA MECHANICAL BABY SCALE




- **27299 SECA 725 - mechanical**
 - **27294 SECA 745 - mechanical - CE M Class III**
- The sturdy classic: with its softly curved weighing tray, the stable all-steel casing and the precise, smooth running weights, this model is the indestructible original baby scale.
User manual: GB, FR, IT, ES, DE, PT, PL.


TECHNICAL SPECIFICATIONS

Capacity: 16 kg
Dimensions: 550x180x290 mm
Made in Asia
Sensitivity: 5 g (27299) 10 g (27294)
Weight: 6.4 kg

MEASURING MAT FOR BABIES



- **27331 GIMA BABY MEASURING MAT**
- Washable and foldable measuring mat for babies and toddlers, easy to carry. For easy and precise measurement of baby length. Size: 121x30 cm. Measuring range: 10-99 cm. Graduation: 5 mm.



- **27329 SECA 210 - BABY MEASURING MAT**
- Mobile measuring mat for babies and toddlers. Kind to the skin, washable, foldable.
The space saving solution for easy and precise measurement of baby length. Range 10-99 cm. Graduation: 5 mm.