

FITNESS SCALE


• 27286 SOEHNLE 7850 FITNESS SCALE

Fitness scale for measurement of muscle, fat and water components, BMI and caloric needs per day - impedance measurement through hands, no need to take off the shoes

- simple self-explanatory menu guide on the glass surface
- no limit on measurement for over- and underweight people
- clothing weight value seasonally adjustable
- print-out of data/time, weight, data analysis and recommendations for nutrition and fitness activities
- software to edit the print-out and settings in DE, GB, FR, NL, ES, IT, HU, CZ, FI, SE, NO, DK, GR, PL, SI, RU, TR can be downloaded from www.soehnle.com

- kg/lbs/st selection function
Manual in: GB, FR, IT, DE, CZ.

- 35674 PAPER ROLL 57mm x 25m
- box of 20



Extra large weighing platform

Glass surface easy to clean

Large Backlit LCD display with 17 mm high digits

Aluminium stand

USB port opposite to LCD for direct connection to PC for customize set-up (charging logo, ...)

Print-out can be with your logo in 17 languages

Integrated printer




TECHNICAL SPECIFICATIONS

Capacity:	180 kg
Sensitivity:	100 g / 0.2 lb
Size - weight:	365x490xh 1,290 mm - 15.8 kg
Platform size:	365x480xh 45 mm
Power supply:	100-240 V, 50/60 Hz

SMART WATCHES WITH 11 PARAMETERS



27232



- 11 PARAMETERS:**
 BODY TEMPERATURE
 HEART RATE
 ECG
 BPM
 SpO₂
 RESPIRATION RATE
 SLEEP
 SPORTS
 CALORIES
 STEPS
 WEATHER



27233

- 27232 QUADRANT SMART WATCH
- 27233 SPORTPIX SMART WATCH

Smart watches designed to monitor different health and sport parameters. Two types of watches with same functionalities: Sportpix is more suitable for sport activities; Quadrant for all occasions, from work to leisure. In addition to the above parameters, both watches include stopwatch, music player and photos. Wireless connection to SIMPLO App, compatible with Android 8.1 and above, iOS 13 and above.

Supplied with power cable and manual in GB, IT, FR, ES, DE available on www.gimaitaly.com. Not Medical Device.

TECHNICAL SPECIFICATIONS

	27232	27233
Accelerometer:	3D sensor	3D sensor
Display:	1.70" IPS TFT	1.3" OLED
Data sync:	4.2 BLE	4.0 BT
Temperature range:	-10° ~ 50°	-10° ~ 50°
Battery:	Lithium polymer	Lithium polymer
Battery capacity:	200 MAH	240 MAH
Battery life	stand-by: 30 days normal use: 15 days frequent use: 7 days	15 days 5 days 3 days
Charging time:	2-3 hours	2-3 hours
Weight:	55 g	38 g
Water resistance:	IP68	IP67

FITBAND ACTIVITY HEALTH TRACKERS

- 27234 FITBAND ACTIVITY HEALTH TRACKER - black
- 27237 FITBAND PLUS ACTIVITY HEALTH TRACKER - black

Waterproof intelligent bracelets for tracking different parameters.

They monitor: temperature (only 27237), heart rate, calories, blood pressure, SpO₂, sleep (only 27237), steps + distance and different sports.

Three kinds of sport modes (running, cycling, rope skipping). They also show call reminders and rejection, information push, weather, stopwatch, alarm clock and photos.

Wireless connection to Dayband App for data management. Devices requirements: iOS 8.0 and above; Android 4.4 and above. Multilingual manual: GB, FR, IT, ES, DE. Not Medical Device.

TECHNICAL SPECIFICATIONS

Display:	4 level brightness adjustment
Size - weight:	42x23x9.5 mm - 25 g
Lithium battery:	90 mAh - charging time 1.5/2 hours
Power supply:	USB rechargeable battery Liquid protection: IP67



27234

4 Health parameters
+ step & distance + 3 sport activities



27237

6 Health parameters
+ step & distance + 5 sport activities