

GIMA BIG DIAL SCALES

- **27239 GIMA BIG DIAL SCALE**
Mechanical scale with classic design: a large and clear dial, a robust housing made of steel. Supplied in multilingual gift box.



MADE IN HUNGARY

- **27240 GIMA DIGITAL BIG DIAL SCALE**
Digital scale with extra large, easy to read LCD display (65 mm). Strong and durable steel construction, comfortable platform. Lb/kg/st switchable by user. Supplied in multilingual gift box.



Multilingual manual and box:
GB, FR, IT, ES,
DE, PT, GR, Arabic

**CAPACITY
225 KG**



27240

27239

TECHNICAL SPECIFICATIONS

Capacity: 160 kg Sensitivity: 1,000 g
Minimum measuring weight: 5 kg
Dial Ø: 180 mm Weight: 3 kg
Dimension: 430x380x70 mm

TECHNICAL SPECIFICATIONS

Capacity: 225 kg Sensitivity: 100 g
Minimum measuring weight: 2 kg
Power: lithium battery (included)
Dimension: 430x380x70 mm



MECHANICAL SCALE



27230

- **27230 MECHANICAL SCALE**
Analogue bathroom scale with stable construction and plastic cover. Multilingual box and manual: GB, FR, IT, ES, PT, DE, GR, HR, HU, SI, RO, EE, LT, LV, CZ, PL, RU, Arabic.

TECHNICAL SPECIFICATIONS

Capacity: 125 kg (275 lb)
Sensitivity: 1,000 g (2 lb) Warranty: 2 years

GIMA BATHROOM DIGITAL SCALES



27248

- **27248 DIGITAL GLASS SCALE - black**
- **27249 DIGITAL GLASS SCALE - grey**
Electronic bathroom scale with 6 mm thick tempered glass platform. Ultra slim design: only 19 mm height. 4 high precision G sensors. Automatic on and off function. Lb/kg/st switchable by user.

Multilingual colour box and manual: GB, FR, IT, ES, DE, PL, PT, HU, RO, FI, SE, HR, GR, Arabic



27249

TECHNICAL SPECIFICATIONS

Capacity: max 150 kg (330 lb)
min. 5 kg (11 lb)
Sensitivity: 100 g (0.2 lb)
Dimension: 300x300x19 mm
Power: 1xCR2032 lithium battery - not included
Warranty: 2 years

EXACTA BODY FAT SCALE



Multilingual manual and box:
GB, FR, IT, ES, PT, DE, PL, NL, HU,
GR, CZ, HR, DK, SE, FI, SI, TR, RU.

27229

**UP TO
180 kg**



- **27229 EXACTA DELUXE BODY FAT SCALE**
Measuring body weight, fat, water, muscle content and daily caloric consumption. Sensor touch. Easy to read LCD display. Large platform made of safety glass. 8 user memories. Kg/lb/st switch. Auto ON/OFF function.

TECHNICAL SPECIFICATIONS

Capacity: max 180 kg (400 lb, 28 st) min 3 kg (6.6 lb) Sensitivity: 100 g (0.2 lb)
Size: 310x310x22 mm LCD Display: 31 mm
Battery: 2x CR2032 - included Warranty: 3 years

OMRON DIGITAL SCALE



27087

**UP TO
180 kg**



- **27087 OMRON HN286 DIGITAL SCALE**
Electronic bathroom scale with four integrated sensors for fast and reliable reading, up to 180 kg. Automatic shutdown after 6 seconds. Switch on by contact on the tray. Multilingual manual in all European languages and NO, TR, RU, Arabic.

TECHNICAL SPECIFICATIONS

Capacity: 5 to 180 kg Sensitivity: 100 g
Power: CR2032 lithium battery - included
Size: 300x300x18 mm Weight: 1.9 kg