

## MANUAL AND ELECTRIC LICE COMB

Natural, ecologic method to eliminate completely lice and nits without any harm. Ideal for allergy sufferers and children.

### • 24400 MANUAL LICE COMB - blister

Ergonomic comb with long metal teeth and a detachable cleaning brush with built-in magnifier. The long, delicately closed teeth of the comb are designed specifically to detach lice and their nits (tiny eggs). The round form of the touching point respects the growth direction of the hair without harming its texture.

Multilingual blister: GB, FR, IT, ES.

Free display box  
any 10 blister



24400

Detachable cleaning brush  
with built-in magnifier

### • 24403 ELECTRIC LICE COMB - IT, ES

### • 24405 ELECTRIC LICE COMB - GB, FR

This innovative product eradicates and kills lice and nits with a quick increasing electrical discharge of 100-250 Volts from the teeth, imperceptible and harmless to humans. Auto switch off. Provided of buzzer indicator and LED light. Supplied with cleaning brush, 1 AAA battery, multilingual manual in a multilingual blister.




24403, 24405

Size:

10.36x7.57x2.66 cm

## DEVICES AGAINST INSECTS



25771

CE 0051

### • 25770 ZANZA-CLICK

AFTERBITE - IT - green

### • 25771 ZANZA-CLICK

AFTERBITE - GB - blue

A few clicks to soothe itching and swelling caused by mosquito bites. MDR Certified


• 25772 BABYFRIEND  
Ultrasonic repelling device against mosquitoes. Suitable for a 15 m<sup>2</sup> room.



25772



## TICK TWEEZER



• 24390 TICK TWEEZER - 8,5 cm  
Intended for holding and removal purpose. Made of stainless steel.

## SPLINTER



• 24391 SPLINTER WITH MAGNIFYING GLASS 2X - 8 cm  
Magnifying glass splinter tweezer made of stainless steel.

## UNSTINGER®

### • 24395 UNSTINGER®

Insect poison extraction tool. An easy and painless way to avoid the discomfort that follows from bites and stings by bees, wasps, hornets and even scorpions or snakes. Reusable one-hand operated vacuum pump which easily and painlessly extracts poison and venom from a sting or bite. Multilingual box and manual: GB, FR, IT, ES.



24395

## DENTAL HYGIENE

36655

### • 36655 GERMORAL MOUTHWASH - bactericidal mouthwash One display contains 12 bottles of 300 ml

For daily oral hygiene, Germoral prevents and acts against the formation of bad breath because it stops fermentation in the oral cavity. It performs an antiplaque and hygienic action. Suitable for smokers, prosthesis or paradontic teeth users or after dental surgery for disinfecting the oral cavity. If regularly applied it prevents and acts against inflammations, gingival disorders, dental plaque formation.



## ECOSAVE FOR POISONOUS BITES


### • 25773 ECOSAVE AGAINST POISONOUS BITES

Pocket-size electronic device that produces - through two electrodes - short discharges (lasting ten microseconds approx. at one second regular intervals) able to reduce local symptoms - pain, burning, redness, swelling - caused by bites and/or stings from insects and other animals, such as hymenopterans (e.g. bees, wasps, bumble-bees, etc.), arachnids (e.g. scorpions, spiders, etc.), marine animals (e.g. fishes, jellyfishes, stinging coral), snakes venom. First aid tool, used while waiting for intervention of a doctor. Not intended for allergic reactions in allergic people. It is not effective against neurotoxic venom (some kinds of snakes, like cobra, mamba and some tropical animals). MDR Certified



25773

CE 0051



ELECTRIC STIMULATION  
TO REDUCE POISONOUS  
BITES EFFECTS