

## OMRON GENTLE TEMP® 720

### • 25611 OMRON GENTLE TEMP® 720 NON CONTACT THERMOMETER

Gentle Temp® 720 is a superfast accurate infrared non contact thermometer. It is intended to measure body, surface and room temperature. It can memorise up to 25 readings, track temperature changes over a period of time and monitor trends. In addition to the forehead mode, it can also be switched to surface measurement mode, which shows the actual and unadjusted surface temperature. It helps to determine if the object temperature (for example a milk bottle or food dispenser) is suitable for the patient or the baby. Room temperature is always displayed. It is intended for home use. Reading scale °C/°F is switchable. Supplied with 1 battery, probe cap and multilingual instructions manual: GB, FR, IT, ES, DE, NL, RU, TR, Arabic.



**OMRON**

**3 APPLICATIONS:  
FOREHEAD, SURFACE, ROOM**

#### TECHNICAL SPECIFICATIONS

Measurement type:	contactless 3 in 1 (forehead, surface, room)
Measurement range/accuracy:	forehead 34.0 - 42.2°C (93.2-108.0°F) / ± 0.2°C surface -22.0 - 80.0°C (-7.6° - 176.0°F) / ± 0.3°C room 10 - 40°C (50°-104°F) / ± 2°C
Auto power off time:	1 min
Measuring time:	1 second
Measuring distance:	1-3 cm (0.4-1.2 inch)
Temperature mode:	°C/°F
Memory:	25 readings
Display:	backlight
Battery:	3.0V DC, 1 CR2032 (code 32212) lithium button battery battery life: approx. 2,500 measurements
Classification - protection:	Type BF - IP22
Size - weight:	155x39x45 mm - 90g

## GIMATEMP NON CONTACT THERMOMETER

### • 25583 GIMATEMP NON CONTACT THERMOMETER GB, FR, IT, ES

### • 25584 GIMATEMP NON CONTACT THERMOMETER IT, DE, PT, SE

Multifunctional infrared non contact thermometer with 4 applications: body, liquid, food, ambient temperature.

- high accuracy and quick measurement in ≤2-3 sec
- 3 colour backlight alarm (red, yellow, green) for checking the temperature status

- measuring distance 1-5 cm
- room temperature is always displayed.
- 50 memories

#### THREE-COLOUR BACKLIGHT ALARM

**NORMAL**  
35.5-37.3°C (95.9-99.1°F)

**SLIGHT FEVER**  
37.4-38.0°C (99.3-100.4°F)

**HIGH FEVER**  
38.1-43.0°C (100.6-109.4°F)



**4 APPLICATIONS:  
BODY, LIQUID, FOOD,  
AMBIENT TEMPERATURE**

#### TECHNICAL SPECIFICATIONS

Measurement method:	Temporal artery scan or forehead spot
Display model - resolution:	°C/°F - 0.1°C/0.1°F
Operating mode:	Adjusted mode (Body mode) Direct mode (Surface mode) Body mode: 34.0~43.0°C/93.2-109.4°F Surface mode: 0~100°C/32-212°F
Output range:	Body mode: 35.0~42.0°C ±0.2°C / 95.0-107.6°F ±0.4°F ±0.3°C/±0.5°F for other range Surface mode: ±2°C/±3.6°F
Accuracy:	
Auto power off time:	≤18 s
Measuring time:	2-3 s
Measuring distance:	1-5 cm (0.4-2 inch)
Batteries (included):	1.5 V (AAA) alkaline batteryX2 (IEC type LR03)
Size - Weight:	122x96x42 mm - 80 g

## INFRARED NON CONTACT THERMOMETERS FOR BODY AND OBJECTS



### • 25591 INFRARED NON CONTACT THERMOMETER

Multifunctional infrared non contact thermometer for body and objects.

Effective distance 2.5 cm.

**Body temperature with fever teller**  
Wide range 32.0-42.9 °C (89.6-109 °F) and green, yellow, red backlight.

**Object temperature 0-100 °C (32-212 °F)**

Body object function for ambient temperature, milk or bath surface temperature.

#### Fever alarm

**32 memories**

**Instant response time: 1 second**

Auto/manual power off.

Supplied with 2AA batteries.

Size 160x95x46 mm.

Display 30x31 mm.

Reading scale °C/°F switchable

Multilingual box and manual: GB, FR, IT, ES.

#### FEVER TELLER



Green backlight  
from 32° to 37.5



Yellow backlight  
from 37.6° to 37.9



Red backlight  
from 38° to 42.9°



### • 25582 A200 INFRARED NON CONTACT THERMOMETER - GB, FR, IT, ES

### • 25590 A200 INFRARED NON CONTACT THERMOMETER - IT, GR, RO, PL

Multifunctional infrared non contact thermometer for body and objects.

Effective distance 2-3 cm.

- result in 1 second
- LCD backlight display

- quality sensor for high accuracy measure
- °C and °F are easily convertible

- 3 colours screen red, orange and green indicate temperature level

- automatically power off after 30s if left it on still

- 25 groups memories

- audible alarm

- low power warning

Supplied with multilingual box and manual (GB, FR, IT, ES or IT, GR, PL, RO) and batteries.

<b>Green backlight</b> MEM 36.5°C 36.5°C < 37.5°C (99.5 °F)	<b>Orange backlight</b> MEM 37.5°C 37.5°C ≥ 37.5°C ≤ 37.9°C (≥ 99.5°F ≤ 100.2°F)	<b>Red backlight</b> MEM 38.0°C 38.0°C ≥ 38.0°C (≥ 100.4 °F)
-------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------

#### TECHNICAL SPECIFICATIONS

Measurement Range:	Body 32 - 42.9°C (89.6 - 109.3 °F) ±0.2 °C Objects 0-99°C (32.0-211.82 °F)
Response time:	1 second
Memory:	25 groups
Display Mode:	°C/°F
Battery Type:	2 AAA 1.5V
Dimension (LxWxH) - Weight:	168x54x44 mm - 75 g

Body or  
objects  
switch



**25582  
25590**