


## X-CHECK B.P. MONITORS

- 49874 X-CHECK BLOOD PRESSURE MONITOR GB, FR, ES, DE, NL, SE, LT, PL, CZ, HR, RO, GR
- 49875 X-CHECK BLOOD PRESSURE MONITOR IT, PT, DK, NO, FI, LV, EE, SK, HU, SI, BG, Arabic

Accurate blood pressure monitor for home use.

- 2 users, 99 memories for each user
  - IHB detection and WHO indicator
  - average reading of last 3 measurements
  - time and date display
  - low battery indicator
  - it can work with type-C USB cable (not included)
- Supplied with cuff (22-42 cm), dustproof plastic bag, 4xAAA batteries and user manual (see above languages)
- 49877 ADULT CUFF 22-42 cm - spare

TECHNICAL SPECIFICATIONS	
Measuring method:	oscillometric
Blood pressure rate:	sys 57~255 mmHg; $\pm 3$ mmHg dia 25~195 mmHg; $\pm 3$ mmHg
Pulse rate:	40~199 bpm; $\pm 5\%$
Operating mode:	continuous operation
Display:	LCD screen, 58.5x44.5x2.9 mm
Cuff size:	22-42 cm
Size-weight:	98x118x62.5 mm - 220 g without batteries
Power supply:	4xAAA batteries DC 6V or DC 5V type C cable
Automatic switch off:	1 minute
Protection:	internally powered supply, type BF, IP21




## OMRON COMPLETE BLOOD PRESSURE AND ECG MONITOR

- 49885 OMRON COMPLETE B.P. AND ECG MONITOR HEM-7530T-E3

The clinically validated OMRON Complete is an accurate 2-in-1 smart blood pressure monitor with single lead ECG recording that delivers easy, instant and accurate results. It can measure blood pressure only, ECG only or B.P. and ECG simultaneously. Quickly screen for a range of heart conditions at once, including: AFib, tachycardia, bradycardia and hypertension. In combination with a smartphone (not included), it records, stores and transfers results. Omron Complete can be connected to OMRON connect App to consult results. The device is intended for use by healthcare professionals, patients with known or suspected heart conditions and health-conscious individuals in a general household situation. Supplied with easy cuff 22-42 cm, 4xAA batteries, storage case and multilingual manual in GB, FR, IT, ES, DE, NL, DK, SE, PL, RU.



TECHNICAL SPECIFICATIONS	
Device	
BP display:	LCD digital display
Operating mode:	continuous operation
Power supply:	4 AA alkaline batteries, 1.5 V
Battery life:	approx. 300 measurements
Size - weight:	232x123x98 mm - 550 g
Cuff size:	145x594 mm, (air tube: 750 mm)
Safety:	Type BF, Type CF applied part
Protection:	IP21, internally powered ME equipment
ECG	
Lead:	single lead
Input dynamic range:	10mV peak-to-peak
Frequency response:	0.67 Hz to 40 Hz
CMRR:	>60 dB
Input impedance:	>10 M $\Omega$
Blood pressure measurement	
Method:	oscillometric
Measurement range:	pressure: 0-299 mmHg; $\pm 3$ mmHg pulse rate: 40-180 beats/min.; $\pm 5\%$
Inflation:	automatic by electric pump
Deflation:	automatic pressure release valve
Internal memory:	90 readings



## OMRON PROFESSIONAL BLOOD PRESSURE MONITOR

- 49890 OMRON PROFESSIONAL B.P. MONITOR HBP-1120-E

The HBP-1120 is designed for use in professional settings and is clinically proven to produce fast and reliable results. Irregular heartbeats detection. It is intended for intensive daily use by healthcare professional in hospitals and medical facilities.

Zero Indicator function shows that the device is ready for next measurement. It works with AC adaptor or batteries. Class IIa. Supplied with cuff 22-32 cm, AC adapter, guarantee card and multilingual manual in GB, FR, IT, ES, DE, NL, TR, RU, Arabic.

- 49891 OMRON GS CUFF2 XS 12-18 cm - optional
- 49892 OMRON GS CUFF2 S 17-22 cm - optional
- 49893 OMRON GS CUFF2 M 22-32 cm - spare
- 49894 OMRON GS CUFF2 L 32-42 cm - optional
- 49895 OMRON GS CUFF2 XL 42-50 cm - optional

TECHNICAL SPECIFICATIONS	
Method:	oscillometric, auscultatory
Measurement range:	pressure: 0-300 mmHg; $\pm 3$ mmHg pulse rate: 40-200 beats/min.; $\pm 5\%$
Operating mode:	continuous operation
Display:	7 segment LCD
Power supply:	AC adapter
Input voltage:	AC 100-240 V; 50/60 Hz
Output voltage:	DC 6V $\pm 0.5$ V
Size - weight:	130x120x175 mm - 510 g
Protection:	Type BF applied part

