

PATIENT TRANSFER CHAIRS

- 44840 ISCHIA PATIENT TRANSFER CHAIR - blue
- 44842 ISCHIA PATIENT TRANSFER CHAIR - any colour

Transfer chair for patient transport within hospital or community environment, suited for admission, accident & emergency department (A&E) and ward transfer purposes.

It is also well suited for external transfers or when there is a need to

jump over obstacles or uneven ground smoothly.

Solid front wheels with Ø 30 cm and 2 swivelling rear castors Ø 12.5 cm. The carer can easily move the chair with patient on board, even in extremely small spaces.

With footrest fully touching the ground and backrest up, it enables the patient to safely perform any activities involving sitting or standing up.

Armrests can be flipped back to allow side access.

The backrest can be reclined via gas spring control creating a seat tilt (-33°/+6°) that allows for a better position to rest and relax.

Equipped with pushing handle, back pocket and pedal for the activation of "central braking system" and "steering-lock system" on the two rear wheels. Max weight capacity 160 kg. Class I M fire resistant.



Push handle with brake

AVAILABLE IN 24 COLOURS


Weight: 42 kg

Seat size: 470x545 mm

Footrest size: 243x380 mm

Armrest flip-back

Seat and backrest with padded surface for comfort



44840



Seat tilt -33°/+6°



Central braking system and steering-lock system

HYDRAULIC TRANSFER CHAIR

- 43430 HYDRAULIC PATIENT TRANSFER CHAIR

Hydraulic transfer chair made of steel, ideal for transferring the patient from bed to toilet, shower or other sites.

It includes 3 interchangeable seat panels:

- fixed toilet panel that can be used with or without a removable toilet pan
- removable plastic shower panel
- removable soft cushion seat panel

Shower and soft cushion seat panels can be disassembled whenever necessary to accommodate varied conditions, leaving the area more hygienic and cleaner.

- 180° easy opening and closing, flexible for movement

- hydraulic lifting adjustment up to 15 cm

- load capacity: 125 kg

- front swiveling castors with brake, one-way rear castors which prevents the machine from drifting due to unstable center of gravity when used by paralyzed people

Packing size: 79x63x33 cm.

Seat size: 48x37xh 42-57cm

3 IN 1
TRANSFER CHAIR
SHOWER CHAIR
TOILET CHAIR



43430



Toilet panel with removable toilet pan

COMMODE WHEELCHAIR

- 43202 FOLDING COMMODE WHEELCHAIR

- one-button foldable commode
- EVA seat and backrest with soft surface and odour free

- detachable seat and backrest pads

- rust-proof frame, made of aluminum alloy

- equipped with removable toilet pan

FOLDABLE

Packing size: 87x23.5x63 cm

Max. load: 125 kg

Weight: 9.4 kg



Folding footrest



43202

Castors with brake, Ø 10 cm

FOLDING BEDSIDE RAIL

- 27704 FOLDING BEDSIDE RAIL

It can be attached to the bed to assist and help those people who have difficulty moving, getting up and lying down.

One-click folding design frame: handle can be easily folded on both left and right side.

Easy to use, single hand operation.

- stepless height adjustment

- 3 levels of handle width adjustment

- comfortable PU handle lowers and raises by 180 degrees

- equipped with light for night-time use

- made of iron

Size: 64-84x63.5xh 73.5-91.5 cm. Weight: 3.95 kg

Packaging size: 60x46x36 cm

Max load: 125 kg

FOR ANY
KIND OF BED

