

ARM SLINGS

Arm sling designed for arm rest and immobilization. Ideal for immobilization after rotator cuff repair or shoulder girdle fracture, shoulder dislocation and subluxation. Breathable and durable material with adjustable padded shoulder strap. Envelope style with thumb support. Machine washable. Fits left and right. Dark blue colour.



GIMA code	ARM SLING	Length
42000	Arm sling-fabric - S	25-28 cm
42001	Arm sling-fabric - M	28-30 cm
42002	Arm sling-fabric - L	30-33 cm
42003	Arm sling-fabric - XL	33-36 cm

ELBOW BRACE WITH SLING

- 42010 UNIVERSAL ELBOW BRACE WITH SLING - right
- 42011 UNIVERSAL ELBOW BRACE WITH SLING - left

Elbow brace designed for post-operative immobilization of the elbow in 15° increments and limitation of ROM (range of motion) by adjustable dial lock hinges. Ideal for:

- conservative treatment of elbow dislocation or luxation
- stable or internally fixed fractures of the distal humerus or proximal radius or ulna.

Features:

- length adjustment (34-43 cm), one size fits all
- extension limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°
- flexion limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
- immobilization limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°
- breathable laminated soft foam material offers air permeability and moisture vapor transport

Supplied with shoulder sling for extra support and relief of arm.



KNEE BRACE

• 42026 UNIVERSAL POST OPERATION KNEE BRACE

Knee brace ideal for post operative immobilization or ROM control of ligament and fracture-related surgeries, injured or operated collateral or cruciate ligaments.

It provides excellent collateral stabilization via long side bars preventing valgus or varus deviation. Adjustable ROM by setting hinges.

Features:

- "One push button" can be released for length adjustment
 - extension limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
 - flexion limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
 - immobilization limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
 - length adjustments from 48 to 70 cm
 - transcendent overall adjustability
 - simply press and rotate settings
- Circumference: lower leg 48 cm, thigh 70 cm. Can be worn over pants or on skin.



42026 Push button for adjustment

KNEE ORTHOSIS

• 42030 CRYO PNEUMATIC KNEE ORTHOSIS

Universal cryo pneumatic knee orthosis indicated for following ACL repair, post total joint surgery, meniscal tears, patella realignment or during rehabilitation. ROM hinges with malleable aluminium uprights offering full flexion and extension control

Made of plush foam laminate with hook and loop closures, reusable gel packs feature soft fabric lining and nontoxic gel. Gel pack temperature: -20°C + 50°C. Size fits up to 61 cm circumference. Supplied with 2 ROM hinges with plastic covers and extension straps, 2 hot/cold gel packs, air compression body wrap, hand pump.



42030

COMPRESSION AND HOT-COLD THERAPY

KNEE SUPPORT

Knee support indicated for:

- mild to moderate instabilities of knee joint and strain of the collateral ligaments
 - meniscus injuries
 - chronic inflammation after operation or trauma
- It provides collateral stabilization of knee joint preventing valgus or varus deviation and hyperextension. Closed cell neoprene ensures compressive and proprioceptive effect. Blood circulation increases through warm effect.

Pull-on model with soft patellar ring that gives stabilization and prevents patellar displacement; soft-cushioned felt forms the contour of knee offering comfort; upper and lower coronal fasteners offer further fit.



42020-3

GIMA code	KNEE SUPPORT	Thigh circumf. cm	Calf circumf. cm
42020	Knee support - S	39-47	29-37
42021	Knee support - M	47-53.5	37-43
42022	Knee support - L	53.5-60	43-50
42023	Knee support - XL	60-67	50-57

ELBOW SUPPORT

• 42040 ELBOW SUPPORT WITH SILICONE PAD

Elbow support indicated for:

- articular effusion and swelling
- chronic inflammation after operation or trauma
- collateral ligament injury
- lateral and medial epicondylitis (tennis elbow and golfer's elbow)

Designed to reduce swelling from minor injuries. It provides relief at muscle insertion point through silicone pad, flexible compression and proprioception improvement. Made of seamless, four-way stretch elastic material, breathable material with moisture absorbed. Hook and loop closure for pressure adjustment.



42040

ANKLE BRACE

Ankle brace with Q-Fit closure system for excellent fit, convenience and durability. Indicated for:

- mild ankle distortion
 - chronic ankle instability
 - prophylactic and therapeutic use in sport activities
- Made of durable and tough nylon fabric with bilateral metal stays for additional support. Extra non-elastic cross-stirrup straps provide maximal stability.



42050-3

GIMA code	ANKLE BRACE
42050	Ankle brace - S
42051	Ankle brace - M
42052	Ankle brace - L
42053	Ankle brace - XL

WRIST BRACE

Q-fit wrist brace with reel locker designed for immobilization of the wrist. Ideal for:

- carpal tunnel syndrome
 - TFCC (triangular fibrocartilage complex) injury
 - radial nerve paralysis.
- Q-fit dial closure system guarantees excellent fit, convenience and durability. Comfortable, breathable skin-friendly material with anatomically-contoured aluminium palmar splint.



42060-3

Left	Right	WRIST BRACE
42060	42065	Wrist brace - S
42061	42066	Wrist brace - M
42062	42067	Wrist brace - L
42063	42068	Wrist brace - XL