

SHELL MEDICINE CABINETS



27960

• 27960 SHELL MEDICINE CABINET - white

Structure in powder-coated aluminium profiles with resistant and durable melamine-faced wood particle panels.

- 8 internal adjustable shelves and lockable cabinet
- 20 unbreakable and recyclable transparent compartment trays (10 on each door).

Edged doors revolving on 270°, hinges complete with glass-filled nylon lock and handles.

4 powder-coated steel feet at the base with adjustable plastic insoles. Size: 910x600xh 1980 mm



27960-1
Closed cabinet



27961

• 27961 SHELL PLUS MEDICINE CABINET - white

Structure made of powder-coated aluminium with melamine panels.

- 2 doors made of 20 mm thick wood particle panels edged in white ABS, revolving on 270°, hinges complete with glass-filled nylon lock and handles
- 20 unbreakable and recyclable transparent compartment trays (10 on each door)
- 5 adjustable shelves with a narcotic cabinet with lock
- 2 revolving panels with 40 unbreakable and recyclable transparent trays (20 on each panel).

4 powder coated steel feet Ø 60 mm with adjustable plastic insoles. Size: 930x600xh 1980 mm

NEO PLUS TREATMENT TROLLEY

• 45716 NEO PLUS TREATMENT TROLLEY

Treatment trolley made of ABS and antibacterial materials externally, while the internal structure is composed of steel and aluminium.

It is equipped with four castors Ø 10 cm (two with brake) and two drawers with ABS inserts, which can be freely organized. Size: 75x47.5x92 cm.



STANDARD ACCESSORIES

2 sharp containers 2 l
2 ABS waste bins with lid

Dividers
Hand sanitizer holder



45716


NEO PLUS ANESTHESIA TROLLEY

• 45718 NEO PLUS ANESTHESIA TROLLEY

Anesthesia trolley made of ABS and antibacterial materials externally, while the internal structure is composed of steel and aluminium. It is equipped with five big and deep drawers with dividers and four anti-static castors Ø 10 cm (two with brake).

Size: 75x47.5x93 cm
Drawers inner size:
2 drawers 42.4x37.5x6.8 cm
2 drawers 42.4x37.5x11 cm
1 drawer 42.4x37.5x22 cm

Central locking
Pull-out worktop



45718

STANDARD ACCESSORIES

Central lock 2 ABS waste bins with single lid
Pull-out worktop 10 anesthesia boxes
Dividers Hand sanitizer holder
Utility container Sharp container 2l



Table board covered by transparent glass
Hand sanitizer holder



NEO PLUS EMERGENCY TROLLEY

• 45720 NEO PLUS EMERGENCY TROLLEY

Emergency trolley made of ABS and antibacterial materials externally, while the internal structure is composed of steel and aluminium.

It is equipped with five big and deep drawers with dividers and 4 castors Ø 10 cm (two with brake). Size: 75x47.5x93 cm

Drawers inner size:
2 drawers 42.4x37.5x6.8 cm
2 drawers 42.4x37.5x11 cm
1 drawer 42.4x37.5x22 cm



STANDARD ACCESSORIES

IV pole Oxygen bottle holder
Defibrillator plate Utility container
Pull-out worktop CPR board
Dividers Hand sanitizer holder
Power outlet 2 sharp containers 2l
Central lock 2 ABS waste bins with lid



45720